



Technology Transition Workshop

Analysis Software for the mtDNA Tiling Assay

Tom Hall, Ph.D.

Ibis Biosciences, Inc.




Outline

- General software interface**
- Features**
- Data processing**
- Data analysis**

Technology Transition Workshop

mtDNA Analysis Software




Software Features

- ❑ **Leverages existing IbisTrack database software for sample tracking**
 - **Integrated system**
 - **MitoDNA analysis module is an independent module**
 - **Updates to MitoDNA → Simply place updated module with most recent IbisTrack; Integration is dynamic**


- ❑ **Provides functions for:**
 - **Analysis / QC of mtDNA assay data and registration of profiles**
 - **Adding / removing databases, populations, sequence profiles, and base composition profiles**
 - **Automated conversion of sequence to base count profiles**
 - **Population-based searching and reporting**
 - **Plate-based analysis reporting**

Technology Transition Workshop

Data Processing



mtDNA Processing Monitor





Technology Transition Workshop

MassCollapse Deconvolution

Raw spectrum







IMAGE COURTESY OF TOM HALL, PH.D.

MassCollapse Deconvolution

Raw spectrum





Deconvolved Spectrum





MassCollapse Deconvolution

Raw spectrum



MassCollapse Deconvolution


Raw spectrum




List of masses and intensities

Collapsed data can also be interrogated for fit of base compositions given a model of our mass tag

- 32855.1 Da**
- 34143.8 Da**
- 37058.2 Da**
- 37254.2 Da**
- 42162.1 Da**
- 42710.3 Da**



Technology Transition Workshop mtDNA Analysis Software



Basic Application Interface

The screenshot shows the IbisTrack application window. The title bar reads "IbisTrack" and the menu bar includes "File", "View", and "Help". Below the menu bar are navigation buttons for "Home", "Casework", "Databasing", "Import", and "Analysis".

Tasks

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Welcome to *IbisTrack*

IbisTrack tracks samples and plates for the Ibis T5000 Biosensor System, from initiating an assay through viewing results.

Please select from one of the following common operations, or select from the tasks in the navigation pane to the left.

- Casework**
Register an assay for casework samples. The wizard associates user defined samples with assay plates as well as provides reports for sample layout on assay plates and robot instructions for setup from a 24 place tube rack. Prior to using the wizard, the user should associate the desired control layout with the appropriate plan.
- Database Population**
Register 9 assay plates for database samples. The wizard expects 91 samples arranged in a specified plate layout for setup. The wizard automatically distributes the 91 samples to 9 assay plates with associated controls.
- Import**
Import information for assay plates and assay configurations into the local kit inventory from files provided by Ibis Biosciences.
- Analysis**
Analyze and review processed results after running assay plates on the T5000.

Tasks
Advanced Tasks
In-House Tasks

Tasks-Home 2.3.0410 MTDNA MTDNA MTDNA

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface

The screenshot shows the IbisTrack application window. The title bar reads 'IbisTrack'. The menu bar includes 'File', 'View', and 'Help'. Below the menu bar are navigation tabs for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. On the left side, there is a 'Tasks' navigation pane with a list of tasks: 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', 'STR Analysis', and 'Analyze STRs'. A blue box highlights the 'Tasks' navigation pane with the text 'Main IbisTrack program menu'. The main content area displays a 'Welcome to IbisTrack' message and a list of common operations: 'Casework', 'Database Population', 'Import', and 'Analysis', each with a brief description. At the bottom of the window, the status bar shows 'Tasks-Home' and version information '2.3.0410 MTDNA MTDNA MTDNA'.

Main IbisTrack
program menu

IMAGE COURTESY OF TOM HALL, PH.D.


Basic Application Interface

The screenshot shows the IbisTrack application window. The title bar reads 'IbisTrack'. The menu bar includes 'File', 'View', and 'Help'. Below the menu bar are navigation buttons for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. On the left is a 'Tasks' navigation pane with a list of tasks: 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The 'mtDNA Analysis' task is expanded, showing sub-tasks: 'Analyze mtDNA', 'Import Mass Data', 'Add Database Items', 'Remove Database Items', 'Move Samples to Populations', 'Rebuild Unique Products', and 'Build Amplified BC Database'. The main content area displays a 'Welcome to IbisTrack' message and instructions to select from common operations: 'Casework', 'Database Population', and 'Import'. A blue box with the text 'mtDNA Analysis menu' is overlaid on the 'Analyze mtDNA' sub-task, with an arrow pointing to it from the left. At the bottom of the window, the status bar shows 'Tasks-Home' on the left and '2.3.0410 | MTDNA | MTDNA | MTDNA' on the right.

mtDNA Analysis menu

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface



IbisTrack

File View Help Home Casework Databasing Import Analysis

Tasks

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Welcome to *IbisTrack*

IbisTrack tracks samples and plates for the Ibis T5000 Biosensor System, from initiating an assay through viewing results.

click here to select from one of the following common operations, or select from the tasks in the navigation pane to the left.

Casework
Register an assay for casework samples. The wizard associates user defined samples with assay plates as well as provides reports for sample layout on assay plates and robot instructions for setup from a 24 place tube rack. Prior to using the wizard, the user should associate the desired control layout with the appropriate plan.

Database Population
Register 9 assay plates for database samples. The wizard expects 91 samples arranged in a specified plate layout for setup. The wizard automatically distributes the 91 samples to 9 assay plates with associated controls.

Import
Import information for assay plates and assay configurations into the local kit inventory from files provided by Ibis Biosciences.

Analysis
Analyze and review processed results after running assay plates on the T5000.

Tasks
Advanced Tasks
In-House Tasks

Tasks-Home 2.3.0410 MTDNA MTDNA MTDNA

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface

The screenshot displays the IbisTrack application window. The title bar reads "IbisTrack". The menu bar includes "File", "View", and "Help". The toolbar contains icons for "Home", "Casework", "Databasing", "Import", and "Analysis".

Tasks Panel (Left):

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Main Interface:

At the top, there are filters: "Plate" (choose plate), "From" (30-Oct-2007), "To" (13-Nov-2007), and "Analyze" (New only). Below this is a table with columns numbered 1 through 9 and rows labeled A and B. Row A, column 1 is highlighted in blue.

Below the table is a navigation bar with tabs: Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup.

At the bottom of the window, the status bar shows "Tasks: Analyze mtDNA" on the left and "2.3.0410 | MTDNA | MTDNA | MTDNA" on the right.

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface

The screenshot displays the IbisTrack application interface. On the left is a 'Tasks' sidebar with categories like 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The main window features a menu bar (File, View, Help) and a toolbar with icons for Home, Casework, Databasing, Import, and Analysis. Below the toolbar, there are search and filter fields for 'Plate' and 'Sample', along with date pickers for 'From' (30-Oct-2007) and 'To' (13-Nov-2007), and a 'New only' checkbox. A table with columns 1-9 and rows A-B is visible. A callout box with the text 'Set time stamp filters' points to the 'From' date field. Below the table is a navigation bar with tabs for Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup. The status bar at the bottom shows 'Tasks: Analyze mtDNA' and version information '2.3.0410 MTDNA MTDNA MTDNA'.

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface




IMAGE COURTESY OF TOM HALL, PH.D.

Technology Transition Workshop

Basic Application Interface

The screenshot displays the IbisTrack software interface. At the top, there is a menu bar with options: File, View, Help, Home, Casework, Databasing, Import, and Analysis. Below the menu bar is a search and filter area with a 'Plate' dropdown set to 'P00005555', a 'From' date of '30-Oct-2005', a 'To' date of '13-Nov-2007', and an 'Analyze' button. A 'Sample' dropdown is set to 'All' and a 'New only' checkbox is present.

The main data area is a table with columns numbered 1 through 9. The rows are labeled 'A' and 'B'. The data in the table is as follows:

	1	2	3	4	5	6	7	8	9	
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004	HUM2925 Scenario: 65 M11819004

Below the table is a navigation bar with tabs: Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup. The 'Analysis' tab is currently selected.

On the left side, there is a 'Tasks' panel with expandable sections: Register Assays, Manage Inventory, View Plates, mtDNA Analysis (with sub-options: Analyze mtDNA, Import Mass Data, Add Database Items, Remove Database Items, Move Samples to Populations, Rebuild Unique Products, Build Amplified BC Database), and STR Analysis (with sub-option: Analyze STRs).

At the bottom of the interface, there is a status bar showing 'Tasks-Analyze mtDNA' on the left and '2.3.0410 MTDNA MTDNA MTDNA' on the right.

IMAGE COURTESY OF TOM HALL, PH.D.

Basic Application Interface

The screenshot shows the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar is a navigation bar with 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The main area is divided into a left sidebar and a central grid.

Left Sidebar (Tasks):

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Central Grid:

At the top of the grid, there are filters: 'Plate' (P00005555), 'From' (30-Oct-2005), 'To' (13-Nov-2007), and 'Sample' (All). There is also an 'Analyze' button and a 'New only' checkbox.

The grid contains 10 columns and 2 rows of sample data. The first cell in the first row is highlighted in blue. A callout box with the text 'Select sample' points to this cell.

	1	2	3	4	5	6	7	8	9	10
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804	HUM2928 Scenario: 86 M11819804

At the bottom of the grid, there are tabs: 'Analysis', 'Profiles', 'Database', 'Mass Data', 'Preferences', 'Scenarios', 'Monthly report', 'Report archive', and 'Databasing plate setup'.

At the bottom of the window, there is a status bar with the text 'Tasks-Analyze mtDNA' on the left and '2.3.0410 MTDNA MTDNA MTDNA' on the right.

IMAGE COURTESY OF TOM HALL, PH.D.

Technology Transition Workshop

Basic Application Interface

The screenshot displays the IbisTrack application window. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below it is a toolbar with icons for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The main interface is divided into a left-hand 'Tasks' sidebar and a central data area.

Tasks Sidebar:

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Main Data Area:

At the top of the main area, there are input fields for 'Plate' (P00005555), 'From' (30-Oct-2005), and 'To' (13-Nov-2007). There is an 'Analyze' button and a 'New only' checkbox. Below these is a table with 10 columns and 2 rows (A and B). The first two cells in each row are highlighted in red.

	1	2	3	4	5	6	7	8	9	10
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70	HUM2925 Scenario: 66 HUM2901 Scenario: 70

Below the table is a horizontal menu with tabs: Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup.

At the bottom of the window, there is a status bar showing 'Tasks-Analyze mtDNA' on the left and '2.3.0410 | MTDNA | MTDNA | MTDNA' on the right.

IMAGE COURTESY OF TOM HALL, PH.D.

Technology Transition Workshop

Basic Application Interface

The screenshot shows the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar is a navigation bar with 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. A callout box with the text 'Press Analyze' points to the 'Analyze' button in the navigation bar.

The main area of the interface is a data table. The table has columns numbered 1 through 9. The rows are labeled 'A' and 'B'. The data in the table is as follows:

	1	2	3	4	5	6	7	8	9	
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924	HUM2925 Scenario: 66 HUM2924

Below the table, there are several tabs: 'Analysis', 'Profiles', 'Database', 'Mass Data', 'Preferences', 'Scenarios', 'Monthly report', 'Report archive', and 'Databasing plate setup'. The 'Analysis' tab is currently selected.

At the bottom of the interface, there is a status bar with the text 'Tasks-Analyze mtDNA' on the left and '2.3.0410 | MTDNA | MTDNA | MTDNA' on the right.

IMAGE COURTESY OF TOM HALL, PH.D.

Technology Transition Workshop

Basic Analysis Interface

The screenshot shows the IbisTrack software interface with several key components highlighted by blue callout boxes:


- Plate selection:** Points to the 'Plate' dropdown menu at the top left, currently set to 'P00005556'.
- Sample selection:** Points to the 'Sample' dropdown menu below the plate, currently set to '1'.
- Plate time-stamp selection:** Points to the 'From' and 'To' date pickers, showing '30-Oct-2005' and '13-Nov-2007' respectively.
- Interface navigation:** Points to the top menu bar containing 'File', 'View', 'Help', 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'.
- Sample / well navigation:** Points to the 'Well 1 (A01)' dropdown menu on the right side of the interface.
- Product assignment info:** Points to the 'Well 13 (B01)' mass spectrum plot, which shows peaks at 2926 + 2891 + 2908.
- Base count profile:** Points to the 'Well 13 (B01)' mass spectrum plot, which shows peaks at 2899 + 2890 + 2907.
- Mass assignment detail:** Points to the 'num error exp mass obs mass' table at the bottom right, showing a match for mass 32854.625.
- Coverage map:** Points to the green bar chart on the left side of the interface, representing signal intensity across different wells.
- Spectral schematic views:** Points to the 'Well 61 (F01)' mass spectrum plot at the bottom, showing peaks at 2917 + 2903 + 2916.

Technology Transition Workshop

Basic Utilities Interface

The screenshot shows the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar are navigation tabs for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The main content area is titled 'Welcome to IbisTrack' and contains the following text: 'IbisTrack tracks samples and plates for the Ibis T5000 Biosensor System, from initiating an assay through to processing results. Please select from one of the following common operations, or select from the tasks in the navigation pane to the left.' Below this text are four sections: 'Casework' (with a yellow folder icon), 'Database Population' (with a blue database icon), 'Import' (with a CD icon), and 'Analysis' (with a magnifying glass icon). On the left side, there is a 'Tasks' navigation pane with a list of tasks: 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis' (with a sub-list: 'Analyze mtDNA', 'Import MtDNA', 'Add Database', 'Remove Database', 'Move Samples to Conditions', 'Rebuild Unique Profiles', 'Build Amplified BC Database'), and 'STR Analysis' (with a sub-list: 'Analyze STRs'). A blue box with the word 'Click' is overlaid on the 'Welcome' text, with several arrows pointing to the 'Analyze mtDNA' and 'Analyze STRs' tasks in the navigation pane. At the bottom of the interface, there are tabs for 'Tasks', 'Advanced Tasks', and 'In-House Tasks'. The status bar at the very bottom shows 'Tasks-Home' on the left and '2.3.0410 | MTDNA | MTDNA | MTDNA' on the right.

IMAGE COURTESY OF TOM HALL, PH.D.



Technology Transition Workshop

Basic Utilities Interface

The screenshot displays the IbisTrack software interface. At the top, five yellow boxes with arrows point to specific menu items: 'Item addition' (pointing to 'Add items to database'), 'Item deletion' (pointing to 'Remove items'), 'Sample arrangement' (pointing to 'Arrange samples'), 'Database maintenance' (pointing to 'Rebuild Unique Products Table' and 'Build Amplified BC database'), and 'Mass data import' (pointing to 'Import mass file'). Below these, a 'Tasks' pane shows a list of actions including 'Sequence profile', 'Individual IDs', 'Base count profile', 'Add database', 'Add population', and 'Import XML'. A large yellow text box in the center lists four methods for sequence imports: 1. Manual profile entry (or copy/paste), 2. Import of raw sequence (software aligns and creates profile automatically), 3. Import LISA-exported profiles (e.g. AFDIL system), and 4. Import CMF-formatted profiles (can be done in batch). Below this list, it states 'Conversion to base composition profiles is automatic for all.' In the bottom left, another yellow box says 'Example: Single addition of a sequence profile'. The bottom right shows a dialog box for adding a sequence profile with a table of entries and instructions: 'Use uppercase for UIPAC codes: 263 N, 263 R, etc. Press 'Format ->' and double-check row entries. Press 'Register into database''. The status bar at the bottom reads 'Tasks-Add Database Items', 'IMAGE COURTESY OF TOM HALL, PH.D.', and '2.5.0520 MTDNA MTDNA MTDNA'.

Item addition **Item deletion** **Sample arrangement** **Database maintenance** **Mass data import**

File View Help Home Casework Databasing Import Analysis

Tasks Add items to database Remove items Arrange samples Rebuild Unique Products Table Build Amplified BC database Import mass file

Register Assays Sequence profile Individual IDs Base count profile Add database Add population Import XML

Sequence imports can be done by:

1. Manual profile entry (or copy/paste)
2. Import of raw sequence (software aligns and creates profile automatically)
3. Import LISA-exported profiles (e.g. AFDIL system)
4. Import CMF-formatted profiles (can be done in batch)

Conversion to base composition profiles is automatic for all.

Example: Single addition of a sequence profile

-	524.3	A
-	524.4	C
T	16224	C
T	16311	C
T	16519	C
A	16524	C

Use uppercase for UIPAC codes:
263 N
263 R
etc.
Press 'Format ->' and double-check row entries.
Press 'Register into database'

Tasks-Add Database Items IMAGE COURTESY OF TOM HALL, PH.D. 2.5.0520 MTDNA MTDNA MTDNA

The screenshot displays the IbisTrack software interface. On the left is a 'Tasks' sidebar with categories: Register Assays, Manage Inventory, View Plates, mtDNA Analysis (with sub-tasks: Analyze mtDNA, Import Mass Data, Add Database Items, Remove Database Items, Move Samples to Populations, Rebuild Unique Products, Build Amplified BC Database), and STR Analysis (with sub-task: Analyze STRs). The main window has a menu bar (File, View, Help) and a toolbar (Home, Casework, Databasing, Import, Analysis). Below the toolbar are filters for 'Plate' (choose plate), 'Sample' (choose sample), 'From' (30-Oct-2007), 'To' (13-Nov-2007), and 'Analyze' (New only). A table with 9 columns and 2 rows (A, B) is visible, with the first cell in row A highlighted. Below the table is a navigation bar with tabs: Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup. The status bar at the bottom shows 'Tasks-Analyze mtDNA', 'IMAGE COURTESY OF TOM HALL, PH.D.', and version '2.3.0410'.

The screenshot displays the IbisTrack software interface. The main window has a menu bar with 'File', 'View', and 'Help'. Below the menu bar are tabs for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The 'Analysis' tab is active, showing a 'Tasks' sidebar on the left with categories like 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The main workspace shows a 'Plate' and 'Sample' selection area with dropdown menus for 'Plate' and 'Sample', and date pickers for 'From' (30-Oct-2007) and 'To' (13-Nov-2007). A blue callout box with the text 'Set time stamp filters' points to the date pickers. Below the date pickers is a grid with columns 1-9 and rows A and B. At the bottom of the window, there is a status bar with 'Tasks-Analyze mtDNA', 'IMAGE COURTESY OF TOM HALL, PH.D.', and version information '2.3.0410 MTDNA MTDNA MTDNA'.

The screenshot shows the IbisTrack software interface. The main window has a menu bar (File, View, Help) and a toolbar with icons for Home, Casework, Databasing, Import, and Analysis. On the left, there is a 'Tasks' sidebar with categories like Register Assays, Manage Inventory, View Plates, mtDNA Analysis, and STR Analysis. The main area displays a 'Sample' list with columns for Plate, Sample ID, and dates. The sample 'P00005556' is selected, and a callout box with an arrow points to it, containing the text 'Select an assay plate'. Below the list are tabs for Analysis, Profiles, Database, Mass Data, Preferences, Reports, Monthly report, Report archive, and Databasing plate setup. The status bar at the bottom shows 'Tasks-Analyze mtDNA', 'IMAGE COURTESY OF TOM HALL, PH.D.', and version information '2.3.0410 MTDNA MTDNA MTDNA'.

Plate	Sample	From	To	Analysis
P00005556	P00005552	30-Oct-2007	13-Nov-2007	<input type="checkbox"/>
	P00005553			<input type="checkbox"/>
	P00005554			<input type="checkbox"/>
	P00005555			<input type="checkbox"/>
	P00005556			<input checked="" type="checkbox"/>
	P00005557			<input type="checkbox"/>
	P00005558			<input type="checkbox"/>
	P00005559			<input type="checkbox"/>

The screenshot displays the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar are navigation tabs for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The main window is titled 'IbisTrack' and contains a 'Tasks' sidebar on the left and a central data grid.

Tasks Sidebar:

- Register Assays
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Main Window Controls:

- Plate: P00005555
- From: 30-Oct-2005
- To: 13-Nov-2007
- Sample: All
- Buttons: Analyze, New only

Data Grid:

	1	2	3	4	5	6	7	8	9	
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86	HUM2925 Scenario: 86

Navigation Tabs: Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, Databasing plate setup

Footer: Tasks-Analyze mtDNA | IMAGE COURTESY OF TOM HALL, PH.D. | 2.3.0410 MTDNA MTDNA MTDNA


The screenshot shows the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar are tabs for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The main window is titled 'IbisTrack' and contains a 'Tasks' sidebar on the left with options like 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The central area displays a grid of sample analysis results. The grid has columns numbered 1 through 9 and rows labeled 'A' and 'B'. Each cell in the grid contains text such as 'HUM2906 Scenario: 75' or 'HUM2901 Scenario: 70'. A blue box with the text 'Select sample' and an arrow points to the first cell in row 'A', column 1. Below the grid are tabs for 'Analysis', 'Profiles', 'Database', 'Mass Data', 'Preferences', 'Scenarios', 'Monthly report', 'Report archive', and 'Databasing plate setup'. At the bottom of the window, there is a status bar with the text 'Tasks-Analyze mtDNA' on the left and '2.3.0410 MTDNA MTDNA MTDNA' on the right.

The screenshot displays the IbisTrack software interface. The top menu bar includes File, View, and Help. Below it are navigation tabs for Home, Casework, Databasing, Import, and Analysis. A left-hand sidebar contains a 'Tasks' menu with options like Register Assays, Manage Inventory, View Plates, mtDNA Analysis, and STR Analysis. The main workspace shows a search filter for Plate P00005555 and Sample 1, with dates from 30-Oct-2005 to 13-Nov-2007. A grid of analysis results is shown, with columns numbered 1-9 and rows labeled A and B. The first cell in row A is highlighted in red. Below the grid are tabs for Analysis, Profiles, Database, Mass Data, Preferences, Scenarios, Monthly report, Report archive, and Databasing plate setup. The status bar at the bottom indicates 'Tasks-Analyze mtDNA' and version '2.3.0410'.

	1	2	3	4	5	6	7	8	9	
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501	HUM2925 Scenario: 86 HUM19501

IMAGE COURTESY OF TOM HALL, PH.D.

The screenshot shows the IbisTrack software interface. At the top, there is a menu bar with 'File', 'View', and 'Help'. Below the menu bar is a navigation bar with 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The 'Analysis' button is highlighted with a callout box that says 'Press Analyze'. Below the navigation bar is a 'Tasks' sidebar on the left with various options like 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The main area shows a table with columns for 'Plate' and 'Sample'. The 'Plate' is set to 'P00005555' and the 'Sample' is set to '1'. The table contains two rows of data, each with 10 columns. The first row is labeled 'A' and the second row is labeled 'B'. The first column of each row is highlighted in red. Below the table is a navigation bar with 'Analysis', 'Profiles', 'Database', 'Mass Data', 'Preferences', 'Scenarios', 'Monthly report', 'Report archive', and 'Databasing plate setup'. At the bottom of the window, there is a status bar with 'Tasks-Analyze mtDNA', 'IMAGE COURTESY OF TOM HALL, PH.D.', and '2.3.0410 MTDNA MTDNA MTDNA'.



ppm = Parts per million: 1 ppm = expected mass / 1,000,000

ppm error = ((expected - observed)/expected) * 1,000,000

Tasks

- Register Assays
 - Register Casework Plates
 - Register New Plate for Database Population
 - Process Repeats for Database
 - Register Quality Control Plates
 - Generate Plate Setup Reports
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Analysis View

Plate: P00005556 From: 18-Nov-2005 To: 02-Dec-2007

Sample: 1

Analysis Profiles Database Mass Data Preferences Scenarios Monthly report Report archive

P00005556-2003X02368-0001A1D (1)

P00005556_E01 P00005556_F01 P00005556_G01 P00005556_H01 P00005556-2003X02368-0001A1D composite

Filter ambiguous assignments Remove primers from profile

Spectrum viewing mode

Raw Deconvolved

Single assignment - top strand

Product strand corresponds to a mass peak

Single assignment - bottom strand

Well: 73 (G01)
Scenario: 65, PP 2896 (HUMMTDNA_ASN_161)
Positions: 16102..16224
Amplified coordinates: 16124..16201
Base count: A45 G13 C41 T24

num	error	exp. mass	obs. mass
1	0.7492	37567.2673	37567.2391

Comments for analysis report

Tasks Analyze mtDNA

IMAGE COURTESY OF TOM HALL, PH.D.

2.3.0410 MTDNA MTDNA THALL

The screenshot shows the IbisTrack software interface. The top menu bar includes File, View, Help, Home, Casework, Databasing, Import, and Analysis. The left sidebar contains a 'Tasks' panel with sections for Register Assays, Manage Inventory, View Plates, mtDNA Analysis, and STR Analysis. The main workspace displays a plate layout with 8 wells (A1-H8) and a detailed view of a 'Single assignment - top strand' for well A1. The 'Profile record' section shows 'Primer pair', 'CRS coords', and 'base count'. The 'Pair of product strands corresponds to a double-stranded PCR product' section shows a chromatogram with peaks at 2896 and 2912. The 'Pair of masses corresponds to a double-stranded PCR product' section shows a table with columns for 'num', 'error', 'exp. mass', and 'obs. mass', with a row containing '1', '0.749', '12.538', and '11.925'. The bottom of the interface includes a 'Comments for analysis report' field and a status bar with 'Tasks-Analyze mtDNA', 'IMAGE COURTESY OF TOM HALL, PH.D.', and '2.3.0410 | MTDNA | MTDNA | THALL'.

Single assignment – top strand

Single assignment – bottom strand

Profile record

Primer pair

CRS coords

base count

Pair of product strands corresponds to a double-stranded PCR product

Pair of masses corresponds to a double-stranded PCR product

2896: 16102..16224: A45 G13 C41 T24

num error exp. mass obs. mass

1 0.749 12.538 11.925

Comments for analysis report


Tasks-Analyze mtDNA

IMAGE COURTESY OF TOM HALL, PH.D.

2.3.0410 | MTDNA | MTDNA | THALL

Technology Transition Workshop


Single Sample Analysis



Three products in a well correspond to three roughly evenly-spaced products

Use of a Mass Tag


High mass precision and mass tag combine to provide unambiguous base compositions in routine operation.



Use of a Mass Tag

Without mass tag:

Product strands differ by 1 Da for two products that differ by a G→A and C→T SNP at the same time.




Use of a Mass Tag

With mass tag:

With ^{13}C -dGTP, the mass separation increases to ~ 10 Da for each strand
This is about 300 ppm (we generally see errors ~ 10 ppm).


A mass tag increases mass separation for these SNPs

 = Corresponding untagged peaks



Real Data with Mass Tag


Primer pair 2906 product from two representative samples
(Data for primer pair 2906 are overlaid).



What Does Mass Tagging Mean?


For all
combinations
of base
differences
plus or minus
10 A, G, C, T ...

For any base count



What Does Mass Tagging Mean?

For all
combinations
of base
differences
plus or minus
10 A, G, C, T ...



What Does Mass Tagging Mean?

Potential assignment ambiguities increase dramatically with measurement error and product size.




IMAGE COURTESY OF TOM HALL, PH.D.

What Does Mass Tagging Mean?

Ambiguous assignments not a problem in our normal working range of mass measurement.

 Working range for TOF (<20 ppm)





IMAGE COURTESY OF TOM HALL, PH.D.

What Does Mass Tagging Mean?

**Most assignments would be fine with natural nucleotides.
Some assignments fall outside of “normal” range.**

 Working range for TOF (<20 ppm)


Which 6% of measurements are a problem?



Percentages are based upon 123,415 mass assignments

Mass Tag Changes Rules

With ^{13}C -dGTP, the mass separation increases to ~10 Da for each strand
This is about 300 ppm (we generally see errors ~10 ppm).



What Does Mass Tagging Mean?

Using ^{13}C -dGTP eliminates ambiguity with a wide margin around our working zone.

- Working range for TOF (<20 ppm)
- 80 ppm measurement error maximum in this assay






IMAGE COURTESY OF TOM HALL, PH.D.


What Does Mass Tagging Mean?

Using ^{13}C -dGTP eliminates ambiguity with a wide margin around our working zone.

 Working range for TOF (<20 ppm)

 80 ppm measurement error maximum in this assay

Threshold at 70 ppm to be conservative




NIJ Technology Transition Workshop

National Institute of Justice

Spectral Viewer




IMAGE COURTESY OF TOM HALL, PH.D.

Used to examine data quality, signal, noise, or to take a picture of the spectrum.




Spectral Viewer - Deconvolved

Zoom view can give an idea of data quality surrounding a given peak.



Base Composition Browser

IMAGE COURTESY OF TOM HALL, PH.D.



Base Composition Browser

Primer pairs in current well are integrated to this view.

Base composition browser ...

Primer pair
2906
2906
2901
2892

Existing base count from database

A G C T Forward Mass Reverse Mass

Clear overlays -->

Check -->


C-length

Add->


Sodium adducts

Potassium adducts

Base Composition Browser




Base Composition Browser



All known products for selected primer pair are in this list.

Base Composition Browser



Base Composition Browser

Base composition browser ...

Primer pair
2906

Forward Primer A5 G2 C8 T9 TCCTTTATCGCACCTACGTTCAAT
Reverse Primer A2 G9 C1 T11 TGGTTGTTATGATGTCTGTGTGG

Existing base count from database

A G C T Forward Mass Reverse Mass
48 18 31 40 42187.0191 42854.099

Clear overlays -->
Check -->

Check +A --> Check +A -->

Sodium adducts
Potassium adducts

Add->

Adenylation products, sodium adducts, and potassium adducts can also be looked for here.

Technology Transition Workshop

Compare Sample to Database

The screenshot shows the IbisTrack software interface. The main window displays a coverage map for a sample labeled 'P00005556-2003X02368-0001A1D composite'. The coverage map consists of a series of green bars of varying heights, representing signal intensity at different positions. A blue box highlights a region of the coverage map with the text 'Right-click in coverage map region'. A context menu is open over this region, listing various actions such as 'Copy coverage map to Clipboard', 'Copy profile', and 'Compare profile pairwise to a profile in the database'. The interface also includes a left-hand navigation pane with sections like 'Register Assays', 'Manage Inventory', 'View Plates', 'mtDNA Analysis', and 'STR Analysis'. The top of the window shows a menu bar with 'File', 'View', and 'Help', and a toolbar with buttons for 'Home', 'Casework', 'Databasing', 'Import', and 'Analysis'. The status bar at the bottom indicates 'Tasks>Analyze mtDNA' and 'IMAGE COURTESY OF TOM HALL, PH.D.'.

Compare Sample to Database

Side-by-side profile comparison shown product-by-product differences between two profiles.

Profile comparison: P00005556-2003X02368-0001A1D (1)

Compare to another stored profile | Compare to another profile on the plate | Search database

Compare Ignore C-length in poly-C regions Tab delimited Force match Display missing entries at bottom

Add relative DS abundances Collapse matching / nonmatching entries for each primer pair

Database: MEASURED
 Population: UNKNOWN
 Individual: 2003X02368-0001A1D
 Reaction set: Tiling set

pp	coords	2003X02368-0001A1D	2003X02368-0001A1D	Diff
2901:	15893..16012:	A47 G18 C25 T30	----> A47 G18 C25 T30	0
2925:	15937..16041:	A35 G14 C24 T32	----> A35 G14 C24 T32	0
2899:	15985..16073:	A26 G15 C21 T27	----> A26 G15 C21 T27	0
2898:	16025..16119:	A26 G17 C27 T25	----> A26 G17 C27 T25	0
2897:	16055..16155:	A31 G13 C29 T28	----> A31 G13 C29 T28	0
2896:	16102..16224:	A45 G13 C41 T24	----> A45 G13 C41 T24	0
2895:	16130..16224:	A36 G7 C33 T19	----> A36 G7 C33 T19	0
2893:	16154..16268:	A44 G7 C46 T18	----> A44 G7 C46 T18	0
2892:	16231..16338:	A40 G9 C40 T19	----> A40 G9 C40 T19	0
2891:	16256..16366:	A37 G9 C41 T24	----> A37 G9 C41 T24	0
2890:	16318..16402:	A20 G14 C30 T21	----> A20 G14 C30 T21	0
2889:	16357..16451:	A21 G17 C36 T21	----> A21 G17 C36 T21	0
2902:	5..97:	A20 G23 C24 T26	----> A20 G23 C24 T26	0
2903:	20..139:	A24 G34 C29 T33	----> A24 G34 C29 T33	0
2904:	83..187:	A23 G21 C30 T31	----> A23 G21 C30 T31	0
2905:	113..245:	A39 G18 C29 T47	----> A39 G18 C29 T47	0
2906:	154..290:	A48 G18 C31 T40	----> A48 G18 C31 T40	0
2908:	204..330:	A42 G16 C38 T32	----> A42 G16 C38 T32	0
2907:	239..363:	A43 G11 C49 T23	----> A43 G11 C49 T23	0
2923:	262..390:	A47 G10 C53 T20	----> A47 G10 C53 T20	0
2910:	331..425:	A33 G9 C27 T26	----> A33 G9 C27 T26	0
2916:	367..463:	A27 G8 C32 T30	----> A27 G8 C32 T30	0
2912:	409..521:	A32 G7 C48 T26	----> A32 G7 C48 T26	0
2913:	464..603:	A44 G10 C63 T23	----> A44 G10 C63 T23	0
Total differences: 0				

Identical profile

Compare Sample to Database

Side-by-side profile comparison shown product-by-product differences between two profiles.

Profile comparison: P00005556-2003X02368-0001A1D (1)

Compare to another stored profile | Compare to another profile on the plate | Search database

Compare Ignore C-length in poly-C regions Tab delimited Force match Display missing entries at bottom

Add relative DS abundances Collapse matching / nonmatching entries for each primer pair

Database: MEASURED
 Population: UNKNOWN
 Individual: 2003X02369-0001A1D
 Reaction set: Tiling set

pp	coords	2003X02368-0001A1D	2003X02369-0001A1D	Diff
2901:	15893..16012:	A47 G18 C25 T30	----> A47 G18 C25 T30	0
2925:	15937..16041:	A35 G14 C24 T32	----> A35 G14 C24 T32	0
2899:	15985..16073:	A26 G15 C21 T27	----> A25 G16 C21 T27	1
2898:	16025..16119:	A26 G17 C27 T25	----> A25 G18 C27 T25	1
2897:	16055..16155:	A31 G13 C29 T28	----> A31 G13 C29 T28	0
2896:	16102..16224:	A45 G13 C41 T24	----> A45 G13 C41 T24	0
2895:	16130..16224:	A36 G7 C33 T19	----> A36 G7 C33 T19	0
2893:	16154..16268:	A44 G7 C46 T18	----> A44 G7 C45 T19	1
2892:	16231..16338:	A40 G9 C40 T19	----> A40 G9 C40 T19	0
2891:	16256..16366:	A37 G9 C41 T24	----> A37 G9 C41 T24	0
2890:	16318..16402:	A20 G14 C30 T21	----> A20 G14 C30 T21	0
2889:	16357..16451:	A21 G17 C36 T21	----> A21 G17 C36 T21	0
2902:	5..97:	A20 G23 C24 T26	----> A19 G24 C24 T26	1
2903:	20..139:	A24 G34 C29 T33	----> A24 G34 C29 T33	0
2904:	83..187:	A23 G21 C30 T31	----> A23 G21 C29 T32	1
2905:	113..245:	A39 G18 C29 T47	----> A40 G18 C28 T47	1
2906:	154..290:	A48 G18 C31 T40	----> A49 G18 C31 T39	1
2908:	204..330:	A42 G16 C38 T32	----> A42 G16 C38 T32	0
2907:	239..363:	A43 G11 C49 T23	----> A43 G11 C49 T23	0
2923:	262..390:	A47 G10 C53 T20	----> A47 G10 C53 T20	0
2910:	331..425:	A33 G9 C27 T26	----> A33 G9 C27 T26	0
2916:	367..463:	A27 G8 C32 T30	----> A27 G8 C32 T30	0
2912:	409..521:	A32 G7 C48 T26	----> A32 G7 C49 T25	1
2913:	454..588:	A44 G10 C63 T23	----> A43 G10 C62 T23	1

Total differences: 9

Different profile

NIJ Technology Transition Workshop

National Institute of Justice Search Database

IbisTrack

File View Help Home Casework Databasing Import Analysis

Tasks

- Register Assays
 - Register Casework Plates
 - Register New Plate for Database Population
 - Process Repeats for Database
 - Register Quality Control Plates
 - Generate Plate Setup Reports
- Manage Inventory
- View Plates
- mtDNA Analysis
 - Analyze mtDNA
 - Import Mass Data
 - Add Database Items
 - Remove Database Items
 - Move Samples to Populations
 - Rebuild Unique Products
 - Build Amplified BC Database
- STR Analysis
 - Analyze STRs

Tasks: Analyze mtDNA

Plate: P00005556 From: 18-Nov-2005 To: 02-Dec-2007 Analyze

Sample: 1 To: 02-Dec-2007 New only

Analysis Profiles Database Mass Data Preferences Scenarios Monthly report Report archive

P00005556-2003X02368-0001A1D (1) P00050052-USA.FBI.000072.SAMPLE (11)

P00005556, E01 P00005556, F01 P00005556, G01 P00005556, H01 P00005556-2003X02368-0001A1D composite

Filter ambiguous assignments Remove primers from profile

Spectrum viewing mode

Raw Deconvolved

P00005556-2003X02368-0001A1D composite

Well: 73 (G01)
Scenario: 65, PP 2896 (HUMMTDNA_ASN_1610)
Positions: 16102..16224
Amplified coordinates: 16124..16201
Base count: A45 G13 C41 T24

2901: 15893..16012: A47 G18 C25 T30
2925: 15937..16041: A35 G14 C24 T32
2899: 15985..16073: A26 G15 C21 T27
2898: 16025..16119: A26 G17 C27 T25
2897: 16055..16155: A31 G13 C29 T28
2896: 16102..16224: A45 G13 C41 T24


num error exp. mass obs. mass

- Copy coverage map to Clipboard
- Copy profile
- Copy profile with abundances
- Copy profile as table with average abundances
- Reset hidden fragments
- Update data paths for raw spectra
- Compare profile pairwise to a profile in the database
- Search database with current profile**
- Compare profile to another profile on this plate
- Register current sample as a contaminant profile
- Set sample status to PASS
- Set sample status to FAIL
- Set sample status to REDO
- Add sample to repeat list
- Remove from repeat list
- Re-evaluate controls and samples on plate
- Override sample status and allow profile to be registered
- Reanalyze current sample with all single SNP variants allowed
- Copy relative abundances of products in sample (per primer pair) to clipboard

Comments for analysis report

2.3.0410 MTDNA MTDNA THALL

IMAGE COURTESY OF TOM HALL, PH.D.



Profile comparison: P00005556-2003X02368-0001A1D (1)

Compare to another stored profile | Compare to another profile on the plate | Search database

Search Show profiles up to 5 differences

Ignore C-length in poly-C regions Use uniform query range in Preferences

Database	Population	Inc
AFDIL	(Select all)	<input type="checkbox"/>
AFDIL_QC	(Select all)	<input type="checkbox"/>
CONTAMINANTS	(Select all)	<input type="checkbox"/>
Fake Database	(Select all)	<input type="checkbox"/>
FBI	(Select all)	<input type="checkbox"/>
GENOMES	(Select all)	<input type="checkbox"/>
IBIS	(Select all)	<input type="checkbox"/>
MEASURED	(Select all)	<input type="checkbox"/>
NIST	(Select all)	<input type="checkbox"/>

Max number of differences to show in interface

Target databases

C-length heteroplasmy treatment

Profile comparison: P00005556-2003X02368-0001A1D (1)

Compare to another stored profile | Compare to another profile on the plate | Search database

Search Show profiles up to 5 differences

Ignore C-length in poly-C regions Use uniform query range in Preferences

1 *	Database	Population	Inc
<input type="checkbox"/>	AFDIL	(Select all)	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Turkmenistan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Uzbekistan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Vietnam	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	African American	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Hispanics	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Greece	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Budapest Hungary	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Russia	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Kyrgyzstan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Kazakhstan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Cyprus	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Asian	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Hungarian Romanies	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Caucasians	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Hong Kong Chinese	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Afghanistan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Egypt	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Nairobi Kenya	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Jordan	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Indonesia	<input checked="" type="checkbox"/>
<input type="checkbox"/>	AFDIL	Tajikstan	<input checked="" type="checkbox"/>
<input checked="" type="checkbox"/>	AFDIL_QC	(Select all)	<input type="checkbox"/>
<input checked="" type="checkbox"/>	CONTAMINANTS	(Select all)	<input type="checkbox"/>
<input checked="" type="checkbox"/>	Fake Database	(Select all)	<input type="checkbox"/>
<input checked="" type="checkbox"/>	FBI	(Select all)	<input type="checkbox"/>
<input checked="" type="checkbox"/>	GENOMES	(Select all)	<input type="checkbox"/>
<input checked="" type="checkbox"/>	IBIS	(Select all)	<input type="checkbox"/>

DB Population ID

Target database

Target populations

Target populations choice

AF-2

USA.FBI.000070


diff

2901: 15893..16012: no data
 2925: 15937..16041: no data
 2899: 15985..16073: no data
 2898: 16048..16098: A15 G7 C17 T12
 2897: 16078..16129: A15 G8 C14 T15
 2896: 16124..16201: A27 G7 C27 T17
 2895: 16157..16201: A17 G1 C19 T8
 2893: 16182..16250: A24 G5 C28 T12
 2892: 16254..16305: A17 G4 C22 T8 N
 2891: 16283..16344: A23 G5 C18 T15 N
 2890: 16318..16402: no data
 2889: 16357..16451: no data
 2902: 5..97: no data
 2903: 20..139: no data
 2904: 103..162: A10 G10 C20 T20
 2905: 138..217: A23 G9 C21 T27
 2906: 178..267: A32 G15 C15 T28
 2908: 234..313: A30 G6 C30 T16
 2907: 263..340: A25 G6 C35 T14
 2923: 262..390: no data
 2910: 331..425: no data
 2916: 367..463: no data
 2912: 409..521: no data
 2913: 464..603: no data



2901: 15893..16012: no data
 2925: 15937..16041: no data
 2899: 16015..16051: A7 G10 C5 T15
 2898: 16048..16098: A15 G7 C16 T13
 2897: 16078..16129: A15 G8 C14 T15
 2896: 16124..16201: A26 G7 C29 T16
 2895: 16157..16201: A16 G1 C21 T7
 2893: 16182..16250: A23 G5 C29 T12
 2892: 16254..16305: A18 G4 C22 T8
 2891: 16283..16344: A24 G5 C19 T14
 2890: 16342..16381: A8 G5 C17 T10
 2889: 16377..16428: A12 G9 C20 T11
 2902: 31..76: A5 G15 C12 T15
 2903: 41..114: A12 G24 C21 T18
 2904: 103..162: A9 G11 C18 T22
 2905: 138..217: A22 G10 C19 T29
 2906: 178..267: A32 G15 C15 T28
 2908: 234..313: A30 G6 C31 T16
 2907: 263..340: A25 G6 C36 T14
 2923: 289..367: A27 G6 C37 T12
 2910: 355..401: A21 G3 C14 T9
 2916: 389..437: A11 G7 C12 T19
 2912: 431..501: A20 G1 C34 T16
 2913: 493..576: A27 G6 C46 T5

1
 2
 2
 1
 3
 3




Total minimum differences

6

Technology Transition Workshop Profile Registration

Analysis of a full plate exposes profile registration interface.

- Samples that pass auto-processing are checked by default
- Controls are not checked by default
- Samples that failed will be disabled
- If positive or negative control fails, registration for entire plate will be disabled; disabled state can be overridden



The screenshot displays the profile registration interface. On the left, a grid shows samples organized by row (A-F) and column (1-4). Each cell contains a sample ID and its associated scenario. A red border highlights the first cell in row A, column 1. To the right of the grid is a 'Register' button and a list of controls with checkboxes. The controls listed are FBI-3, FBI-9, FBI-22, FBI-28, FBI-32, FBI-33, FBI-37, FBI-47, FBI-48, FBI-49, WATERCTRL, and WATERCTRL. The first nine controls have their checkboxes checked. A status dialog box is open in the foreground, titled 'Status' with a close button (X). It contains an information icon (i) and the following text: '532 products successfully registered.', '0 failures.', and '0 products already registered.' An 'OK' button is at the bottom of the dialog.

	1	2	3	4
A	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70	HUM2906 Scenario: 75 HUM2901 Scenario: 70
B	HUM2925 Scenario: 86 HUM2891 Scenario: 60	HUM2925 Scenario: 86 HUM2891 Scenario: 60	HUM2925 Scenario: 86 HUM2891 Scenario: 60	HUM2925 Scenario: 86 HUM2891 Scenario: 60
C	HUM2899 Scenario: 68 HUM2890 Scenario: 59	HUM2899 Scenario: 68 HUM2890 Scenario: 59	HUM2899 Scenario: 68 HUM2890 Scenario: 59	HUM2899 Scenario: 68 HUM2890 Scenario: 59
D	HUM2898 Scenario: 67 HUM2889 Scenario: 58	HUM2898 Scenario: 67 HUM2889 Scenario: 58	HUM2898 Scenario: 67 HUM2889 Scenario: 58	HUM2898 Scenario: 67 HUM2889 Scenario: 58
E	HUM2893 Scenario: 62 HUM2910 Scenario: 79	HUM2893 Scenario: 62 HUM2910 Scenario: 79	HUM2893 Scenario: 62 HUM2910 Scenario: 79	HUM2893 Scenario: 62 HUM2910 Scenario: 79
F	HUM2897 Scenario: 66 HUM2903 Scenario: 72	HUM2897 Scenario: 66 HUM2903 Scenario: 72	HUM2897 Scenario: 66 HUM2903 Scenario: 72	HUM2897 Scenario: 66 HUM2903 Scenario: 72

Register

- FBI-3
- FBI-9
- FBI-22
- FBI-28
- FBI-32
- FBI-33
- FBI-37
- FBI-47
- FBI-48
- FBI-49
- WATERCTRL
- WATERCTRL

Status X

532 products successfully registered.
0 failures.
0 products already registered.

OK

IMAGE COURTESY OF TOM HALL, PH.D.

NIJ Technology Transition Workshop

Report Generation

National Institute of Justice

Click to generate a report

Generate report for P20000034

Log general plate comments here

Log sample-specific comments here

Comments for analysis report

This plate was actually run on 2/7/2006. Tom somehow managed to mess up his calculation of DNA concentration after running the Quantifiler assay. The target DNA concentration in all reactions was 500 pg. Because of the calculation mishap, the actual DNA concentrations used were:

Sample	pg/rxn
AF-12	85.6
AF-13	62.7
AF-14	46.4
AF-15	2.0
AF-16	120.5
AF-17	9.6

Comments for analysis report

85.6 pg per reaction used instead of 500 pg per reaction.

Comparison to other profiles on plate P20000034:
 1: AF-12: 13 min differences.
 2: AF-13: 11 min differences.
 3: AF-14: 0 min differences.

Tasks: Register Assays, Manage Inventory, View Plates, mtDNA Analysis, STR Analysis

Analysis Profiles Database Mass Data Preferences Scenarios Monthly report Report archive

P20000034-AF-12 (1) P20000034-AF-13 (2) P20000034-AF-14 (3) P20000034-AF-15 (4) P20000034-AF-16 (5) P20000034-AF-17 (6) P20000034-AF-18 (7) P20000034-AF-19 (8)

P20000034, B03 P20000034, C03 P20000034, D03 P20000034, E03 P20000034, F03 P20000034, G03 P20000034, H03 P20000034-AF-14 composite

Filter ambiguous assignments Remove primers from profile

Spectrum viewing mode Raw Deconvolved

P20000034-AF-14 composite

Well 3 (A03)
2906 + 2901 + 2892

Well 15 (B03)
2925 + 2891 + 2908

num	error	exp. mass	obs. mass	mod
2901:	15893..16012:	A46 G19 C25 T30		
2925:	15937..16041:	A35 G14 C24 T32		
2899:	15985..16073:	A26 G15 C21 T27		
2898:	16025..16119:	A26 G17 C28 T24		
2897:	16055..16155:	A32 G12 C30 T27		
2896:	16102..16224:	A46 G12 C41 T24		

Tasks: Advanced Tasks, In-House Tasks

Tasks-Analyze mtDNA

IMAGE COURTESY OF TOM HALL, PH.D.

2.3.0410 MTDNA MTDNA THALL



Technology Transition Workshop

Mitochondrial DNA Plate Analysis Report

Plate: P00008859 (Mag = M00004214, PCR = C00014215)

Date: Sunday, May 27, 2007 9:19:56 PM

User: thall

Experiment: NIST -SAMPLES

TTrack database: mtdna

mtDNA database: mtdna

First page of report summarizes samples and controls on a plate.

- Quantity (or dilution) used
- Status
- Ave abundance (for controls)
- Ave error (for positive)
- Control status
- Samples that have registered profiles
- User comments

12 samples on plate P00008859:

	Status
1. NIST-GT37869 (500.00 pg)	PASS
2. NIST-TT51399 (500.00 pg)	PASS
3. NIST-ZT80737 (500.00 pg)	PASS
4. NIST-ZT80870 (500.00 pg)	PASS
5. NIST-UT57317 (500.00 pg)	PASS
6. NIST-WT51359 (500.00 pg)	PASS
7. NIST-WA29584 (500.00 pg)	PASS
8. NIST-MT94875 (500.00 pg)	PASS
9. NIST-OT05896 (500.00 pg)	PASS
10. NIST-PT84222 (500.00 pg)	PASS
11. PDB (1.00 DF)	PASS
12. SC35495-10-POS (500.00 pg)	PASS

Positive control: 1

12: SC35495-10-POS (500.00 pg)

Average product abundance: 2698.771

Mean measurement error magnitude: 11.09 ppm

Status: PASS

Positive control status: PASS

Negative control: 1

11: PDB

Number of detected products: 0

Average product abundance: NaN

Ratio of average product abundance to corresponding positive control products: NaN

Status: PASS

Negative control status: PASS

10 samples, 1 positive, 1 negative.

Base counts are registered for NIST-GT37869, NIST-TT51399, NIST-ZT80737, NIST-ZT80870, NIST-UT57317, NIST-WT51359, NIST-WA29584, NIST-MT94875, NIST-OT05896, NIST-PT84222

User comments:

Plate #9 of NIST samples sent by John Butler.

Sample NIST-ZT80737 should be repeated to confirm the odd G-length heteroplasmy observed in overlapping primer pairs 2902 and 2903.



Technology Transition Workshop

Second page of report summarizes analysis parameters used.

Plate: P00008859 (Mag = M00004214, PCR = C00014215)
Date: Sunday, May 27, 2007 9:19:56 PM
User: thall
Experiment: NIST-SAMPLES

TTrack database: mtdna
mtDNA database: mtdna

Analysis parameters for plate P00008859

Matching threshold: 70 ppm
Matching error delta threshold for double-stranded product: 80 ppm
Look for adenylations: OFF
Max ratio of strand abundance in double-stranded product: 2.5
Max negative to positive abundance ratio: 0.1
Min sample to positive abundance ratio: 0.25
Min fraction of ave peak abundance to report unmatched mass: 0.25
There is a filter in place to filter out minor SNP heteroplasmy assignments where the minor variant is less than 15% of the abundance of the major product.
There is a filter in place to filter out apparent +A artifacts where the minor variant is less than 30% of the abundance of the major product.

¹³C-enriched dGTP was used on plate P00008859:
Lot 4: 99.07% ¹³C enrichment

Criteria for passing the positive control:

All expected products for the positive control sample must be detected and the sample must match 100% to the expected positive control profile.

Criteria for passing the negative control:

If no positive control is present on the plate, the average abundance of any detected products in the negative control cannot exceed 10% of the average abundance of the corresponding products in all passed samples on the plate. If there is a positive control on the plate, the average abundance of any detected products in the negative control cannot exceed 10% of the average abundance of the corresponding products in the positive control.

Regardless of average abundance in the negative control, if identifiable products are observed in more than 1/2 of the primer pair regions, the negative control will be failed.

If either the negative or positive control fails, the plate must be rerun.

Criteria for passing a sample:

1. A double-stranded product must be detected for all primer pairs.
2. The average product abundance for the sample must equal or exceed 25% of the average for the positive control.

Each sample gets a summary page.


- Coverage map
- Average abundance
- Average error
- Status
- Comparison to other samples on the plate
- Log of manual assignments and filtered assignments (if any)
- User comments

Plate: P00008859 (Mag = M00004214, PCR = C00014215)
 Date: Sunday, May 27, 2007 9:19:56 PM
 User: thall
 Experiment: NIST-SAMPLES

TTrack database: mtdna
 mtDNA database: mtdna

Sample 1: NIST-GT37869 (500.00 pg)

P00008859-NIST-GT37869 composite



2901: 15924..15985: A25 G9 C14 T14
 2925: 15963..16017: A22 G4 C13 T16
 2899: 16015..16051: A7 G10 C5 T15
 2898: 16048..16098: A15 G7 C15 T14
 2897: 16078..16129: A15 G8 C15 T14
 2896: 16124..16201: A28 G6 C29 T15
 2895: 16157..16201: A17 G1 C20 T7
 2893: 16182..16250: A24 G5 C27 T13
 2892: 16254..16305: A18 G4 C22 T8
 2891: 16283..16344: A24 G5 C19 T14
 2890: 16342..16381: A8 G5 C17 T10
 2889: 16377..16428: A12 G9 C20 T11
 2902: 31..76: A5 G16 C10 T15
 2903: 41..114: A12 G25 C19 T18
 2904: 103..162: A10 G10 C18 T22
 2905: 138..217: A23 G9 C18 T30
 2906: 178..267: A32 G15 C13 T30
 2908: 234..313: A30 G6 C27 T18
 2907: 263..340: A25 G6 C33 T15
 2923: 289..367: A27 G6 C34 T13
 2910: 355..401: A21 G3 C14 T9
 2916: 389..437: A11 G7 C12 T19
 2912: 431..501: A20 G1 C34 T16
 2913: 493..576: A27 G6 C46 T5

Well 1 (A01)
2906 + 2901 + 2892

Well 13 (B01)
2925 + 2891 + 2908

Well 25 (C01)
2899 + 2890 + 2907

Well 37 (D01)
2898 + 2889 + 2923

Well 49 (E01)
2893 + 2910 + 2902

Well 61 (F01)
2897 + 2903 + 2916

Well 73 (G01)
2896 + 2913 + 2904

Well 85 (H01)
2905 + 2895 + 2912

System Comments:
 Average abundance: 2129.774
 Mean measurement error: 4.86 ppm
 Mean measurement error magnitude: 12.1 ppm
 24 detected products


Analysis Status: PASS
 Reg status: COMPLETE: last updated 5/27/2007 9:17:09 PM

Comparison to other samples on plate:
 Comparison to other profiles on plate P00008859:
 1: NIST-GT37869: 0 min differences.
 8: NIST-MT94875: 9 min differences.
 9: NIST-OT05896: 11 min differences.
 10: NIST-PT84222: 12 min differences.
 2: NIST-TT51399: 12 min differences.
 5: NIST-UT57317: 10 min differences.
 7: NIST-WA29584: 9 min differences.
 6: NIST-WT51359: 9 min differences.
 3: NIST-ZT80737: 8 min differences.
 4: NIST-ZT80870: 12 min differences.
 12: SC35495-10-POS: 10 min differences.

User Comments:

Summary of raw data signal corresponding to primers in each well of negative control(s).

If there are any products in the negative control, a comparison to each sample (one sample per page) will follow.





Technology Transition Workshop

Summary

- Data processing after MS is automated**
- Processed data persist in database with registered sample and plate information**
- Sample analysis is automated (triggered manually)**
- Interface provides dynamic graphical functions to evaluate, QC, and adjust (if necessary) profile assignments**
- Profiles can be registered a plate at a time**
- Automated comparison of profiles to database entries**
- Search database**
- Add sequences and convert to base compositions in the database**



Technology Transition Workshop

Contact Information:

Tom Hall, Ph.D.

Ibis Biosciences, Inc.

thall@ibisbio.com