




Technology Transition Workshop | *Manohar Furtado, Ph.D.*

GenPlexTM : A CE-based SNP Detection Platform


Summary of SNP Detection Methods



Technology
Transition Workshop



SNP / Mutation Detection



Sample Types Work-Flow

Sample Types

Swabs
Blood
Serum
Semen
Saliva
Tissue
Feces
Urine
Nasal wash
Soil
Food matrices
Meat
Milk
Rinses
Vegetables
Beverage
Creams
Vaccines
Injectables
Water
Air

Sample Pre-processing

Sample
Collection
Devices

Sample
Clean-up

Sample
Conc.
Devices

Nucleic Acid Isolation Station

6100
GenePod
Enigma
IQuum, HandyLab®
Network Bio, Microchip
Robotics, KF, Tecan
Irlekan, IGene

Amplification Detection Systems

Real-time Systems
7500, m2000, Aztec
Bax® Q7, Biotrove,
Fluidigm®, HandyLab®,
Enigma, IQuum

Thermal Cyclers
Fragment Analysis
Sequencing

CE Systems
31XX
Network Biosystems
Micro Chip/Mathies

Luminex®

E-TaqMan®
Mass Spec

Software Integration: Sample ID & tracking, LIMS, Data Analysis & Reporting, Archiving

Technology
Transition Workshop



Mutation Detection Assay Formats

Customer Needs

- **Low sample input (0.5 – 1 ng range), sensitivity**
- **Heterozygote detection / balance**
- **Accuracy, specificity, reproducibility**
- **Inhibitor tolerance**
- **Degraded samples / small amplicons**
- **Contamination controls**
- **Work-flow (ease of use)**
- **Detection software**
- **Assay cost**
 - **Ease of manufacturing**

Mutation Detection Assay Formats

- **Single-plex PCR and sequencing (dye-labeled terminator) CE**
 - **MicroSeq[®] Fungal ID kits (330-500 bp)**
 - **ViroSeq[®] Genotyping System (FDA cleared) (1200 bp)**
 - **Detailed analysis of heterozygote / mixture detection capability for over 70 mutations**
 - **Detailed analysis of inhibitor effects**
 - **VariantSeQr[™] (Cancer-related genes; > 65,000 amplicons)**
 - **MitoSeQr[™] (Mito-sequencing)**

Mutation Detection Assay Formats

- **Multiplex-PCR and SBE**
 - SNaPShot[®] (dye-labeled terminator) CE
 - PinPoint (ddNTPs / no dye) MALDI
- **Fundamental issue with dye-labeled terminators are incorporation efficiencies of polymerases on hetero / mixture detection capabilities**
 - **Mutation K103N**
 - Detected at 5 to 10% mutant with a > 95% accuracy
 - **Mutation D67N**
 - Detected at 60% mutant with a < 90% accuracy

Mutation Detection Assay Formats

- **Multiplex-PCR and OLA on CE / Luminex™ / Qdot® Nanocrystals / Microarrays**
 - Cystic fibrosis kit (FDA Cleared)
 - GenPlex® format for human ID etc.
- **Multiplex-OLA and PCR**
 - SNPLex® (variable SNP sets for disease association) CE
 - Universal TaqMan® (real-time allelic discrimination)

Mutation Detection Assay Formats

- **Allelic discrimination (multiplexed on arrays)**
 - **TaqMan[®] Arrays (48 to 384 single-plexes)**
 - **Open array formats (Mustang, Biotrove, Fluidigm[®])**
- **Whole genome sequencing**
 - **The ultimate solution: Alec Jeffreys SOLiD[™] System**

Other Formats

- **Real-time PCR (Taqman[®] and SYBR Green)**
 - Gene expression profiling
 - Pathogen detection
 - Biothreat applications
- **Allelic discrimination**
 - Genotyping on TaqMan[®]
 - HCV genotyping (sequencing versus OLA versus TaqMan[®])
- **Allele-specific PCR formats**
- **Endonuclease cleavage**
- **Microarray-based formats**
- **LAMP**
- **FRET probes and melt curves**
- **Invader assay**
- **Scorpion probes (ARMS)**


Addressing Sample Requirements

- **1 ng of human DNA = 300 copies**
- **Sequencing formats**
 - Need about 300 to 900 target genome copies
 - HIV 1000 copies (larger amplicon)
- **Multiplex PCR formats (SBE, PCR – OLA formats)**
40- to 50-plex
 - Short amplicons multiplexed: 150 to 300 target genome copies
 - < 1 ng for a 50-plex


Addressing Sample Requirements

- **Methods requiring much higher sample input**
 - **SNPLex[®] (OLA-PCR) 30 to 40 ng / 48-plex reaction**
 - **Universal detector (OLA-TaqMan[®]) 3 ng / SNP or about 150 ng for a 50-plex**
 - **Allelic discrimination formats (high sample input)**
 - **TaqMan[®] 600 to 900 copies / SNP: 100 ng for a 50-plex**
 - **Open arrays (300 copies / SNP): 50 ng for a 50-plex**

SeqScape[®] Analysis for PCR Template Titration Study




Effect of Genomic DNA Amount on Heterozygote SNP Calling and Quality of DNA Sequence



DNA Titration: Universal Detector


DNA NA17202

3 ng / SNP




Total: 144 ng

0.3 ng / SNP



14.4 ng



0.03 ng / SNP



1.44 ng

Technology
Transition Workshop 

Protocol Development for 49-Plex Panel: A Sensitive Assay at 1 ng Input



Technology
Transition Workshop



Pre-amplification

Tested for TaqMan[®] assays:

ITN	1000-plex	
Body fluid ID	28-plex	<i>heterozygote balance</i>
Sequencing	24-plex	
Identifiler	17-plex	



Split sample for 10 assays



Limited sample: about 20 copies



10 cycle pre-amplification



Split sample for 1000 assays



TaqMan[®] detection limit of 5 to 10 copies / assay



2 copies / assay



20 copies / assay

As many as 100 to 1000 targets


Sample about 20,000 copies

Carry-over Contamination in PCR

- **Exponential amplification followed by exponential amplification**
 - Nested PCR **needs containment**
 - Pre-amplification by PCR followed by TaqMan[®]
 - Pre-amplification by PCR followed by PCR-OLA or OLA PCR
- **Exponential amplification followed non-exponential amplifications / other detection formats**
 - PCR amplification followed by exo-SAP then cycle sequencing
 - PCR amplification followed by exo-SAP the OLA
 - PCR amplification followed by single-based extensions
- **Detection instrument dynamic range**
 - CE about 2 logs
 - TaqMan[®] about 6 logs
 - Luminex[®] about 4 logs

Complete Solution: Sample to Answer


SNPs / Indels and Target Detection



Sample Preparation Systems

- **Direct lysis in PCR compatible buffers**
- **Column-based**
 - Rapid spin column
 - LySep column
- **PrepSEQ™ & PrepFiler™ Kits**
 - Manual
 - KingFisher®
 - Tecan
 - Automate Express







2007

Technology
Transition Workshop



Chemistry Overview – The GenPlex™ Genotyping System



Technology
Transition Workshop



Pros and Cons

- **Pros**

- High sensitivity
- Rapid time-to-result turn around time, < 7 hours
- OLA reaction targets PCR products at high concentration; not gDNA; less complexity
- Very useful for fixed / known sets of SNPs / mutations (eg. CF kit)
- Highly specific; specificity at the PCR and OLA steps
- Flexible; multiple applications possible


- **Cons**

- < 1 ng gDNA
- Multiplexed PCR needs some development
- Not as useful for screening variable sets of SNPs needed for disease association studies

Three Unique OLA Probes per SNP Assay



OLA Principle




Technology
Transition Workshop







Strep hybridization




Strep elution



Capillary Electrophoresis

GenPlex™ Chemistry and Workflow: Current Format



System features

- High sensitivity
 - 100 pg gDNA
- Highly specific
 - Specificity at the PCR and OLA steps
- Rapid time-to-result turn around time
- Automated
 - About 5 hours
- Multiple applications possible
 - 48 SNPs
 - 96 DNA targets
 - 48 targets and 24 SNPs
 - Any combination

TTR ~5 hrs



Technology
Transition Workshop

GenPlex™ Chemistry and Workflow: Other Formats

Eliminate the capture and wash steps

Multiplex PCR

OLA

Run on CE

Existing CF format

Run on Prodigy™
microarrays


Run on Luminex®
bead-based



About 4 hrs

Technology
Transition Workshop 

GenPlex™ High Throughput Solution



Biomek-3000 and GenPlex SNP Genotyping in Forensic Genetics

Michael Stangegaard,* Carmen Tomas, Anders J. Hansen, Rune Frank-Hansen, Claus Børsting, and Niels Morling
University of Copenhagen, Copenhagen, Denmark

Stangegaard, M. et al., *JALA* (2008) 13(5) 297-303

Technology Transition Workshop 

SNPforID Consortium

Electrophoresis 2006, 27, 1713–1724

1713

Juan J. Sanchez¹
Chris Phillips²
Claus Børsting¹
Kinga Balogh³
Magdalena Bogus³
Manuel Fondevila²
Cheryl D. Harrison⁴
Esther Musgrave-Brown⁴
Antonio Salas²
Denise Syndercombe-Court⁴
Peter M. Schneider³
Angel Carracedo²
Niels Morling¹

¹Department of Forensic Genetics,
Institute of Forensic Medicine,
University of Copenhagen,
Copenhagen, Denmark

²Institute of Legal Medicine,
University of Santiago de
Compostela,
Santiago de Compostela, Spain

³Institute of Legal Medicine,
Lebanese American University

Research Article

A multiplex assay with 52 single nucleotide polymorphisms for human identification

A total of 52 SNPs reported to be polymorphic in European, Asian and African populations were selected. Of these, 42 were from the distal regions of each autosome (except chromosome 19). Nearly all selected SNPs were located at least 100 kb distant from known genes and commonly used STRs. We established a highly sensitive and reproducible SNP-typing method with amplification of all 52 DNA fragments in one PCR reaction followed by detection of the SNPs with two single base extension reactions analysed using CE. The amplicons ranged from 59 to 115 bp in length. Complete SNP profiles were obtained from 500 pg DNA. The 52 loci were efficiently amplified from degraded samples where previously only partial STR profiles had been obtained. A total of 700 individuals from Denmark, Greenland, Somalia, Turkey, China, Germany, Taiwan, Thailand and Japan were typed, and the allele frequencies estimated. All 52 SNPs were polymorphic in the three major population groups. The mean match probability was at least 5.0×10^{-19} in the populations studied. Typical paternity indices ranged from 336 000 in Asians to 549 000 in Europeans. Details of

Sanchez, J. et al., *Electrophoresis* (2006) **27** 1713-1724

Initial Strategy

**Select p-arm and
q-arm SNPs on
each autosome**



**Pilot multiplex designs
with 22 - 25plex SNaPshot™**

**Re-validate SNP allele
frequencies**




**Review and modify the first 50
(on the basis of PCR performance and frequency)**

Custom GenPlex™ Kit for SNP and Indel Genotyping: 49-Plex




Custom Kit

Working prototype in the forensic area: SNPforID panel



SNP Genotyping Using GeneMapper® Software 4.0




Allelic ladder

Singleplex

Multiplex
(49-plex)

Technology
Transition Workshop 

Automated Allele Calling – ZipChute[®] Probe Identification



e.g. Hom. 1
AA


Het.
AG

Hom. 2
GG


Marker = SNP
Marker = 2 bins
Bin = Allele

Gender Determination

DNA1




DNA2



Polyacrylamide Gel QC for Multiplex PCR Products

Targets 48 Multiplex Targets

1 2 3 4 5 6 7 8 M



DNA:

1, 3, 5, 7 = 14529

2, 4, 6, 8 = 09947A

Primers:

1, 2, 5, 6 = SNPforID

3, 4, 7, 8 = 1.2 μ M

Annealing temperature:



1, 2, 3, 4 = 57° C

5, 6, 7, 8 = 60° C

60 bp

Low molecular weight PCR products were observed

Annealing Temperature Testing: 57° C and 60° C



DNA: NA14529

Accuracy Evaluation Studies Across Multiple Platforms



Available online at www.sciencedirect.com



Forensic Science International: Genetics 1 (2007) 180–185



www.elsevier.com/locate/bsig

Evaluation of the Genplex SNP typing system and a 49plex forensic marker panel

C. Phillips^{a,*}, R. Fang^b, D. Ballard^c, M. Fondevila^a, C. Harrison^c, F. Hyland^b,
E. Musgrave-Brown^c, C. Proff^d, E. Ramos-Luis^a, B. Sobrino^a, A. Carracedo^a,
M.R. Furtado^b, D. Syndercombe Court^c, P.M. Schneider^d

The SNPforID Consortium

^aForensic Genetics Department, Genomic Medicine Group, University of Santiago de Compostela, Galicia, Spain

^bApplied Markets Group, Applied Biosystems, Foster City, CA, USA

^cDepartment of Haematology, ICMS, Queen Mary's School of Medicine & Dentistry, London E1 2AT, UK

^dInstitute of Legal Medicine, University of Cologne, Germany

Received 29 January 2007; accepted 3 February 2007

100% Concordance:

- TaqMan[®]
- SNaPshot[™]
- Sequenom[®] iPLEX[™]

Phillips, C. et al., *For Sci International Genetics* (2007) **1(2)** 180-185



High Accuracy and Sensitivity

Forensic Science International: Genetics 3 (2008) 1–6



Contents lists available at ScienceDirect

Forensic Science International: Genetics

journal homepage: www.elsevier.com/locate/bsif



Sensitivity: 250 pg
Error Rate: 1/13728
(caused by mutation)

Typing of 48 autosomal SNPs and amelogenin with GenPlex SNP genotyping system in forensic genetics

Carmen Tomas*, Michael Stangegaard, Claus Børsting, Anders Johannes Hansen, Niels Morling

The SNPforID Consortium

Section of Forensic Genetics, Department of Forensic Medicine, Faculty of Health Sciences, University of Copenhagen, 11 Frederik V's Vej, DK-2100 Copenhagen, Denmark

ARTICLE INFO

Article history:

Received 11 April 2008

Received in revised form 19 June 2008

Accepted 30 June 2008

Keywords:

GenPlex
Autosomal SNPs
Amelogenin
Forensic genetics

ABSTRACT

GenPlex (Applied Biosystems) is a new SNP genotyping system based on an initial PCR amplification followed by an oligoligation assay (OLA). The OLA consists of the hybridization of allele and locus specific oligonucleotides (ASOs and LSOs) to PCR products and posterior ligation of ASOs and LSOs. The ligation products are immobilized to microtitre plates and reporter oligonucleotides (ZipChute® probes) are hybridized to the ligation products. ZipChute® probes are subsequently eluted and detected using capillary electrophoresis. Applied Biosystems developed the GenPlex SNP genotyping system with amelogenin and 48 of the 52 SNPs used in the 52 SNP-plex assay developed by the SNPforID consortium. The system requires equipment that is usually found in forensic genetic laboratories. The use of a robot for performance of the pipetting steps is highly recommendable.

A total of 286 individuals from Denmark, Somalia and Greenland were investigated with GenPlex using a Biomek® 3000 (Beckman Coulter) robot. The results were compared to results obtained with an ISO 17025 accredited SNP typing assay based on single base extension (SBE). With the GenPlex SNP genotyping system, full SNP profiles were obtained in 97.6% of the investigations. Perfect concordance was obtained in duplicate investigations and the SNP genotypes obtained with the GenPlex system were concordant with those of the accredited SBE based SNP typing system except for one result in rs901398 in one of 286 individuals most likely due to a mutation 6 bp downstream of the SNP. Reproducible SNP genotypes were obtained from as little as 250 pg of DNA.

© 2008 Elsevier Ireland Ltd. All rights reserved.

Tomas, C. et al., *For Sci International Genetics* (2008) **3(1)** 1–6

Technology
Transition Workshop



Mixture and Degraded DNA Sample SNP Genotyping



ELSEVIER

Available online at www.sciencedirect.com



Forensic Science International: Genetics Supplement Series 1 (2008) 389–393



www.elsevier.com/locate/FSIGSS

Research article

Forensic validation of the Genplex SNP typing system—Results of an inter-laboratory study

Esther Musgrave-Brown^{a,*}, David Ballard^a, Manuel Fondevila Álvarez^b, Rixun Fang^c, Cheryl Harrison^a, Chris Phillips^b, Yogesh Prasad^d, Bea Sobrino Rey^b, Catherine Thacker^a, Joerg Wiluhn^d, Angel Carracedo^b, Peter M. Schneider^e, Denise Syndercombe Court^a

The SNPforID Consortium¹

^aCentre for Haematology, ICMS, Barts & the London, Queen Mary's School of Medicine & Dentistry, London, United Kingdom

^bInstitute of Legal Medicine, University of Santiago de Compostela, Santiago de Compostela, Spain

^cApplied Biosystems, Foster City, USA

^dApplied Biosystems, Spain

^eInstitute of Legal Medicine, University of Cologne, Cologne, Germany

Received 6 September 2007; accepted 11 October 2007

Abstract

We present data from a multi-laboratory validation study of the Genplex typing system [C. Phillips, et al., Evaluation of the Genplex SNP typing system and a 49plex forensic marker panel, *Forensic Sci. Int.: Genet.* 1 (2) (2007) 180–185.] (Applied Biosystems), which interrogates a subset of 48 SNPs selected from the panel of 52 previously developed for forensic analysis by the SNPforID consortium [J.J. Sanchez, et al., A multiplex assay with 52 single nucleotide polymorphisms for human identification, *Electrophoresis* 27 (9) (2006) 1713–1724], plus amelogenin. The Genplex technology was developed through modification of the SNPlex™ system (also Applied Biosystems) and utilises oligo-ligation of PCR products followed by probe hybridisation to generate dye-labeled, allele-specific oligonucleotides that are detected with capillary electrophoresis.

We compare the success rate of Genplex in typing 55 samples in three laboratories with results obtained from STR typing of the same samples using Powerplex® 16 (Promega) and SGMPlus® (Applied Biosystems). The sample set was chosen to mimic extracts encountered in forensic situations and includes low concentration and degraded material as well as mixtures. We demonstrate that the Genplex technique provides a significantly higher success rate than STR-based methods when typing degraded DNA and is also capable of mixture detection up to a level of 1 in 10.

© 2008 Elsevier Ireland Ltd. All rights reserved.

Keywords: SNP; Oligo-ligation; Genplex

Musgrave-Brown, E. et al., *For Sci International Genetics Suppl. Series 1* (2008) 389–393

Demonstrate:

- Higher success rate than STR-based methods when typing degraded DNA
- Mixture detection up to a level of 1 in 10

Technology
Transition Workshop



Automating GenPlex™




Tomas, C. et al., *For Sci International Genetics* (2008) **3(1)** 1-6

Stangegaard, M. et al., *JALA* (2008) **13(5)** 297-303


- **286 individuals typed**
- **Biomek 3000 robot used for automation**
- **Sensitivity at 250 pg**
- **No carry-over contamination**
- **High Accuracy: only 1/14014 allele discrepancy due to a mutation 6 bp downstream of the SNP**

Profiles of 3130xl and 3730xl Platforms

3130xl




3730xl



NTCs Were Clean



AB

SNP
forID







PHC

Protocol Development for 49-Plex Panel




Forensic Samples




Degraded DNA Testing

2U




4U

3U



5U

Degraded DNA: STR versus SNP Comparison





Identifiler[®] and SNPs: Empirical Data Genotypes on CEPH families

- **Three ‘CEPH’ families, European origin**
- **3-generation families with genotyped parents and grandparents, many children**
- **Genotyped for Identifiler[®] and 49 SNPs**
 - **Power of discrimination**
 - **Paternity exclusion**

STRs: Number of Loci Different Between Samples


Overall mean # loci different = 13.8




SNPs: Number of Genotype Differences Between Samples

Overall mean # loci different about 27



SNPs with different genotypes in 2 samples: full sibs in Red



SNPs with different genotypes between 2 samples:
Full Sibs in Red



Number of Genotype Differences Between Related and Unrelated Samples




Identifiler[®] and SNPforID Panels: Average Probability of Paternity Exclusion

	Identifiler [®]	49 SNPs	30 SNPs	20 SNPs
CEPH 20 samples	$8.986 * 10^{-7}$	$2.08 * 10^{-5}$	$8.22 * 10^{-4}$	$1.36 * 10^{-2}$

Paternity exclusion analysis: CERVUS 2.0

Identifiler[®] and SNPs


- Identifiler[®] has higher LOD score for paternity detection in 20 CEPH children



LOD score analysis: CERVUS 2.0

Identifiler[®] and SNPs

- SNPs can better distinguish between close relatives for paternity identification



LOD score analysis: CERVUS 2.0

Conclusions

- **A 49-plex single tube reaction was designed**
- **A gender-specific marker was included**
- **Able to reliably detect all SNPs on the 3130 and 3730**
- **Comparison with TaqMan[®] allelic discrimination assays indicated a concordance rate of 100%**
- **The method was used to successfully genotype 41 samples from three CEPH families**
- **The method worked with with several sample types (blood, buccal swabs, semen, etc.)**
- **Testing with DNA degraded using sonication and DNase I treatment indicated that the method works well with degraded DNA**

Conclusions: Identifiler[®] and SNPforID SNP Panel

- Overall informativeness is comparable
- SNPforID panel is better for identity
 - About 45 SNPs are comparable to Identifiler[®]
- Identifiler[®] is better for paternity exclusion
 - But SNPs give greater distinction between closely related candidates
- SNPs give larger absolute number of loci different

SNP Panels From Ken Kidd

G Model
FSIGSS-442; No of Pages 2

ARTICLE IN PRESS

Forensic Science International: Genetics Supplement Series xxx (2009) xxx–xxx



Contents lists available at ScienceDirect

Forensic Science International: Genetics Supplement Series

journal homepage: www.elsevier.com/locate/FSIGSS



Research article

Multiplexed SNP detection panels for human identification

Rixun Fang^a, Andrew J. Pakstis^b, Fiona Hyland^a, David Wang^a, Jaiprakash Shewale^a,
Judith R. Kidd^b, Kenneth K. Kidd^b, Manohar R. Furtado^{a,*}

^a Applied Markets, Applied Biosystems/Life Technologies, 850 Lincoln Centre Dr, m/s 402, Foster City, CA 94404, United States

^b Genetics Dept. Yale University School of Medicine, New Haven, CT 06520, United States

Fang, R. et al., *For Sci International Genetics Suppl. Series*
(2009) doi:10.1016/j.fsigss.2009.08.161

Technology
Transition Workshop 

Features

- A panel of 92 SNPs was assembled for human identification after testing more than 500 potential candidate SNPs from the list available since January 2009 at the site
[\(<http://info.med.yale.edu/genetics/kkidd/contents.html>\)](http://info.med.yale.edu/genetics/kkidd/contents.html)
- SNPs were selected based on testing samples from 44 populations across the globe using TaqMan[®] allelic discrimination formats in a collaboration
- High average heterozygosity (> 0.4)
- Low global F_{st} values (< 0.06) as a set of highly discriminative SNPs suitable for human identification in all geographic populations

Features


- **The SNP list includes 45 essentially unlinked SNPs distributed across the human autosomes**
 - **These SNPs are an excellent panel for individual identification providing high match probabilities in the range of values comparable to those for the 13 CODIS STR markers**
- **The unlinked status of these 45 SNPs also make them useful for resolving close biological relationships**
- **Additional SNPs on the list can be used to provide a higher power of discrimination, if required**
- **No meaningful departures from Hardy-Weinberg ratios were seen for any SNP in the populations tested**

Features

- **Separating out six SNPs exhibiting some linkage disequilibrium (LD) left 86 SNPs in the panel with no significant pair-wise LD**
- **Two panels of 49 and 41 SNPs were multiplexed for initial testing**
- **Coding SNPs and SNPs with known functional or phenotypic manifestations were excluded**
- **To test the performance characteristics of the two panels of SNPs and to compare their utility to STR analysis for human identification and paternity testing, 41 individuals from three different CEPH families spanning three generations were genotyped**

Informativeness: SNPs versus STRs

- **Probability of Identity: 42 SNPs = 15 STRs**



Data Analysis

- **41 samples**
- **80 markers gave good results**
- **Mean # genotype differences between any pair of samples = 44.76**
- **Standard deviation of # genotype differences between any pair of samples = 6.7**
- **Minimum # genotype differences between pair of samples (many full-sibs) = 18**
- **Maximum # genotype differences between any pair of samples = 64**

Identity

# of SNPs	Log ₁₀ Probability of Identity
30	-12.7774
35	-14.9069
40	-17.0367
45	-19.1661
50	-21.2955
55	-23.4252
60	-25.5550
65	-27.6845
70	-29.8138
75	-31.9433
80	-34.0728

Leverage CE Platform Flexibility for Fixed Sets of SNPs / Mutations

SNPforID
49-plex

KenKidd
SNPs ID
3 pools of 45

AIM markers
38-plex

Y-SNPs
42-plex

Mito SNPs
2-panels

*Custom kit on market.
Dtl for a new product*

Clock gene
mutation set

Bipolar
mutation
set

Ras mutation
set (12)

50-plex
ZipChutes®

MRSA
typing panel

Coccidiosis
typing panels

Multiplex design pipeline


EGFR mutation
set (28)

Microbial panel with detection of over
36 targets, 4 drug resistance mutations
4 virulence factors in one tube

**50 zip chutes with 2-dyes, expandable to 100 with
additional dyes**

Technology
Transition Workshop 

Platform Flexibility: Manufacturing



Population Indicative Panel Study: Chris Phillips – Feasibility Study


- **SNP targets:**
 - 32 SNPs
 - 30 bi-allelic
 - 2 tri-allelic
- **Tri-allelic OLA oligo design**
 - Two bins were used for each locus

SNP ID	SNP_NAME	GENOTYPE	ASO1	ASO2
rs5030240	NBS-1	A/C/G	M5.1_017_1_A40037=G	M5.1_017_2_A40118=C
			M5.1_005_1_A40033=G	M5.1_005_2_A30007=A
rs4540055	NBS-9	A/C/T	M5.1_032_1_A50055=T	M5.1_032_2_A50056=G
			M5.1_013_1_A40034=T	M5.1_013_2_A40027=A


Creating the Final Combination of SNPs

- **Select best loci for between-population F_{st}**
- **Optimize the chromosome spread**


32 SNP Panel Study For a Male's DNA



32 SNP Analysis For a Female's DNA



Tri-allelic Markers for DNA Mixture Study




Y-SNP Study: Potential Collaboration with Genographic Project

- **Genographic Project:**
 - A five-year effort to understand the human journey
 - Launched by The National Geographic Society, IBM, geneticist Spencer Wells, and the Waitt Family Foundation have launched the Genographic Project
 - M humanity's genetic journey through the ages
- **131 targets were selected to cover all Y haplogroups**
 - SNPs were selected based on information from Peter Underhill, YCC, and publications
 - Multiplex PCR primers and OLA oligos designed

Y-SNPs

- **Tree**

Y-SNPs – 29 Loci




Multiplexed SNP Detection System for Mitochondrial DNA


- **Rixun Fang PhD***
- **Pius Brzoska PhD***
- **Peter Ma PhD***
- **Peter Vallone PhD****
- **Manohar Furtado PhD***
 - ***Applied Markets and Molecular and Cell Biology Division, Applied Biosystems, Foster City CA 94404**
 - **** NIST Building 227, Room B24, 100 Bureau Drive Gaithersburg, MD 20899-8312 USA**
- **AAFS 2007**

GeneMapper[®] Software Cluster Plots for SNP Genotyping

MtDNA




Nuclear DNA





Caucasian Multiplex Panels for Distinguishing the Most Common HV1 / HV2 Types (Vallone et al., Coble et al.)

A	B	C	D	E	F	G	H
477 T/C	477	72	482	4808	64	3826	64
3010 G/A	3010	513	5198	5147	4745	3834	4688
4580 G/A	3915	4580	6260	9380	10211	4688	11377
4793 A/G	5004	5250	9548	9899	10394	6293	12795
5004 T/C	6776	11719	9635	11914	10685	7891	13293
7028 C/T	8592	12438	11485	15067	11377	11533	14305
7202 A/G	10394	12810	11914	16519	14470	12007	16519
10211 C/T	10754	14770	15355	14560	12795		
12858 C/T	11864	15833	15884	16390	15043		
14470 T/A	15340	15884	16368	14869	16390		
16519 T/C	16519	16519	16519				
H:1	H:2,H:3,H:6	V:1,H:5	J:1,J:2,K:2,K:3	J:4,T:2,T:3,H:4	V:1,H:1,H:2,H:3	J:1,J:3,T:1	K:1

SNP Detection for Panel A Coding SNPs




Mitochondrial Genotype Verification by *mitoSEQr*TM Sequencing System




Human Genomic DNA

PCR Reaction



Sequencing Reaction



Electrophoresis



Data Analysis

Mitochondrial SNP Detection by mitoSEQr™ System

Various_Target RC	Various_Target_exon	Various_Target_intron_1
Summary	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
NT Variants	- - - - -	- - - - -
Index	6441	6451
Reference	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NAO2254	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NAO5920	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NAO8587	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NAO9947	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA10924	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA10959	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA12593	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA14448	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA14529	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA14548	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ NA14672	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A
▶ CEPH	T A C G T T G T A G C T	C A C T T C C A C T A T G T C C T A T C A A

Correlation Between GenPlex™ and Sequencing Methods: 100%

	DNA	7202AG	5004TC	4793AG	3010GA	10211CT	14470TA	4580GA	477TC	12858CT	7028CT	16519TC
1	NA00131	A	T	A	A	C	T	G	T	C	C	C
2	NA00333	A	T	A	G	C	T	G	T	C	T	C
3	NA00546	A	T	A	G	C	T	G	T	C	C	T
4	NA00607	A	T	A	G	C	T	G	T	C	T	C
5	NA00893	A	T	A	G	C	T	G	T	C	C	T
6	NA00946	A	T	A	G	C	T	G	T	C	T	C
7	NA01805	A	T	A	A	C	T	G	T	C	C	C
8	NA01814	A	T	A	A	C	T	G	T	C	T	T
9	NA01953	A	T	A	G	C	T	G	T	C	C	T
10	NA01954	A	T	A	G	C	T	G	T	C	T	T
11	NA01990	A	T	G	G	C	T	G	T	C	C	C
12	NA02254	A	T	A	G	C	T	G	T	C	T	C
13	NA05920	A	T	A	G	C	T	G	T	C	C	T
14	NA08587	A	T	A	G	C	T	G	T	C	T	T
15	NA09947	A	T	A	G	C	T	G	T	C	C	C
16	NA10924	A	T	A	G	C	T	G	T	C	C	T
17	NA14448	A	T	A	G	C	T	G	T	C	C	T
18	NA14529	A	T	A	G	C	T	G	T	C	T	T
19	NA14548	A	T	A	G	C	T	G	T	C	T	C
20	NA14672	A	T	A	G	C	T	G	T	C	T	C


Correlation with SNaPshot[®] Study by NIST

- 44 U.S. Caucasian samples sharing a common HVI / HVII type
- 100% correlation


NIST calls	3010	4793	10211	5004	7028	7202	16519	12858	4580	477	14470
mt1	GT38069	G	G	C	T	C	A	C	G	T	T
mt10	MT97163	G	G	C	T	C	A	C	G	T	T
mt11	TT50697	G	A	C	T	C	A	C	G	T	T
mt12	TT50715	A	A	C	T	C	A	C	G	T	T
mt13	TT50720	G	A	C	T	C	A	C	G	T	T
mt14	TT51702	G	A	C	T	C	A	C	G	T	T
mt15	TT51703	G	A	C	T	T	A	C	A	T	T
mt16	UA16894	A	A	C	T	C	A	C	G	C	T
mt17	UT57291	G	A	C	T	A	T	C	A	T	T
mt18	UT57293	A	A	C	T	C	A	C	G	T	T
mt19	UT57301	G	G	C	T	C	A	C	G	T	T
mt2	GT38075	G	G	C	T	C	A	C	G	T	T
mt20	UT57310	G	A	C	T	C	A	C	G	T	T
mt21	UT58314	G	A	C	T	C	A	C	G	T	A
mt22	UT58333	G	A	C	C	C	A	T	C	G	T
mt23	UT58334	A	A	C	T	C	A	C	G	T	T
mt24	UT58335	G	A	C	T	C	A	C	G	T	T
mt25	WT51345	A	A	C	T	C	A	C	G	T	T
mt26	WT51359	A	A	C	T	C	A	C	G	T	T
mt27	WT51378	G	A	C	T	C	A	C	G	T	T
mt28	WT52457	A	A	C	T	C	A	C	G	T	T
mt29	WT52472	A	A	C	T	C	A	C	G	T	T
mt3	GT38078	G	A	C	T	C	A	C	G	T	T
mt30	WT52474	A	A	C	T	C	A	C	G	T	T
mt31	MT94861	G	A	C	T	C	A	C	G	T	T
mt32	GT38094	G	A	C	T	C	A	C	G	T	T
mt33	JT51178	G	A	C	T	C	A	C	G	T	T
mt34	MT94846	G	A	C	T	C	A	T	C	G	T
mt35	MT94847	G	A	C	T	C	A	C	G	T	T
mt36	MT94850	A	A	C	T	C	A	C	G	T	T
mt37	MT94892	G	A	C	T	C	A	T	C	G	T
mt38	MT97152	G	A	C	T	C	A	C	G	T	T
mt39	MT97180	A	A	C	T	C	A	C	G	T	T
mt4	GT38081	G	A	C	T	T	A	T	C	A	T
mt40	TT51669	A	A	C	T	C	A	C	T	G	C
mt41	GT38091	G	A	C	T	T	A	T	C	A	T
mt42	MT94812	G	A	C	T	T	A	T	C	A	T
mt43	OT07753	G	A	C	T	T	A	C	A	T	T
mt44	PT86478	A	A	C	T	C	A	C	C	G	T
mt5	GT38086	A	A	C	T	C	A	C	C	G	C
mt6	GT38100	A	A	C	T	C	A	C	C	G	T
mt7	GT38106	G	A	C	T	C	A	T	C	G	T
mt8	GT38108	A	A	C	T	C	A	C	C	G	T
mt9	MT94842	A	A	C	T	C	A	C	C	G	T

Sensitivity of 11-Plex SNP Detection System

1 pg
gDNA



0.1 pg
gDNA




Technology
Transition Workshop




System Validation Using Forensic Samples


Semen stain




Blood stain




Saliva stain




Milk stain




Urine stain



Vaginal swab



Fingerprint




Technology
Transition Workshop



Potential to Genotype 58 SNPs

A	B	C	D	E	F	G	H
477 T/C	477	72	482	4808	64	3826	64
3010 G/A	3010	513	5198	5147	4745	3834	4688
4580 G/A	3915	4580	6260	9380	10211	4688	11377
4793 A/G	5004	5250	9548	9899	10394	6293	12795
5004 T/C	6776	11719	9635	11914	10685	7891	13293
7028 C/T	8592	12438	11485	15067	11377	11533	14305
7202 A/G	10394	12810	11914	16519	14470	12007	16519
10211 C/T	10754	14770	15355	14560	12795		
12858 C/T	11864	15833	15884	16390	15043		
14470 T/A	15340	15884	16368	14869	16390		
16519 T/C	16519	16519	16519				
H:1	H:2,H:3,H:6	V:1,H:5	J:1,J:2,K:2,K:3	J:4,T:2,T:3,H:4	V:1,H:1,H:2,H:3	J:1,J:3,T:1	K:1

Potential of Scaling up to 58 SNPs Using Two Multiplex Designs



Conclusions

- **An easy and robust method has been developed for mitochondrial SNP genotyping**
- **The results from this approach are completely concordant with sequencing and SNaPshot[®]**
- **This method fulfills the following requirements for forensic studies:**
 1. **Low consumption of gDNA samples (picogram of genomic DNA)**
 2. **High throughput (up to 48 SNPs / reaction)**
 3. **High accuracy**
 4. **Fast and ease of use (< 7 hours)**
- **This method can be used for varieties of forensic samples**

References

- Peter M. Vallone, Rebecca S. Just, Michael D. Coble, John M. Butler, Thomas J. Parsons. “A multiplex allele-specific primer extension assay for forensically informative SNPs distributed throughout the mitochondrial genome” *Int J Legal Med* (2004) 118(3) 147–157
- Coble, M.D., Just, R.S., O'Callaghan, J.E., Letmanyi, I.H., Peterson, C.T., Irwin, J.A., Parsons, T.J. “Single nucleotide polymorphisms over the entire mtDNA genome that increase the power of forensic testing in Caucasians” *Int J Legal Med* (2004) 118(3) 137–146

Clinical Applications: Pathogen and Drug Resistance Detection

- **Application for detection of multiple pathogens, drug resistant mutations and virulence factors**
- **Equal or better than existing formats**

SIMULTANEOUS DETECTION OF 17 RESPIRATORY VIRUSES AND SUBTYPES WITH THE LUMINEX DIAGNOSTICS® RESPIRATORY VIRUS PANEL KIT

Christine Robinson, Kristin Murray, Kristi Lookner, and Qi Wei
The Children's Hospital, Denver CO

Respiratory Viral Panel




Virus	Type/Subtype	Gene Target
Influenza A		Matrix gene
	H1	Hemagglutinin gene
	H3	Hemagglutinin gene
	H5	Hemagglutinin gene
Influenza B		Matrix gene
Respiratory Syncytial Virus	Type A	Fusion gene
	Type B	Fusion gene
Coronavirus	229E	Nucleocapsid gene
	OC43	Nucleocapsid gene
	SARS	Nucleocapsid gene
	NL63	Nucleocapsid gene
	NHKU1	Nucleocapsid gene
Parainfluenza virus	Parainfluenza 1	HN gene
	Parainfluenza 2	HN gene
	Parainfluenza 3	HN gene
	Parainfluenza 4	HN gene
Metapneumovirus		Polymerase (L) gene
Enterovirus/ Rhinovirus		5'NCR
Adenovirus		Hexon gene

U.S. Naval Medical Research Unit 2

Total time about 6 to 7 hours


ID Tag™ RVP Overview




II. Setup & Multiplex PCR (3 hrs)

Viruses	Primer	Comment	Results
Influenza A	H1, H3	H5 not scored	3 (4)
Influenza B			1
RSV	A, B		2
Coronaviruses	NL63, HKU1, 229E, OC43	SARS available outside the US	4 (5)
Parainfluenza	1-4		4
HRP1		Dutch sequence	1
Enterovirus		No differentiation	1
Adenovirus			1
Control	M1C bacteriophage	Internal extraction control	1
Lambda DNA		Run control	1
Total Analytes: 19 (21)			


III: Amplicon Treatment & Multiplex TSPE (2.5 hrs)



IV: Bead Hybridization + Reporter Addition (1 hr)



V & VI: Detection and Analysis (30 minutes)



Technology Transition Workshop



Proof of Concept: GenPlex™ Assay Design for Multiple Pathogen Targets (41 Targets and Six Mutations)

Design	Assay Targets	Targets Tested
1	1	Humanparainfluenza
2	1	InfluenzaA(H1N1)KlA82D(-N)
3	1	EnterovirusA
4	1	SARScoronavirusF1Mut1
5	1	InfluenzaBvirus
6	1	Humanparainfluenza1
7	1	Humanparainfluenza2
8	1	Humanparainfluenza3
9	1	InfluenzaA(H1N1)G2G0(-N)
10	1	SARScoronavirusF1Mut1
11	1	Humanrotavirus4
12	1	InfluenzaA(H3N2)
13	1	Humanrespiratorysyncytialvirus
14	1	InfluenzaA(H1N1)G1G1(-N)
15	2	InfluenzaA(H1N1)G1G1(-N)
16	2	Brillusathais
17	2	BrillusspCBD1B
18	2	Francislatularemisisubpularis
19	2	Francislatularemisis
20	2	Humanparainfluenza3
21	2	Humanpapillomavirus16
22	2	HumanRhinovirus1
23	2	HumanRhinovirus2

Design	Assay Targets	Targets Tested
24	2	EBV
25	2	Humanhepatitis4
26	2	Humanhepatitis5
27	2	EKpolyomavirus
28	2	JCpolyomavirus
29	2	Sriantivirus40
30	2	Humanhepatitis7
31	2	HumanT-lymphotropicvirus2
32	2	HIV1
33	2	HepatitisBvirus
34	2	HepatitisBvirus
35	2	HepatitisCvirus
36	2	Humanimmunodeficiencyvirus1
37	2	Humanadenovirus6
38	2	Humanadenovirus6
39	2	Humanimmunodeficiencyvirus2
40	2	Humanimmunodeficiencyvirus2
41	2	OW_UL9_C82G
42	2	OW_UL9_C82G
43	2	OW_UL9_d1981-2
44	2	OW_UL9_d1981-2
45	2	OW_UL9_MEDV
46	2	OW_UL9_MEDV
47	2	Humanpapillomavirus


All in one tube

Technology Transition Workshop




Multiplexing Study for Pathogen Detection

Singleplex Assay



Multiplex Assay




Cytomegalovirus Drug Resistance Mutations

- Mutations tested: UL97, M460V, C592G, A594V, L595S and pol del981-2

CMV_UL97_C592G	C[T/G]G
CMV_UL97_M460V	C[A/G]T
CMV_UL97_A594V	G[C/T]G
CMV_UL97_L595S	T[T/C]G
CMV_UL97_del981-2	C[ATCGAC]-C

- OLA assays: two alleles for mutation



Companion Diagnostics

- **Application for detection of mutation in EGFR and KRAS pathway**

EGFR Pathway Mutations Influencing Cancer Therapy

- **Mutations in EGFR and Ras confer drug resistance**

Current Technologies Used for EGFR Mutation Detection In Lung Cancer

Table 1. Methods for detecting *EGFR* mutations in lung cancer specimens

Technique	Reference	Sensitivity (% mutant DNA)	Mutations identified	Comprehensive detection of deletions and insertions?
Direct sequencing	Multiple	25	Known and new	Yes
PCR-SSCP	(10)	10	Known and new	Yes
TaqMan PCR	(11, 12)	10	Known only	No
Loop-hybrid mobility shift assay	(13)	7.5	Known only	Yes
Cycleave PCR	(14)	5	Known only	Yes
PCR-RFLP and length analysis	(15)	5	Known only	Yes
MALDI-TOF MS – based genotyping	(16)	5	Known only	No
PNA-LNA PCR clamp	(17)	1	Known only	No
Scorpions ARMS	(18)	1	Known only	No
dHPLC	(19 – 21)	1	Known and new	Yes
Single-molecule sequencing	(22)	0.2	Known and new	Yes
Mutant-enriched PCR	(23)	0.2	Known only	No
SMAP	(8)	0.1	Known only	No

Abbreviations: SSCP, single-strand conformation polymorphism; PNA-LNA: peptide nucleic acid – locked nucleic acid; MALDI-TOF MS, matrix-assisted laser desorption/ionization time-of-flight mass spectrometry; ARMS, amplified refractory mutation system; dHPLC, denaturing high performance liquid chromatography.

Clin Cancer Res 2007;13(17) September 1, 2007

DxS Developed EGF29 and K-RAS Mutation Test Kits

- **DxS's EGF29 and K-RAS Mutation Test Kits detect key mutations from EGFR and K-RAS that correlate with responsiveness to tyrosine kinase inhibitors, such as gefitinib and erlotinib, with some mutations having a sensitizing effect and others being linked to resistance**
- **Total of (7 + 29) = 36 tubes for assay**
- **2 to 3 ng DNA input / tube**
- **Total about 70 ng**

GenPlex™ Assay to Detect 34 Mutations In One Reaction


UID	Gene	Exon	SNP/Indel Mutation
EGFR_A750_E758del	EGFR	exon19	A[AGCAACATCTCCGAAAGCCAACAAGGA/-]A
EGFR_A767_S768_TLAins	EGFR	exon20	A[-/CGCTGGCCA]G
EGFR_E746_A748del	EGFR	exon19	G[GAATTAAGA/-]G
EGFR_E746_A750del	EGFR	exon19	A[GAATTAAGAGAAGCAACA/-]T
EGFR_E746_E749del	EGFR	exon19	G[GAATTAAGAGAA/-]G
EGFR_E746_T751del	EGFR	exon19	G[GAATTAAGAGAAGCAACA/-]T
EGFR_E746del	EGFR	exon19	A[AGG/-]A
EGFR_G719A	EGFR	exon18	G[G/C]C
EGFR_G719C	EGFR	exon18	G[G/T]G
EGFR_G719S	EGFR	exon18	G[G/A]G
EGFR_K745_E749del	EGFR	exon19	C[AAGGAATTAAGAGAA/-]
EGFR_K745_L747del	EGFR	exon19	C[AAGGAATTA/-]A
EGFR_L747_A750del	EGFR	exon19	A[TTAAGAGAAGCA/-]A
EGFR_L747_E749del	EGFR	exon19	A[TTAAGAGAA/-]G
EGFR_L747_K754del	EGFR	exon19	A[TTAAGAGAAGCAACATCTCCGAAA/-]G
EGFR_L747_P753del	EGFR	exon19	A[TTAAGAGAAGCAACATCTCCG/-]A
EGFR_L747_S752del	EGFR	exon19	A[TTAAGAGAAGCAACATCT/-]C
EGFR_L747_T751del	EGFR	exon19	A[TTAAGAGAAGCAACA/-]T
EGFR_L858R	EGFR	exon21	C[T/G]G
EGFR_L861Q	EGFR	exon21	C[T/A]G
EGFR_M766_A767_Alns	EGFR	exon20	G[-/GCCATA]G
EGFR_N756_L760del	EGFR	exon19	C[CAACAAGGAAATCCT/-]C
EGFR_P753_I759del	EGFR	exon19	T[CCGAAAGCCAACAAGGAAATC/-]C
EGFR_R748_P753del	EGFR	exon19	T[TTAAGAGAAGCAACATCTC/-]C
EGFR_S752_I759del	EGFR	exon19	T[ATCTCCGAAAGCCAACAAGGAAAT/-]C
EGFR_S768I	EGFR	exon20	A[G/T]C
EGFR_T751_E758del	EGFR	exon19	C[AACATCTCCGAAAGCCAACAAGGA/-]A
EGFR_T790M	EGFR	exon20	A[C/T]G
EGFR_V769_D770_ASVins	EGFR	exon20	G[-/GCCAGCGTG]G
KRAS_Ala12Val	K-RAS	exon2	G[C/T]T
KRAS_Arg12Cys	K-RAS	exon2	T[C/T]G
KRAS_Gly12Asp	K-RAS	exon2	G[G/A]T
KRAS_Gly12Ser	K-RAS	exon2	T[G/A]G
KRAS_Gly13Asp	K-RAS	exon2	G[G/A]C

DNA Input = 0.5 ng

Technology
Transition Workshop




Cancer Cell Line DNA SNP and INDEL Detection Using GenPlex® Approach (NCI-60 Panel)




KRas Gly13Asp


MDA-MB-435:
M14 melanoma
cells




MDA-MB-231:
Breast
adenocarcinoma



HCT-15:
colon
carcinoma



HCT-116
cells:colorect
al carcinoma




SK-MEL 28:
Human
melanoma cell

SK-MEL 2:
Human
melanoma cell


SK-MEL 5:
Human
melanoma cell

EGFR T751-E758 del




100% Concordance of Mutation Detection Between DNA Sequencing and GenPlex™ Platform

KRas _Ala12Val



KRas Gly12Ser





Minor Component Detection Sensitivity

Cell line DNA: ACHN/UACC1

Competitive OLA format to increase sensitivity

0:1
1:1
1:10
1:20
1:100



Technology Transition Workshop

Conclusions

- **The GenPlex™ system provides a solution for multiple mutation and nucleic acid detection:**
 - 98 nucleic acid targets
 - 49 SNP targets
 - Any combination of both per reaction
- **GenPlex™ system has high sensitivity: suitable for low amount genomic DNA (0.1 ng to 1 ng of NA)**
- **High specificity is ensured at two levels; both at the PCR and OLA reactions**
 - > 99.5% accuracy reported in customer studies
- **The system is amenable to automation**
- **Streamlined protocol with time to result time of about 4 to 5 hours**
- **Highly Integrated GenPlex™ System will be an excellent platform for molecular diagnostic applications**

Questions?

Contact Information

Manohar Furtado Ph.D.
Applied Biosystems/Life Technologies
(650) 554 2670
Manohar.furtado@lifetech.com

© 2009 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. All images are courtesy of Life Technologies unless otherwise indicated.