


GeneMapper® *ID-X* Software

Version 1.0/1.1

Note: To improve the clarity of graphics in this PDF file, use the zoom tool to increase magnification to 150% or greater.

GeneMapper®
ID-X



GeneMapper® *ID-X* Software Version 1.0/1.1

Getting Started

1

Installation
Requirements

2

Performing GMID-
X v1.0
Pre-Installation
and GMID-X v1.1
Upgrade
Procedures

3

Performing a Full
Software
Installation

4

Performing a
Client Software
Installation

5

Performing Post-
Installation
Procedures

6

Setting Up
Autoanalysis

7

© Copyright 2008, Applied Biosystems. All rights reserved.

For Research Use Only. Not for use in diagnostic procedures.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

APPLIED BIOSYSTEMS DISCLAIMS ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. IN NO EVENT SHALL APPLIED BIOSYSTEMS BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Notice to Purchaser: License Disclaimer.

Purchase of this software product alone does not imply any license under any process, instrument or other apparatus, system, composition, reagent or kit rights under patent claims owned or otherwise controlled by Applera Corporation, either expressly, or by estoppel.

Applera, Applied Biosystems, AB (Design), ABI PRISM, AmpFLSTR, COfiler, GeneMapper, Profiler Plus, SGM Plus and Yfiler are registered trademarks and MiniFiler and SEfiler are trademarks of Applera Corporation or its subsidiaries in the US and/or certain other countries.

This product includes software developed by the Apache Software Foundation.

This product includes software developed by the ExoLab Project.

JNIRegistry is copyrighted © by ICE Engineering, Inc.

Microsoft and Windows are registered trademarks of Microsoft Corporation.

Oracle is a registered trademark of Oracle Corporation.

All other trademarks are the sole property of their respective owners.

Part Number 4375953 Rev. B

07/2008

Contents

| | | |
|------------------|---|------------|
| | Preface | vii |
| | How to Use This Guide | vii |
| | How to Obtain More Information | x |
| | How to Obtain Support | xi |
| Chapter 1 | Getting Started | 1 |
| | Start Here: Getting the Most Out of this Guide | 2 |
| | Installation Options | 3 |
| Chapter 2 | Installation Requirements | 7 |
| | Installation CDs | 8 |
| | Full Installation Computer Requirements | 9 |
| | Client Installation Computer Requirements | 11 |
| | Verifying Your Computer Setup | 13 |
| | Analyzers and Compatible Software | 14 |
| Chapter 3 | Performing GMID-X v1.0 Pre-Installation and GMID-X v1.1 Upgrade Procedures | 15 |
| | Overview | 16 |
| | Requirements for Upgrading GMID v3.1 to GMID v3.2 | 18 |
| | Exporting GMID v3.2.x Data | 18 |
| | Using the Export Utility (GMID v3.2.x Software Only) | 19 |
| | Manually Exporting GMID v3.2 Data Objects | 23 |
| | Uninstalling the GMID v3.x Software | 24 |
| | Enabling Autoanalysis | 26 |
| | Starting the Data Collection Software on the 31xx Data Collection Computer | 26 |

| | | |
|------------------|--|-----------|
| | Obtaining the Data Collection Computer Name | 29 |
| | Obtaining the GeneMapper® <i>ID-X</i> Software Host Computer Name | 30 |
| | Requirements for Upgrading GMID-X v1.0/1.0.1 to GMID-X v1.1 | 31 |
| Chapter 4 | Performing a Full Software Installation | 33 |
| | Overview | 34 |
| | User Account Requirements for Full Installations | 34 |
| | Performing a Full Installation | 36 |
| | Full Installation Requirements | 36 |
| | Installing the Full GeneMapper® <i>ID-X</i> Software | 37 |
| Chapter 5 | Performing a Client Software Installation | 43 |
| | Overview | 44 |
| | User Account Requirements for Client Installations | 44 |
| | Performing a Client Installation | 46 |
| | Client Installation Requirements | 46 |
| | Installing the Client GeneMapper® <i>ID-X</i> Software | 47 |
| | Resetting Permissions for the Installation | 52 |
| Chapter 6 | Performing Post-Installation Procedures | 53 |
| | Registering the GeneMapper® <i>ID-X</i> Software | 54 |
| | Logging Into GMID-X After Installation | 55 |
| | Connecting to a New Host | 56 |
| | Logging Out of the GeneMapper® <i>ID-X</i> Software | 58 |
| | Additional Post-Installation Tasks | 59 |
| Chapter 7 | Setting Up Autoanalysis | 61 |
| | Overview of Autoanalysis Process | 62 |
| | Determining if Your Computer Can Perform Autoanalysis | 64 |
| | Autoanalysis Instrument and Software Configurations | 65 |
| | Overview of Autoanalysis Setup | 66 |
| | Local and Remote Autoanalysis | 66 |

| | | |
|-------------------|--|------------|
| | Data Processing Options | 67 |
| | Setting Up Autoanalysis | 68 |
| | Procedure 1: Creating a User Account for Autoanalysis | 69 |
| | Procedure 2: Creating an Instrument Protocol for Autoanalysis | 70 |
| | Procedure 3: Creating a Results Group for Autoanalysis | 72 |
| | Procedure 4: Creating a Plate Record for Autoanalysis | 80 |
| | Recommended Data Processing Tasks | 83 |
| | Optimizing Data Sharing | 83 |
| | Synchronizing GMID-X Software with the Data Collection Software | 90 |
| Appendix A | Installing the Generic Updater Software | 91 |
| | Installing the Generic Updater Software | 92 |
| | Removing the Generic Updater Software | 94 |
| Appendix B | Troubleshooting the Installation | 97 |
| | Troubleshooting Tasks | 98 |
| | Disabling Firewall Protections | 98 |
| | Canceling an Installation | 99 |
| | General Installation Solutions | 99 |
| | Troubleshooting Tools | 105 |
| | Checklist | 105 |
| | Interdepartmental Forms | 107 |
| Appendix C | Glossary of Terms | 111 |
| | Definition of Terms | 112 |
| | Index | 115 |

Preface

How to Use This Guide

| | |
|------------------------------|--|
| Purpose of This Guide | <p>This guide describes the procedures for installing the GeneMapper® <i>ID-X</i> Software Version 1.0 or Version 1.1 (GMID-X v1.0 or v1.1).</p> <ul style="list-style-type: none">• Chapter 1, Getting Started – Provides a roadmap to help you determine the best approach for installing the GMID-X v1.0 or v1.1 on the target computer.• Chapter 2, Installation Requirements – Contains the specifications for the computers that are targeted for a new GMID-X v1.0 or v1.1 software installation.• Chapter 3, Performing GMID-X v1.0 Pre-Installation and GMID-X v1.1 Upgrade Procedures – Provides procedures on backing up and exporting application data and objects, and on uninstalling GeneMapper® <i>ID</i> Software Version 3.x (GMID v3.x).• Chapter 4, Performing a Full Software Installation – Provides procedures for installing the full software on the Data Collection or non-Data Collection computer, to serve as the host.• Chapter 5, Performing a Client Software Installation – Provides procedures for installing the client software on the Data Collection or non-Data Collection computer.• Chapter 6, Performing Post-Installation Procedures – Describes the basic functions of registering the software, logging into and out of the computer, and general information following an installation.• Chapter 7, Setting Up Autoanalysis – Describes how to set up your system to perform automatic remote or local analysis of .fsa sample files generated on the Data Collection computer. |
|------------------------------|--|

- **Appendix A, Installing the Generic Updater Software** – Provides the procedure for installing the Generic Updater Software.
- **Appendix B, Troubleshooting the Installation** – Provides possible avenues for troubleshooting and resolving software installation problems.
- **Appendix C, Glossary of Terms** – Defines the terms commonly used throughout the installation guide.

Note: For instructions on configuring the administrative features of the software, see the *GeneMapper® ID-X Software v1.0 Administrator's Guide*. For instructions on configuring the software to prepare for analysis, see the *GeneMapper® ID-X Software v1.0 Getting Started Guide* and the *GeneMapper® ID-X Software Online Help*.

Audience This guide is written for laboratory personnel responsible for installing the GeneMapper *ID-X* Software v1.0 or v1.1.

Assumptions This guide assumes that you have a working knowledge of the Microsoft® Windows® operating system.

Text Conventions This guide uses the following conventions:

- **Bold** indicates user action. For example:
Enter **0**, then press **Enter** for each of the remaining fields.
- *Italic* text indicates new or important words and is also used for emphasis. For example:
Before analyzing, *always* prepare fresh matrix.
- A right arrow bracket (►) separates successive commands you select from a drop-down or shortcut menu. For example:
Select **File ► Open ► Spot Set**.
Right-click the sample row, then select **View Filter ► View All Runs**.

User Attention Words

Two user attention words appear in Applied Biosystems user documentation. Each word implies a particular level of observation or action as described below.

Note: Provides information that may be of interest or help but is not critical to the use of the product.

IMPORTANT! Provides information that is necessary for proper instrument operation, accurate chemistry kit use, or safe use of a chemical.

Examples of the user attention words appear below:

Note: The size of the column affects the run time.

Note: The Calibrate function is also available in the Control Console.

IMPORTANT! Make certain that you uninstall the GMID v3.2 software before you install the GMID-X v1.0 or v1.1 software on the non-Data Collection computer.

IMPORTANT! You must create a separate Sample Entry Spreadsheet for each 96-well plate.

How to Obtain More Information

Related Documentation

The following related documents are shipped with the software:

- ***GeneMapper® ID-X Software Version 1.0 Administrator's Guide*** – Describes how to configure the administrative features of the software, maintain the GMID-X and database dashboard software, and use the command-line interface.
- ***GeneMapper® ID-X Software Version 1.0 Getting Started Guide*** – Explains how to set up the software, set up a project, and analyze the example data provided with the GeneMapper ID-X Software.
- ***GeneMapper® ID-X Software Online Help*** – Describes all features and functions of the software and provides step-by-step procedures for using the software.
- ***GeneMapper® ID-X Software Version 1.0 Reference Guide*** – Describes process quality values (PQVs) and algorithms.
- ***GeneMapper® ID-X Software Version 1.0 Quick Reference Guide*** – Provides an easy-to-follow workflow for using the software with instrument compatibility information.
- ***Dashboard Online Help*** – Provides information about the database dashboard software.

Note: For additional information and support, see “How to Obtain Support” on page xi.

Obtaining Adobe Acrobat Reader

Portable document format (PDF) versions of this guide and all supporting documentation are available on the documentation CD available with your software.

To open the user documentation on the CD drive, you need the Adobe® Acrobat® Reader® software, which is available from www.adobe.com or www.appliedbiosystems.com/support/software.

Send Us Your Comments

Applied Biosystems welcomes your comments and suggestions for improving its user documents. You can e-mail your comments to:

techpubs@appliedbiosystems.com

How to Obtain Support

For HID support send an e-mail to:

HIDTechSupport@appliedbiosystems.com, or call (**from within North America only**) 888.821.4443 and select option **1**. For support **outside North America** and **internationally**, please contact your local support office or visit the Applied Biosystems support page to obtain contact information.


For the latest services and support information for all locations, go to **<http://www.appliedbiosystems.com>**, then click **Support**. At the Support page, you can:

- Access worldwide telephone and fax numbers to contact Applied Biosystems Technical Support and Sales facilities.
- Search through frequently asked questions (FAQs).
- Submit a question directly to Technical Support.
- Order Applied Biosystems user documents, MSDSs, certificates of analysis, and other related documents.
- Download PDF documents.
- Obtain information about customer training.
- Download software updates and patches.

1

Getting Started

1



This chapter covers:

- Start Here: Getting the Most Out of this Guide 2
- Installation Options 3

Start Here: Getting the Most Out of this Guide

Installation Guide Overview

This chapter presents an overview of the installation options available for GeneMapper® ID-X Software Version 1.0 or Version 1.1 (GMID-X v1.0 or v1.1). To benefit from this guide, you need to be familiar with your current computer and GeneMapper® ID Software v3.1 or v3.2 setup (if applicable).

To get started:

1. Read through [Chapter 1](#) to become familiar with your installation options.
2. See [Chapter 2](#) for general installation requirements.
3. See [Chapter 3](#) for pre-installation instructions.
4. See [Chapter 4](#) for full installation instructions.
5. See [Chapter 5](#) for client installation instructions.
6. See [Chapter 6](#) for post-installation instructions.
7. See [Chapter 7](#) for autoanalysis setup instructions.
8. See [Appendix A, Installing the Generic Updater Software](#) for generic updater installation instructions.
9. See [Appendix B, Troubleshooting the Installation](#) for troubleshooting instructions.
10. See [Appendix C, Glossary of Terms](#), for a definition of terms used frequently in this manual.

After you have read the applicable chapters in this guide and are installing your software, refer to the *GeneMapper® ID-X Software Administrator's Guide* for information on configuring the GeneMapper® ID-X Software administrative features and electronic data chain of custody systems.

Installation Options

Installation options for the GeneMapper® *ID-X* Software are based on the type of computer you are on and your database-sharing requirements.

You can either install a "full" version of the software, which consists of both the GMID-X software and the database, or you can install a "client" version of the software, which consists only of the GMID-X software. The full version on a computer can function as the database host for the client. Multiple clients can connect to the same database through a network connection for easy data sharing.

If you choose not to install the full GMID-X v1.0 or v1.1 software (as a co-installation) or the client GMID-X v1.0 or v1.1 software on your 31xx or 3730 Data Collection computer, you can leave the existing GeneMapper® *ID* Software Version 3.x (GMID v3.x) installed, or you can uninstall GMID v3.x and install the Generic Updater to generate .fsa files.

The following table provides information on the different types of installations and where to locate the corresponding information in this guide.


Table 1 Installation Options Summary

| Computer Type | Full Installation Summary | Client Installation Summary [‡] | Generic Updater |
|--|---|---|--|
| Non-Data Collection (Non-DC) Computer or 310 Data Collection (DC) Computer | <ul style="list-style-type: none"> Consists of Full Installation (with database) Non-DC computer can be configured for Remote Autoanalysis with 31xx DC computers Requires exporting of GMID v3.x data (page 18) Requires uninstallation of GMID v3.x (page 24) Review Full Installation Computer Requirements on page 9 prior to installing Installation time: approx. 1 hour <p>Refer to page 37 for the full installation procedure.</p> | <ul style="list-style-type: none"> Consists of Client Installation (no database) Non-DC computer can be configured for Remote Autoanalysis with 31xx DC computers Requires exporting of GMID v3.x data (page 18) Requires uninstallation of GMID v3.x (page 24) Review Client Installation Computer Requirements on page 11 prior to installing Installation time: 5 to 10 minutes <p>Refer to page 46 for the client installation procedure.</p> | N/A |
| 31xx or 3730 Data Collection (DC) Computer | <ul style="list-style-type: none"> Consists of Co-Installation (31xx) (DC software shares database) Can be set up for Local Autoanalysis Requires exporting of GMID v3.x data (page 18) No uninstallation of GMID v3.2 required[§] Installation time: approx. 1 hour <p>Refer to page 37 for the co-installation procedure.</p> | <ul style="list-style-type: none"> Consists of Client Installation (no database) Can be set up for Local Autoanalysis Requires exporting of GMID v3.x data (page 18) Requires uninstallation of GMID v3.x (page 24) Installation time: approx. 5 to 10 minutes <p>Refer to page 46 for the client installation procedure.</p> | <ul style="list-style-type: none"> Used on 31xx/3730 DC computers only, to generate .fsa files without installing GMID-X software Requires exporting of GMID v3.x data (page 18) Requires uninstallation of GMID v3.x (page 24) Installation time: less than 1 minute <p>Refer to page 91 for the generic updater procedure.</p> |

[‡] In order to install a client version of the software, you must first perform at least one successful full installation. The client installation requires you to enter the Computer Name of the database host computer during installation.

[§] You must recreate all the files associated with autoanalysis (e.g. Results Groups) since old files created as part of the GMID v3.2 integration cannot be reused. See [Chapter 7, “Setting Up Autoanalysis.”](#)

See [Figure 1](#) to determine the steps to take based on the desired installation configuration.



‡ Must upgrade to the GeneMapper® ID Software v3.x

Figure 1 Possible GeneMapper® ID-X Software installation configurations



Chapter 1 Getting Started

Installation Options


2

Installation Requirements

2

This chapter covers:

- Installation CDs 8
- Full Installation Computer Requirements 9
- Client Installation Computer Requirements 11
- Verifying Your Computer Setup 13
- Analyzers and Compatible Software 14



Installation CDs

Three different installation CDs are available:

- **New Full Installation** – The two *GeneMapper® ID-X Software Version 1.0* or *Version 1.1 Full Install* CDs (Disks 1 and 2) install the full version (GeneMapper ID-X Software and multi-user database) on the target computer. The full version CDs also include extra Generic Updater software and an installation Export Utility. You can install the software on any compatible computer.
or
- **Client Installation** – The *GeneMapper® ID-X Software Version 1.0* or *Version 1.1 Client Install* CD installs the GeneMapper ID-X Software Client on the target computer. You can install the client on any compatible computer (including one running Data Collection Software) and connect to any multi-user database on the same local area network.
- **Upgrade Installation (v1.1)** – The *GeneMapper® ID-X Software Version 1.1 Upgrade Install* CD installs the GeneMapper ID-X Software v1.1 mixture analysis tool on the target computer (client or full). You can install the software on any target computer already running GeneMapper ID-X Software v1.0 or v1.0.1. However, if you upgrade on a client computer you must also upgrade the computer that hosts the full multi-user database on the same local area network.

An additional CD comes with the software. This CD includes all installation and associated GMID-X documentation.

Note: The *GeneMapper® ID-X Software Version 1.0/1.1 Installation Guide* and supporting documentation are provided on a separate CD. To open the user documentation on the CD, you need Adobe® Acrobat® Reader® software. Download from www.adobe.com or www.appliedbiosystems.com/support/software.

Full Installation Computer Requirements

Minimum Configurations

Table 2 shows the recommended and minimum computer configurations for a *new full installation* (the GeneMapper® ID-X Software v1.0 or v1.1 and database).

Note: To install the GMID-X software, you need a local user account with administrative privileges. You can run the computer on regional settings but you need an English Operating System.

Table 2 Computer requirements for a new full installation

| Component | Recommended Configuration | Minimum Configuration‡ |
|--|--|--|
| IMPORTANT! The GeneMapper® ID-X Software will not function properly on computers with multiple physical processors. The GeneMapper® ID-X Software requires a single-processor or dual-core processor computer. It is not advised to install the software on a computer with other Oracle® or SQLPlus applications and clients, other than the database that is installed with the Data Collection software. If you must run other applications, install GMID-X first. | | |
| Computer | <ul style="list-style-type: none"> Intel Pentium® IV processor, >2.8 GHz 1 GB of RAM Two 120-GB hard drives§ Free disk space: <ul style="list-style-type: none"> 200 MB on the boot drive (drive on which the operating system is installed) 7 GB on the drive on which the GeneMapper ID-X Software is installed 20/48X IDE CD-ROM 10/100 NIC with RWU (internal) | <ul style="list-style-type: none"> Intel Pentium® processor, 733 MHz 512 MB of RAM Free disk space: <ul style="list-style-type: none"> 200 MB on the boot drive (drive on which the operating system is installed) 7 GB on the drive on which the GeneMapper ID-X Software is installed 20/48X IDE CD-ROM 10/100 NIC with RWU (internal) |
| Monitor | <ul style="list-style-type: none"> 1024 × 768 pixel resolution 19-inch color monitor (enables greater viewing of data) | <ul style="list-style-type: none"> 1024 × 768 pixel resolution 17-inch color monitor |

Table 2 Computer requirements for a new full installation (*continued*)

| Component | Recommended Configuration | Minimum Configuration [‡] |
|---------------------|--|------------------------------------|
| Operating System | Either of the following: <ul style="list-style-type: none"> • Microsoft® Windows® 2000 Professional Operating System, Service Pack 4 Update • Microsoft® Windows® XP Professional Operating System, Service Pack 2 or later Note: The GMID-X software is not supported on the Microsoft® Windows® 2003 Operating System and has not been tested on the Windows® Vista® Operating System. Additionally, operating system requirements for co-installation (installation on a Data Collection computer) may differ from those listed above. For information, see “Analyzers and Compatible Software” on page 14 . | |
| Ethernet Capability | <ul style="list-style-type: none"> • Network card for database installation • TCP/IP must be installed before database installation | |

[‡] The minimum configuration may not provide optimal performance.

[§] *Required* if the Data Collection Software and the GeneMapper® ID-X Software are installed on the same computer.

Note: Intel® Core™ 2 Duo processors ship as standard with many computers. The GMID-X software has been validated by Applied Biosystems to run on Intel® Core™ 2 Duo processor computers.

Required Time for Full Installation

The amount of time required to perform a full installation is:

- Approximately 1 hour on a non-data collection computer, depending on the speed of the computer.
- Approximately 25 minutes on a Data Collection computer (340 SCSI).
- Approximately 1 hour 15 minutes on other Data Collection computers.

Client Installation Computer Requirements

Minimum Configurations

Table 3 shows the minimum computer configurations for a *client installation* (GeneMapper® ID-X Software v1.0 or v1.1 only).

Note: To log in to the GMID-X software, you need a local user account with administrative privileges. For sites outside the U.S., you may need to change regional settings to U.S. settings following software installation.

Table 3 Computer requirements for client installation

| Component | Minimum Configuration |
|------------------|--|
| Client Computer | <ul style="list-style-type: none"> Intel Pentium® processor, 733 MHz 512 MB of RAM† 20/48X IDE CD-ROM drive 10/100 NIC with RWU (internal) 250 MB of free disk space 200 MB on the boot drive (drive on which the operating system is installed) |
| Monitor | <ul style="list-style-type: none"> 1024 × 768 pixel resolution 17-inch color monitor (smaller monitor lessens ability to view a full screen of data) |
| Operating System | <p>Either of the following:</p> <ul style="list-style-type: none"> Microsoft® Windows® 2000 Professional Operating System, Service Pack 4 Microsoft® Windows® XP Professional Operating System, Service Pack 2 or later <p>Note: The GMID-X software is not supported on the Microsoft® Windows® 2003 Operating System and has not been tested on the Windows® Vista® Operating System.</p> |

† Although you can install the GeneMapper ID-X Software on a computer with 512 MB of RAM, Applied Biosystems recommends using 1 GB or more of RAM for better performance.


**Required Time for
Client Installation**


The amount of time required to perform a client installation is approximately 5 to 10 minutes.

**Optimizing
Computer
Performance**

For recommendations on the number of client connections to use per host to optimize computer performance, contact Applied Biosystems.

Verifying Your Computer Setup





1. On the desktop, right-click  **My Computer**, then select **Properties**.
2. Click the General tab of the System Properties dialog box to verify that your computer meets the minimum requirements for installation.



Analyzers and Compatible Software

Instrumentation Table 4 lists Applied Biosystems Genetic Analyzer instruments and the corresponding Data Collection and Operating system software that are supported by the GeneMapper® ID-X Software

Table 4 Genetic Analyzers and Compatible Software

| Genetic Analysis Instrument | Data Collection Software and Operating System |
|---|--|
|  <p>ABI PRISM® 310 Genetic Analyzer</p> | <ul style="list-style-type: none"> • 310 Data Collection v3.0 • Windows 2000, SP 3 and 4 <i>and</i> • 310 Data Collection v3.1 • Windows XP, SP 2 or later |
|  <p>Applied Biosystems 3130/3130xl Genetic Analyzer</p> | <ul style="list-style-type: none"> • 3130 Data Collection v3.0/ 3130xl Data Collection v3.0 • Windows XP, SP 2 or later |
|  <p>ABI PRISM® 3100/3100-Avant™ Genetic Analyzer</p> | <ul style="list-style-type: none"> • 3100 Data Collection v2.0/ 3100-Avant™ Data Collection v2.0 • Windows 2000, SP 3 or later |
|  <p>Applied Biosystems 3730 Genetic Analyzer[‡]</p> | <ul style="list-style-type: none"> • 3730 Data Collection v3.0 • Windows XP, SP 2 or later |

[‡] The 3730 Genetic Analyzer does not include the 3730xl 96-capillary configuration. The 48-capillary configuration is validated for analysis of Identifiler® data generated from single-source data samples only.


3

Performing GMID-X v1.0 Pre-Installation and GMID-X v1.1 Upgrade Procedures

This chapter covers:

- Overview 16
- Requirements for Upgrading GMID v3.1 to GMID v3.2 18
- Exporting GMID v3.2.x Data 18
- Uninstalling the GMID v3.x Software 24
- Enabling Autoanalysis. 26
- Obtaining the GeneMapper® ID-X Software Host
Computer Name 30
- Requirements for Upgrading GMID-X v1.0/1.0.1 to GMID-X
v1.1 31

3



Overview

Reminders This chapter provides procedures that you must perform to install GeneMapper® *ID-X* Software v1.0 or v1.1 software on computers that *have* GeneMapper® *ID* Software v3.x currently installed. It also provides the procedure for upgrading from GeneMapper® *ID-X* Software v1.0 or v1.0.1 to v1.1.

IMPORTANT! If you are installing GMID-X v1.0 or v1.1 software on a computer *without* GMID v3.x, proceed to [Chapter 4 on page 33](#) for full installation instructions or [Chapter 5 on page 43](#) for client installation instructions.

IMPORTANT! If you are upgrading from GMID-X v1.0/1.0.1 to v1.1, go to [“Requirements for Upgrading GMID-X v1.0/1.0.1 to GMID-X v1.1” on page 31](#) of this chapter.

Before installing GMID-X software on computers that have GMID v3.x software installed, you must export the data objects (projects, analysis methods, etc.) as a means of backing up the data, then uninstall the existing software. Because data objects are deleted when GMID 3.x is uninstalled, if they are not exported and saved to a different location, they cannot be restored.

The *GeneMapper® ID-X Software Version 1.0* (or *Version 1.1*) *Full Install* CD contains a GMID v3.2.x Software Export Utility, which automatically exports all data objects stored within your GMID v3.2.x database and stores the files in a user-defined location. This utility is only designed to work in conjunction with GMID v3.2 and above. Therefore, if you are still running GMID v3.1 software, upgrade to GMID v3.2 to take advantage of this export utility (see [“Requirements for Upgrading GMID v3.1 to GMID v3.2” on page 18](#)).

Upgrade to GMID v3.2 if you also wish to import the exported files into GMID-X v1.0 or v1.1 once the software is installed, and if you plan to install the full version of GMID-X v1.0 or v1.1 on a 3lxx Data Collection computer.

Note: GMID v3.2.x data objects may be imported and viewed in GMID-X v1.0 or v1.1; however, their use is limited. (See [“Optimizing Data Sharing” on page 83](#)).

Requirements for Upgrading GMID v3.1 to GMID v3.2

To upgrade from GMID v3.1 to GMID v3.2, you must have a *GeneMapper® ID Software v3.2 Installation* CD and a valid GMID v3.1 registration code. If you do not have one or both of these items, please contact your local technical support organization.

Exporting GMID v3.2.x Data

To export data, you can use the Export Utility to automatically export all data objects from the GeneMapper Manager and Panel Manager within the GeneMapper® ID Software v3.2.x application. Alternatively, you can export each data object manually.


Applied Biosystems recommends that you use the Export Utility to more efficiently export data. Once the data is exported and saved to a new location, archive these objects using your lab's standard procedure.

Using the Export Utility (GMID v3.2.x Software Only)


Exporting GMID v3.2.x Data with Utility Tool

1. Insert the *GeneMapper® ID-X Software Version 1.0* (or *Version 1.1*) *Full Install CD* into the CD drive.
2. On the splash screen, click **Export Data from GeneMapper ID v3.2 or 3.2.1**.

Note: If there is no GMID v3.2 on the target system, the export button is disabled.




3. When the Export tool window is displayed, click **Browse** to specify the location to save all exported files, or enter a destination folder in the field provided:




IMPORTANT! Choose an export location outside the **Applied Biosystems ▶ GeneMapper** application folder. The GeneMapper folder is deleted when GMID v3.x is uninstalled.

4. Click **Export**.
 - If the utility detects a conflict with the location specified, an error message is displayed. Follow the instructions given in the error messages before proceeding. These error messages include:
 - No destination is specified for export files. Please select a destination.
 - Specified location does not exist. Please select a new location.
 - Access to specified location denied. Please check folder permissions or select a new location.
 - You have chosen to export data to the **Applied Biosystems ▶ GeneMapper** application folder. This folder is deleted when GMID 3.2 is uninstalled. Select a location outside of the **Applied Biosystems ▶ GeneMapper** application folder.

- If the export is successful the following message is displayed:



- If the export is only partially successful, the following message is displayed:



Proceed to the log file and manually export all data objects that failed to export.

After you have exported all your files, review the following table to determine the steps to take next.

Table 5 Installation Options following Auto Export

| Computer Type | Full Installation | Client Installation |
|----------------------------------|--|--|
| Non-DC or 310 DC Computer | <ol style="list-style-type: none"> 1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24. 2. Install Full GMID-X v1.0 or v1.1. See “Performing a Full Installation” on page 36. | <ol style="list-style-type: none"> 1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24. 2. Install Client GMID-X v1.0 or v1.1. See “Performing a Client Installation” on page 46. |

| Computer Type | Full Installation | Client Installation |
|---------------------------------|--|---|
| 31xx or 3730 DC Computer | <p><i>Do not uninstall.</i></p> <p>Install Full GMID-X v1.0 or v1.1. See “Performing a Full Installation” on page 36.</p> | <ol style="list-style-type: none">1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24.2. Install Client GMID-X v1.0 or v1.1. See “Performing a Client Installation” on page 46. |

Manually Exporting GMID v3.2 Data Objects

For information on manually exporting files out of GMID v3.x, refer to the section on exporting items from GMID Manager in the *GeneMapper® ID Software Version 3.1 User's Guide*.

After you have exported all your files, review the following table to determine the steps to take next.

Table 6 Installation options following manual export

| Computer Type | Full Installation | Client Installation |
|----------------------------------|--|--|
| Non-DC or 310 DC Computer | <ol style="list-style-type: none"> 1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24. 2. Install Full GMID-X v1.0 or v1.1. See “Performing a Full Installation” on page 36. | <ol style="list-style-type: none"> 1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24. 2. Install Client GMID-X v1.0 or v1.1. See “Performing a Client Installation” on page 46. |
| 31xx DC Computer | <p><i>Do not uninstall.</i></p> <p>Install Full GMID-X v1.0 or v1.1. See “Performing a Full Installation” on page 36.</p> | <ol style="list-style-type: none"> 1. Uninstall GMID v3.x. See “Uninstalling the GMID v3.x Software” on page 24. 2. Install Client GMID-X v1.0 or v1.1. See “Performing a Client Installation” on page 46. |

Uninstalling the GMID v3.x Software

IMPORTANT! *For 31xx or 3730 Data Collection Full Installations only:* If you are installing the *full* GeneMapper® ID-X Software v1.0 or v1.1 on a 31xx or 3730 Data Collection Computer, *do not uninstall* the GMID v3.x software. The GeneMapper® ID-X Software v1.0 or v1.1 installer automatically overrides the existing GMID v3.x software and installs the *full* GeneMapper® ID-X Software v1.0 or v1.1 for you.

IMPORTANT! For all other installations, you must restart your computer before uninstalling. If you forget to do so, upon running the uninstall, the computer prompts you the restart your computer (click **Cancel** to exit the uninstall). If you have already restarted your computer, click **OK** to continue the uninstall.

Follow the steps below to *uninstall* GeneMapper® ID Software v3.x:

1. Log in to the computer using a user account with administrative privileges.
2. Close the GeneMapper® ID Software v3.x and all other applications, then restart your computer before proceeding with the uninstall.

IMPORTANT! If you do not restart your computer before proceeding with the uninstall and/or if the uninstall does not complete successfully the first time, perform this procedure a second time. If the uninstall still does not work, have a knowledgeable technician review and proceed with “[Running the Clean-up Utility](#)” on [page 102](#) or contact Applied Biosystems Technical Support to proceed.

3. Select **Start ▶ Control Panel**.

4. Double-click **Add or Remove Programs**.

Note: If a message indicates that you need administrator rights or privileges, log off the computer, then log on again as a user with administrator privileges.

5. Select the GeneMapper® *ID* Software.
6. Click **Change/Remove**. The InstallShield Wizard opens.
7. In the Welcome page, select **Remove**, then click **Next**.
8. At the prompt, click **OK** to verify the uninstall.
9. When the Maintenance Complete window opens, select **Yes, I want to restart my computer now**.

IMPORTANT! For the uninstallation to complete, you must restart the computer. *Do not install* GMID-X v1.0 or v1.1 software without first restarting your computer.

IMPORTANT! After you restart your computer, a command prompt window may be displayed. This window can take up to 40 seconds to close. *Do not close* this window. Wait for the process to complete on its own. If you close the window, the GMID-X software may not install correctly.

Enabling Autoanalysis

Use autoanalysis if you want to communicate with the Data Collection software for the collection and processing of .fsa sample files, and for data analysis.

For installations in which you want autoanalysis enabled, follow the procedures below:

- Start the Data Collection Software to enable autoanalysis during the full and client installation procedure (see [page 26](#)).
- Retrieve the Data Collection Computer name to enable autoanalysis during the full and client installation procedure (see [page 29](#)).
- Retrieve the Host Computer name to enable client access to the host computer during the client installation procedure (see [page 30](#)).

IMPORTANT! The Data Collection software must be running during installation if you wish to run autoanalysis.

Starting the Data Collection Software on the 31xx Data Collection Computer

Start the Data Collection Software to Enable Autoanalysis


Start the Data Collection Software on the 31xx Data Collection computer *before* installing the full GeneMapper® *ID-X* Software v1.0 or v1.1 on your computer for a *remote autoanalysis configuration*, or on the 31xx Data Collection computer for a *co-installation*.

Note: *For 310 Data Collection Computers only:* The following Data Collection Software sections are not applicable. Proceed to “[Performing a Full Installation](#)” on page 36 or “[Performing a Client Installation](#)” on page 46.

Starting the Data Collection Software

Use the following procedure to enable autoanalysis only.

1. On the 31xx Data Collection computer, select **Start ▶ All Programs ▶ Applied Biosystems ▶ Data Collection ▶ Run <Data Collection version>**, where <Data Collection version> is:
 - 3100/3100-*Avant*[™] Data Collection v2.0,
 - or*
 - 3130/3130*xI* Data Collection v3.0.
2. After the Service Console opens, wait until all four symbols change to green squares .




Note: If the software services do not start automatically, click **Start All**.


3. If the Data Collection Software requires a password, a login dialog box opens. Enter the Login Name and Password, then click **OK**.

Note: If you do not know the Data Collection computer Login Name or Password, contact your Data Collection computer administrator.

4. Verify that the Data Collection software started without errors:
 - a. In the Service Console, right-click the square next to Data Service, then select **Show Console** to display the Data Service output message box.




- b. Verify that no errors are displayed in the lower pane of the message box, then close the Data Service output message box.



Obtaining the Data Collection Computer Name

If you are installing the client GeneMapper® *ID-X* Software v1.0 or v1.1 on a non-Data Collection Computer to enable remote autoanalysis, you must first obtain the name of the 31xx Data Collection computer:

1. On the 31xx Data Collection Computer, right-click **My Computer**, select **Properties**, then select the **Computer Name** tab to view the full computer name.




Computer
name

2. Make a note of the computer name and enter this name when prompted during installation of the client GeneMapper® *ID-X* Software.

Obtaining the GeneMapper® ID-X Software Host Computer Name

Before you install the client GeneMapper® ID-X Software, obtain the name of the multi-user database host computer to enable access to the host computer from a network:

1. On the multiuser database host computer, right-click **My Computer**, select **Properties**, then select the **Computer Name** tab.



2. Make a note of the computer name and enter this name when prompted during installation of the client GeneMapper® ID-X Software.

Requirements for Upgrading GMID-X v1.0/1.0.1 to GMID-X v1.1

To upgrade your target computer (client or full) to v1.1, make certain that either GMID-X v1.0 or v1.0.1 is already installed on your system. The upgrade installer works only if it detects either of these previous versions.

IMPORTANT! When you upgrade on a client computer, you must also upgrade the computer that hosts the full software version on the same local area network.

Note: You can reuse your existing password when upgrading from GMID-X v1.0 or v1.0.1 to v1.1.

3

To upgrade from GMID-X v1.0 or v1.0.1 to v1.1:

1. Insert the *GeneMapper® ID-X Software Version 1.1 Upgrade Install CD* into your CD drive.
2. In the Welcome Window, click **Next**.
3. Read the Release notes, then click **Next**.
4. Click **Next** to install the Software to the default location.
or
To install the GeneMapper® ID-X Software to a drive different from the default drive, select the destination folder.
5. When the software is installed, remove the installation CD from the CD drive.


The upgrade installation takes approximately 5 to 10 minutes, depending on the computer being upgraded.

4

Performing a Full Software Installation

This chapter covers:

- Overview 34
- User Account Requirements for Full Installations 34
- Performing a Full Installation 36



Overview

This chapter describes how to perform a full installation of the GeneMapper® *ID-X* Software Version 1.0 or Version 1.1 (GMID-X v1.0 or v1.1) on your target computer.

Note: If you have GeneMapper® *ID* Software v3.x (GMID v3.x) on your computer and haven't already performed pre-installation tasks, go back to [Chapter 3, “Performing GMID-X v1.0 Pre-Installation and GMID-X v1.1 Upgrade Procedures,”](#) on [page 15](#) for additional instructions. The installation may fail if these tasks are not performed first.

User Account Requirements for Full Installations

Log On Requirements


To perform a new full installation of GMID-X v1.0 or v1.1, you must:

- Log on to the local computer (not a network domain).
- Use an Administrator account (unrestricted access).

Note: After the software is installed, you can run the GMID-X v1.0 or v1.1 without using an Administrator account.

Verifying User Accounts

1. On the desktop, select **Start ► Control Panel**.
2. In the Control Panel window, double-click **User Accounts**.
3. In the Users tab, verify that the user account belongs to the Administrators group and the domain name is the same as the computer name.



Belongs to the Administrators

Domain is the same as the computer name

Performing a Full Installation

A full GeneMapper® *ID-X* Software v1.0 or v1.1 installation:

- Installs the GMID-X software and database.
- Enables the GMID-X v1.0 or v1.1 software to serve as the host for multiple client computers needing access to the database.

Note: It is not recommended that the Data Collection computer function as the host computer due to performance constraints.

Full Installation Requirements

Important Reminders

Below is a list of items to take into consideration prior to installing GMID-X v1.0 or v1.1.

- *DO NOT* cancel an installation before the installation process is finished. If you do so, or if the installation otherwise stalls, you may need technical support to perform an uninstall. To uninstall your installation, see [“Uninstalling the GMID-X Software” on page 100](#). If uninstallation fails, contact your local Applied Biosystems technical support organization.
- In a co-installation (Data Collection computer with Data Collection Software), you are limited to one client connection. To allow more client connections to the database, install the full version GeneMapper® *ID-X* Software (software and database) on a computer other than the Data Collection (DC) computer.
- A full installation on a 31xx DC Computer is only possible if the DC Computer is running DC Software v2.0 or higher.
- When installing GMID-X v1.0 or v1.1 on a 31xx DC Computer that has GMID v3.2 software, you do not need to uninstall the v3.2 software; however, you must recreate all the files associated with autoanalysis (e.g. Results Groups) because old files created as part of the GMID v3.2 integration cannot be reused. See [Chapter 7, “Setting Up Autoanalysis.”](#)


- When you install GMID-X v1.0 or v1.1 on a 31xx DC Computer or on a non-DC computer (computer without the Data Collection software) for remote autoanalysis, the data collection software must be running on the DC Computer. For more information on starting the data collection software, see [page 26](#).
- To set up autoanalysis, obtain the computer name of the DC Computer: Right-click **My Computer** on the DC computer, select **Properties**, then select the **Computer Name** tab to view and record the full computer name. You are required to enter this during installation (see [page 29](#)).

Installing the Full GeneMapper® ID-X Software

Full Installation Procedure (DC and Non-DC Computers)

To install the full GMID-X v1.0 or v1.1 software onto your computer:

1. Insert the first disk of the *GeneMapper® ID-X Software Version 1.0* (or *v1.1*) *Full Install* CDs into the CD drive (if not inserted during use of the Export Utility).
2. On the splash screen, double-click **Install GeneMapper ID-X v1.0** (or **v1.1**) **Application** to start the installer.
3. Close all other applications and windows (except the Data Collection application), then click **OK** to close the following message.



4. In the Welcome window, click **Next**.
5. Review the installation requirements status, then click **Next**.

6. Select **Remote Autoanalysis** or **Stand-alone** for type of installation, then click **Next**:


- **Remote Autoanalysis** – Select this option to set up remote autoanalysis. In this configuration, the target computer is networked to the 31xx DC computer. The autoanalysis manager alerts GMID-X when a run is complete. GMID-X then opens, and you can import the .fsa files into a project and analyze them.
- **Stand-Alone** – Select this option if you *do not* want to set up remote autoanalysis or if your DC computer is not networked to the target computer during installation.

IMPORTANT! After the GMID-X v1.0 or v1.1 Software is installed, you cannot switch from stand-alone to remote autoanalysis. To change the configuration, you must uninstall, then reinstall the software.


Note: For the 31xx Data Collection computer, you do not have the option of selecting stand-alone or remote autoanalysis. You can optionally set up local autoanalysis later. See [“Setting Up Autoanalysis” on page 68](#) for instructions.

Note: For the 310 Data Collection Computer, autoanalysis (local or remote) is not available. Proceed to [step 8](#).

Note: If you select **Stand-alone**, skip [step 7](#) and [step 13](#).




7. (Remote Autoanalysis only) If you select Remote Autoanalysis, enter the full name of the Data Collection computer that you intend to have support remote autoanalysis (see [page 29](#)), then click **Next**.




8. Read the release notes, then click **Next**.

9. Click **Next** to install the Software to the default location. To install the GeneMapper® ID-X Software to a drive different from the default drive, select the destination folder:
 - a. Verify that the boot drive (drive on which the Windows operating system is installed) contains at least 200 MB of available free space.
 - b. Select a destination folder with at least 7 GB of free space, then click **Next**.

IMPORTANT! If you choose an invalid destination such as A: drive or CD-DVD ROM and click **Next**, an error message is displayed saying the drive is invalid. When you click **OK** the installation is stopped. To recover, click **Cancel** and start the installation over again, using a valid destination folder.



10. Select a Setup Type if desired: either Desktop Short Cut and/or Quick Launcher or deselect both, then click **Next**.




11. Read the current settings, then click **Next** to start the installation.

IMPORTANT! During installation, a DOS window may open while DOS commands execute. Do not delete, close, or click the DOS window. If you accidentally click the DOS window, press the **Esc** key to exit the window.

12. When prompted, insert Disk 2 and follow instructions.

13. (Remote Autoanalysis only) When GeneMapper® *ID-X* Software installation is complete, if you were running the Data Collection Software, stop the software on the 31xx Data Collection computer by clicking **Stop All** in the Service Console.



14. Remove the installation CD from the CD drive.
15. In the InstallShield Wizard Complete page, select **Yes, I want to restart my computer now**, then click **Finish**.

When you finish the installation, proceed to [Chapter 6, “Performing Post-Installation Procedures,”](#) on [page 53](#).


To set up autoanalysis after performing post-installation procedures, proceed to [Chapter 7, “Setting Up Autoanalysis,”](#) on [page 61](#).

5

Performing a Client Software Installation

This chapter covers:

- Overview 44
- User Account Requirements for Client Installations. 44
- Performing a Client Installation 46
- Resetting Permissions for the Installation. 52



Overview

This chapter describes how to perform a client installation of the GeneMapper® *ID-X* Software Version 1.0 or Version 1.1 on your target computer.

Note: If you have GeneMapper® *ID* Software v3.x on your computer and haven't already performed pre-installation tasks, go back to [Chapter 3, “Performing GMID-X v1.0 Pre-Installation and GMID-X v1.1 Upgrade Procedures,”](#) on [page 15](#) for additional instructions. The installation may fail if these tasks are not performed first.

User Account Requirements for Client Installations

Log In Requirements


To perform a client installation of the GeneMapper® *ID-X* Software:

- Confirm that the target computer is connected to the database host computer when you log in. If there is no login window that connects you to the local area network served by the host, or if you are not automatically connected, contact your lab's on-site technical support department.
- Use an Administrator account (unrestricted access).

Note: After the software is installed, you can run the GeneMapper® *ID-X* Software without using an Administrator account.

Verifying User Accounts

1. On the desktop, select **Start ► Control Panel**.
2. In the Control Panel window, double-click **User Accounts**.
3. In the Users tab, verify that the user account belongs to the Administrators group and the domain name is the same as the computer name.



Belongs to the Administrators

Domain is the same as the computer name

Performing a Client Installation

A client GeneMapper® *ID-X* Software v1.0 or v1.1 installation:

- Allows the user access to the GMID-X v1.0 or v1.1 software without a database installed on the same computer.
- Connects to one or more database hosts for access to the data stored in the database and saves data objects back to the database.

Client Installation Requirements

Important Reminders

Below is a list of items to take into consideration before you perform a client installation:


- *DO NOT* cancel an installation before the installation process is finished. If you do so, or if the installation otherwise stalls, you may need technical support to perform an uninstall. To uninstall your installation, see [“Uninstalling the GMID-X Software” on page 100](#). If uninstallation fails, contact your local Applied Biosystems technical support organization.
- At least one full installation must be performed on another computer for it to act as the database host for the client.
- Users are required to have a valid user name and password with an account on the database host computer.
- You must obtain the name of the database host computer to enable access to the host computer from a network: Right-click **My Computer** on the database host computer, select **Properties**, then select the **Computer Name** tab to view and record the full computer name. You are required to enter this name during installation (see [page 30](#)).
- When you install GMID-X v1.0 or v1.1 on a 31xx Data Collection computer or on a non-DC computer for remote autoanalysis, the data collection software must be running on the Data Collection computer (see [page 26](#).)

- If you choose to set up autoanalysis, obtain the Computer Name of the Data Collection computer: Right-click **My Computer** on the Data Collection computer, select **Properties**, then select the **Computer Name** tab to view and record the full computer name. You are required to enter this during installation (see [page 29](#)).

Installing the Client GeneMapper® ID-X Software


To install the client GeneMapper® ID-X Software Version 1.0 or Version 1.1 on the target computer:

1. Insert the *GeneMapper® ID-X Software Version 1.0 (or Version 1.1) Client Install* CD into the CD drive.
2. On the splash screen, double-click **Install GeneMapper ID-X v1.0 (or v1.1) Client** to start the installer.
3. Be sure that the Data Collection software application is running but close all other applications and windows, then click **OK** to close the following message.




4. In the Welcome window, click **Next**.
5. Review the installation requirements status, then click **Next**.
6. Read the release notes, then click **Next**.

7. When prompted, enter the GeneMapper® ID-X Software (GM) host (server) computer name, then click **Next**.



8. After the installer establishes a connection with the multi-user database computer, click **OK** to continue.



9. Select **Remote Autoanalysis** or **Stand-Alone** for type of installation, then click **Next**.
- **Remote Autoanalysis** – Select this option to set up remote autoanalysis. In this configuration, the target computer is networked to the 31xx Data Collection computer. The autoanalysis manager alerts GMID-X when a run is complete. GMID-X then opens, and you can import the .fsa files into a project and analyze them.


- **Stand-Alone** – Select this option if you *do not* want to set up remote autoanalysis or if your DC computer is not networked to the target computer during installation.

IMPORTANT! After the GMID-X v1.0 or v1.1 Software is installed, you cannot switch from stand-alone to remote autoanalysis. To change the configuration, you must uninstall, then reinstall the software.


Note: For the 31xx Data Collection computer, you do not have the option of selecting stand-alone or remote autoanalysis. You can optionally set up local autoanalysis later. See “[Setting Up Autoanalysis](#)” on page 68 for instructions.

Note: For the 310 Data Collection Computer, autoanalysis (local or remote) is not available. Proceed to [step 11](#).

Note: If you select **Stand-alone**, skip [step 10](#) and [step 13](#).




10. (Remote Autoanalysis) If you selected Remote Autoanalysis, enter the full name of the 31xx Data Collection computer that will support the remote autoanalysis configuration, then click **Next**.



11. To install the GeneMapper® *ID-X* Software to a drive on the client computer other than the default drive, select the destination folder and verify that the boot drive (drive on which the Windows operating system is installed) contains at least 250 MB of available free space.

IMPORTANT! If you choose an invalid destination such as A: drive or CD-DVD ROM and click **Next**, an error message is displayed saying the drive is invalid. When you click **OK** the installation is stopped. To recover, click **Cancel** and start the installation over again, using a valid destination folder.



12. Read the current settings, then click **Next** to start the installation.

IMPORTANT! During installation, a DOS window may open while DOS commands execute. Do not delete, close, or click the DOS window. If you accidentally click the DOS window, press the **Esc** key to exit the window.

13. When installation of the GeneMapper® *ID-X* Software is complete, if you had the 31xx Data Collection Software running, stop the software on the Data Collection computer by clicking **Stop All** in the Service Console.
14. Remove the installation CD from the CD drive.
15. In the InstallShield Wizard Complete window, select **Yes, I want to restart my computer now**, then click **Finish**.
16. Repeat the procedure for additional client installations.

When you finish all client installations, proceed to [Chapter 6, “Performing Post-Installation Procedures,”](#) on page 53.

To set up autoanalysis after performing post-installation procedures, proceed to [Chapter 7, “Setting Up Autoanalysis,”](#) on page 61.

Resetting Permissions for the Installation


Resetting Permissions on Folders and Files

The client installer specifies folder privileges at the time of the GeneMapper® *ID-X* Software installation.

However, the NTFS permissions set by the system administrator may not allow the GeneMapper® *ID-X* Software to access the GeneMapper folder following installation. Before you run the software, check the permissions to make sure that the GeneMapper® *ID-X* Software can access the GeneMapper folder.

If you do not have the required permissions, you may see a “can't access file” error message when using the software. To set the privileges manually, either before running the software or in the event of an error message:

1. Right-click the AppliedBiosystems folder.
2. Select **Properties**, then click the **Security** tab.
3. In the top pane, select **Users**.
4. In the bottom pane, select **Full Control**.




6

Performing Post-Installation Procedures

This chapter covers:

- Registering the GeneMapper® ID-X Software 54
- Logging Into GMID-X After Installation 55
- Connecting to a New Host. 56
- Logging Out of the GeneMapper® ID-X Software 58
- Additional Post-Installation Tasks. 59





Registering the GeneMapper® ID-X Software

Register the GeneMapper® ID-X Software v1.0 or v1.1 the first time you start the software on each computer. The registration code is on the registration card packaged with the software.

Note: If you upgraded from GeneMapper® ID-X Software v1.0/v1.0.1 to v1.1, you do not need to re-register. Your registration code remains the same.

Enter the same registration code on each computer, based on the number of users each registration code supports.


1. On the desktop, double-click  **GeneMapper ID-X v1.0** (or **v1.1**) or select **Start ▶ All Programs ▶ Applied Biosystems ▶ GeneMapper ▶ GeneMapper ID-X v1.0** (or **v1.1**).
2. Complete the Registration dialog box:




The image shows a Windows-style dialog box titled "GeneMapper". Inside the dialog, the title "Product Registration" is centered. There are three text input fields: "Your Name:", "Organization:", and "Registration Code:". Below these fields are two buttons: "OK" and "Quit".

3. Click **OK**.

Logging Into GMID-X After Installation

1. If the Login dialog box is not displayed in the desktop, double-click  **GeneMapper ID-X v1.0** (or **v1.1**) or select **Start** ▶ **All Programs** ▶ **Applied Biosystems** ▶ **GeneMapper** ▶ **GeneMapper ID-X v1.0** (or **v1.1**).
2. In the Login to GeneMapper ID-X dialog box, enter or select “*gmidx*” as the administrator account user name and the temporary password, “**password**”, then click **OK**:



3. If prompted, specify a new password. Keep a record of the new password. Click **OK**.

Note: The *gmidx* account requires you to set a new password when you log in for the first time. If you follow this procedure after first log in, the passwords for the accounts listed may be different from the passwords listed in this guide.

IMPORTANT! Do not lose the password for the *gmidx* account. Applied Biosystems cannot retrieve these passwords.

4. Read the license agreement, then click **Yes** if you agree with the terms.
5. Review the license agreement and product warranty, then click **Accept**.

Connecting to a New Host


If you install more than one multi-user database computer, any client computer can connect to any of the databases if the client is on the same network as the database computers.

Users must set up separate user accounts on each multi-user database computer that they wish to access.


Note: For detailed information on setting up user accounts, see the *GeneMapper® ID-X Software Version 1.0 Administrator's Guide*.

To make the multi-user database computers available to client computers, perform the following steps for each client computer:

1. Click **New Host** in the Login to GeneMapper ID-X dialog box.



2. Complete the information in the New Host dialog box:




- a. In the **Host Name** field, enter the full computer name or IP address of the multi-user database host you want to access. If the GeneMapper® *ID-X* Software cannot connect to the database host, the error message “You have entered an invalid host” is displayed. Click **OK** to exit, then re-enter the database host information.
- b. Select the applicable machine type.
- c. Click **OK**.


The Login window is updated to include the name of the new database host. The User Name list reflects the user accounts on the new database host. See the *GeneMapper® ID-X Software Version 1.0 Administrator's Guide* for additional information.


Logging Out of the GeneMapper® ID-X Software

You can close the GeneMapper® ID-X Software by:

- Selecting **File ▶ Logout**.
- Selecting **File ▶ Exit**.
- Clicking  (Close).

When you log out of the GeneMapper® ID-X Software with a project that has unsaved changes, the software prompts you to save or discard the changes you made since the last time you updated the project.



If you used the **File ▶ Logout** command, the login window re-opens after you select either **Yes** or **No**. The GeneMapper ID-X Software closes completely if you use either **File ▶ Exit** or  (Close).

Additional Post-Installation Tasks

Importing Application Data

If you exported application data from GeneMapper® *ID* Software v3.2.x before you installed GeneMapper® *ID-X* Software v1.0 or v1.1, you may wish to import the data objects (projects, analysis methods, etc.) into the GeneMapper® *ID-X* Software v1.0 or v1.1 database host computer.

Select **Tools ► GeneMapper *ID-X* Manager**. Then in the GeneMapper *ID-X* Manager window, click each tab to import the corresponding objects into GMID-X.

Creating User Accounts and Configuring Administrative Functions

For information on creating user accounts and configuring the Security Manager, Audit Manager, and E-Signature Administrator, see the *GeneMapper® ID-X Software Version 1.0 Administrator's Guide*.

Setting Up Autoanalysis


If you installed the GeneMapper® *ID-X* Software as either a remote autoanalysis configuration, co-installation, or Data Collection client installation, follow the procedures in [Chapter 7, “Setting Up Autoanalysis,”](#) on [page 68](#).



Setting Up Autoanalysis

This chapter covers:

- Overview of Autoanalysis Process 62
- Determining if Your Computer Can Perform Autoanalysis .. 64
- Autoanalysis Instrument and Software Configurations 65
- Overview of Autoanalysis Setup 66
- Setting Up Autoanalysis 68
- Recommended Data Processing Tasks 83



Overview of Autoanalysis Process

Autoanalysis Process

The autoanalysis process occurs in the following sequence:

1. You close the GeneMapper® *ID-X* Software and start the Autoanalysis Manager.
2. On the Data Collection (DC) computer, you select the plate to run and schedule the run.
3. When the DC Software finishes a run as specified in the results group, the data (.fsa sample files) are stored in the destination specified by the DC Software results group.
4. The DC Software sends the message “Run Completed” to the Autoanalysis Manager. The Autoanalysis Manager searches for jobs every two minutes.
5. When the Autoanalysis Manager receives a “Run Completed” message, it adds the job to its job queue.
6. The Autoanalysis Manager:
 - a. Opens the GeneMapper® *ID-X* Software.
 - b. Creates a project using the name of the Data Collection results group.
7. The GeneMapper® *ID-X* Software:
 - a. Imports the .fsa sample files from the location specified in the Data Collection results group.
 - b. Uses the settings specified in the Data Collection plate record (panels and bins, analysis method, size standard) to analyze the .fsa sample files.
 - c. Saves the project, then associates the data in the project with the gmidx security group (for information on security groups, see the *GeneMapper® ID-X Software Version 1.0 Administrator's Guide*).

- d. Stores the results in the default host GeneMapper® *ID-X* Software database.

Note: The default host is the GeneMapper *ID-X* Software computer that is specified in the Login dialog box of the computer that performs autoanalysis. For instructions on changing the default host, see [“Connecting to a New Host” on page 56](#).

8. The Autoanalysis Manager:
 - a. Closes the GeneMapper *ID-X* Software.
 - b. Updates the status of the run in the Autoanalysis Manager job queue.

For information on the Autoanalysis Manager, see the *GeneMapper® ID-X Software Online Help*.

Determining if Your Computer Can Perform Autoanalysis

Autoanalysis Software Components

An autoanalysis configuration uses:

- **Data Collection Software** – Runs the instrument and collects data from samples. See [“Instrument and software options for autoanalysis” on page 65](#) to determine if your Data Collection computer supports autoanalysis.
- **Autoanalysis Manager** – Communicates between the Data Collection Software and the GeneMapper® *ID-X* Software. It queues jobs and tracks the status of their processing. A job can contain a single run or multiple runs as part of a result group.
- **GeneMapper® *ID-X* Software v1.0 (or v1.1)** – Creates a project, imports .fsa sample files, analyzes the data, then saves the results in the multi-user database.

Figure 2 shows the relationship and interaction between the autoanalysis components.




Figure 2 Relationship between autoanalysis components




Autoanalysis, which can be performed from a multi-user database computer or a client, requires the Autoanalysis Manager.

To determine if your computer contains the Autoanalysis Manager, select **Start ▶ Applied Biosystems ▶ Autoanalysis Manager ▶ Autoanalysis Manager 3.0**.

Autoanalysis Instrument and Software Configurations

With autoanalysis, the GeneMapper® *ID-X* Software v1.0 or v1.1 can automatically analyze data generated on the following systems:

Table 7 Instrument and software options for autoanalysis

| Genetic Analyzer Instrument | Compatible Data Collection Software and Operating Systems |
|--|--|
|  <p>Applied Biosystems 3130/3130xl Genetic Analyzer</p> | <ul style="list-style-type: none"> • 3130 Data Collection v3.0/ 3130xl Data Collection v3.0 • Windows XP, SP 2 or later |
|  <p>ABI PRISM® 3100/3100-Avant™ Genetic Analyzer</p> | <ul style="list-style-type: none"> • 3100 Data Collection v2.0/ 3100-Avant™ Data Collection v2.0 • Windows 2000, SP 3 or later |
|  <p>Applied Biosystems 3730 Genetic Analyzer‡</p> | <ul style="list-style-type: none"> • 3730 Data Collection v3.0 • Windows XP, SP 2 or later |

‡ The 3730 Genetic Analyzer does not include the 3730xl 96-capillary configuration. The 48-capillary configuration is validated for analysis of Identifiler® data generated from single-source data samples only.

Overview of Autoanalysis Setup

Local and Remote Autoanalysis

With the analyzers listed in [Table 7](#), you can configure the Data Collection Software and the GeneMapper® *ID-X* Software v1.0 or v1.1 to perform data collection and then analyze the data automatically, without user interaction.

The GeneMapper® *ID-X* Software v1.0 or v1.1 can be a co-installation or a client installation on the DC computer (for local autoanalysis). Alternatively, the software can be installed on a non-DC computer and configured for remote autoanalysis with the DC computer, as shown [Figure 3](#).

Note: If you install the GeneMapper® *ID-X* Software on the same computer as the Data Collection software (co-installation), you must manually import the panel and bin files. See the *GeneMapper® ID-X Software Online Help* for information on this procedure.




Figure 3 Examples of autoanalysis configurations

Data Processing Options

The Autoanalysis feature of the GeneMapper® *ID-X* Software allows you to set up a Data Collection run, then process the data automatically in the GeneMapper® *ID-X* Software. Autoanalysis can occur:

- **Locally** – from the DC computer when the Data Collection software and GeneMapper® *ID-X* Software are installed on the same computer.
- **Remotely** – from the non-DC computer when the Data Collection software and GeneMapper® *ID-X* Software are installed on different computers.

The Autoanalysis Manager is automatically installed when you install the GeneMapper® *ID-X* Software on a:

- DC computer as a co-installation or client installation
- Non-DC computer for remote autoanalysis

Autoanalysis Manager is not installed when you install the GeneMapper® *ID-X* Software on a non-DC computer in a stand-alone configuration.

If the Autoanalysis Manager is not installed on your computer and you want to use the computer for autoanalysis, uninstall the GeneMapper® *ID-X* Software (see [“Uninstalling the GMID-X Software” on page 100](#)), then reinstall for Remote Autoanalysis.

Before Setting Up Autoanalysis

Before you set up autoanalysis:

- Ensure that the Data Collection Software and GeneMapper® *ID-X* Software are installed in either of the following configurations:




Figure 4 Configurations necessary for autoanalysis

- If the software is installed on different computers, ensure that the two computers are connected to the same network.
- Ensure that the remote GeneMapper® *ID-X* Software computer has the Autoanalysis Manager installed: Select **Start ▶ Applied Biosystems ▶ Autoanalysis Manager ▶ Autoanalysis Manager 3.0**.
- If the Autoanalysis Manager is not installed, you must uninstall the GeneMapper® *ID-X* Software, then reinstall the GeneMapper® *ID-X* Software and select the Remote Autoanalysis option (see [“Uninstalling the GMID-X Software” on page 100](#)).
- Import or create analysis methods, panels, bins, and size standards on the multi-user database host computer (see the *GeneMapper® ID-X Software Online Help*).
- If you install the GeneMapper® *ID-X* Software on the same computer as the Data Collection software (co-installation), you must manually import the panel and bin files. See the *GeneMapper® ID-X Software Online Help* for information on this procedure.

Setting Up Autoanalysis

The autoanalysis set-up process consists of four procedures as shown in the following table.

| Procedure | Action | See page |
|-----------|--|----------|
| 1. | Create a User Account. | 69 |
| 2. | In the Data Collection Software, create an instrument protocol for autoanalysis. | 70 |
| 3. | In the Data Collection Software, create a results group for autoanalysis. | 72 |
| 4. | In the Data Collection Software, create a plate record for autoanalysis. | 80 |

Procedure 1: Creating a User Account for Autoanalysis

IMPORTANT! The user account for autoanalysis must exist on the computer that is specified as the host computer, which may be different from the computer on which you run autoanalysis and store data files (see [Figure 5 on page 84](#)).

Create a user account for autoanalysis with the following settings (for information on creating user accounts, see the *GeneMapper® ID-X Software Version 1.0 Administrator's Guide*). In the GeneMapper ID-X Software Security Manager New User screen:

1. Deselect **Pre-expire**.

IMPORTANT! If the account is set to pre-expire, you must log in with the account one time to change the password before you use the account for autoanalysis.

2. Select **Scientist** profile.

Note: If you create custom profiles for your system, create a profile for autoanalysis with the following minimum privileges: Read privileges for Panel Manager, Size Standard, Analysis Method; Read and Update privileges for creating and analyzing projects.

3. For Full Name, enter **Remote Autoanalysis, Data Collection Client**, or **Co-installation**.

4. Deselect **Show EULA**.

5. Associate the user account with any user group. The project that autoanalysis creates, analyzes, and saves is associated with the GeneMapper® ID-X Software security group, which is accessible by all user groups.

When you are finished, go to Procedure 2 and create an instrument protocol.

Procedure 2: Creating an Instrument Protocol for Autoanalysis

An instrument protocol contains all the settings needed to run the instrument, including the protocol name, type of run, run module, and dye set.


IMPORTANT! Create new instrument protocols after installing or reinstalling the GeneMapper® *ID-X* Software v1.0 or v1.1.

For more detailed information about instrument protocols, see your instrument user guide.

Creating Instrument Protocols

1. In the Navigation pane of the Data Collection Software, select **GA Instruments** ► **<Instrument Name>** ► **Protocol Manager**.
2. In the Instruments Protocols section, click **New** to open the Protocol Editor.

2



3. In the Protocol Editor, enter instrument protocol information:
 - a. Enter a name for the protocol.
 - b. (Optional) Enter a description for the protocol.
 - c. Select **REGULAR** in the Type drop-down list.

- d. Select the **HIDFragmentAnalysis36_POP4_1** run module in the Run Module drop-down list.
- e. Select the appropriate dye in the Dye Set drop-down list (see the table that follows to determine the dye).

Protocol Editor

Name: HID 3a

Description: 3b

Type: REGULAR 3c

Run Module: HIDFragmentAnalysis36_POP4_1 3d

Dye Set: G5 3e

OK Cancel

| Kit | Dye Set |
|---|---------|
| AmpF ℓ STR $^{\text{®}}$ COfiler $^{\text{®}}$ Kit AmpF ℓ STR $^{\text{®}}$ Profiler Plus $^{\text{®}}$ Kit AmpF ℓ STR $^{\text{®}}$ Profiler Plus ID Kit AmpF ℓ STR $^{\text{®}}$ SGM Plus $^{\text{®}}$ Kit Other 4-Dye AmpF ℓ STR $^{\text{®}}$ Kits | F |
| AmpF ℓ STR $^{\text{®}}$ SEfiler $^{\text{™}}$ Kit AmpF ℓ STR $^{\text{®}}$ Identifiler $^{\text{®}}$ Kit AmpF ℓ STR $^{\text{®}}$ Yfiler $^{\text{®}}$ Kit AmpF ℓ STR $^{\text{®}}$ Minifiler $^{\text{™}}$ Kit Other 5-Dye AmpF ℓ STR $^{\text{®}}$ Kits | G5 |

- f. Click **OK**.

When you are finished, go to Procedure 3 to create a Results Group for autoanalysis.

Procedure 3: Creating a Results Group for Autoanalysis

Overview of Results Groups

Create a Results Group for each computer that performs autoanalysis.

A Results Group specifies the file names and storage location (shared folder) for .fsa sample files, and the user name and password for the GeneMapper® *ID-X* Software computer that performs autoanalysis. After you create the Results Group, the software can perform autoanalysis.

Create new Results Groups after installing or reinstalling the GeneMapper® *ID-X* Software v1.0 or v1.1.

IMPORTANT! The maximum number of .fsa sample files that the Data Collection Software can assign to a Results Group, when using the Analysis Type GeneMapperIDX-Generic or GeneMapperIDX-*<computer name>*, is 7000. After you collect 7000 .fsa sample files in a Results Group, you must create a new results group.

For more detailed information about results groups, see your specific instrument user guide.

Allelic Ladder Requirements and Results Groups

IMPORTANT! If a run folder does not contain an allelic ladder, the samples are reported as off ladder (OL) in the GeneMapper® *ID-X* Software. For more information, see the *GeneMapper® ID-X Software Online Help*.

To set up the Data Collection Software so at least one allelic ladder is present in each run folder, create a results group that *groups by plate* as described below, instead of by run (which is the default).

Creating a Results Group for Autoanalysis

For each GeneMapper® *ID-X* Software client computer that runs autoanalysis:

1. In the Navigation pane of the Data Collection Software, select **GA Instruments ► Results Group**.
2. Click **New** to open the Results Group Editor, then click the tabs of the results groups and enter the information as described below.

Entering Results Group Information

1. In the General tab, enter the Results Group Name. This name is assigned as the default autoanalysis project name (all results for autoanalysis are stored in this GeneMapper® *ID-X* Software project). The Results Group name must be unique to the Data Collection Software.
2. (Optional) Enter the Results Group Owner.
3. (Optional) Enter a Results Group Comment.

The screenshot shows the 'Results Group Editor' window with the 'General' tab selected. The 'Results Group Name' field is highlighted with a blue selection bar and labeled with a '1'. The 'Results Group Owner' field is labeled with a '2'. The 'Results Group Comment' field is labeled with a '3'. The 'OK' and 'Cancel' buttons are at the bottom right.

Entering Analysis Parameters

1. Click the **Analysis** tab.
2. For the Analysis Type, select **GeneMapperIDX-
<computer name>** (the name of the host or client computer on which you intend to run autoanalysis).

Note: The computer names are added to the Analysis Type list when you install the GeneMapper® *ID-X* Software as a co-installation or client installation on the Data Collection computer, or choose the remote autoanalysis option for the non-DC computer.


Note: You cannot perform autoanalysis if you select **GeneMapperIDX-Generic**.

3. Select **Do Autoanalysis**.
4. Select **Results Group Entry Complete** to start autoanalysis *after* all samples that use the same results group (all samples on the plate) have been run and to ensure that each run folder contains an allelic ladder.

IMPORTANT! If a run folder does not contain an allelic ladder, the samples are reported as off ladder (OL) in the GeneMapper® *ID-X* Software. For more information, see the *GeneMapper® ID-X Software Online Help*.


5. In the Login ID and Password fields, enter a GeneMapper® *ID-X* Software user account username and password for the computer on which you intend to run autoanalysis. This account must exist on the computer specified as the default host. The default host computer may not be the computer on which you perform autoanalysis (the computer running Autoanalysis Manager) (see [Figure 5 on page 84](#)).

IMPORTANT! The GeneMapper *ID-X* Software requires a user account with a Scientist or higher security profile to perform autoanalysis (see [“Procedure 1: Creating a User Account for Autoanalysis” on page 69](#)).




Specifying the Sample File Storage Location

1. Click the **Destination** tab.
2. Specify the location in which to store the .fsa sample files
3. Select **Use Custom Location**.




4. Set the destination as described in the following table.

| For... | Do the Following... |
|---------------------|--|
| Local autoanalysis | <ol style="list-style-type: none"> 1. Click Browse, then select a destination. 2. Click Test. |
| Remote autoanalysis | <ol style="list-style-type: none"> 1. Before you specify the destination, verify that you can connect to the computer on which you intend to store the .fsa sample files and that the shared folder you previously created (see page 72) is accessible from that computer. <p>Note: To optimize data sharing, store files on a non-GeneMapper® ID-X Software computer (see page 83).</p> <ol style="list-style-type: none"> a. On the Data Collection computer, select the Start menu, then select Run. b. Enter the names of the computer on which you intend to store the .fsa sample files and the shared folder, using the following format: <code>\\<server name>\<shared folder name></code> For example: <code>\\GMID_PC\Remote_AutoAnalysis</code>  <ol style="list-style-type: none"> c. Click OK. If Windows Explorer opens to the shared folder on the designated computer, the computer and shared folder are accessible. Close Windows Explorer, then proceed to step 2. If Windows Explorer does not open to the shared folder on the designated computer, either the computer is not accessible, or the shared folder is not set up correctly. Troubleshoot your network, computer, or shared folder problem before proceeding to step 2. 2. In the Destination tab (see page 76), enter the names of the computer on which you intend to store the .fsa sample files and the shared folder, using the format shown in step 1b. 3. Click Test to test the location path name connection. If the test: <ul style="list-style-type: none"> • Passes – The message “Test succeeded” is displayed. • Fails – The message “Could not make the connection. Please check that the Path Name is correct” is displayed. Make sure you correctly entered the computer and shared folder names, then troubleshoot your network or computer before continuing. |

Creating Sample and Name Plates

In the naming tab, provide information about the samples:

1. Click the **Naming** tab.
2. From the Format list in the Sample File Name Format section (top), select a series of fields to create unique sample names, for example, Well Position, Sample Name, and Capillary Position.



3. From the Format list in the Run Folder Name Format section (bottom), select **Plate Name**.

IMPORTANT! To ensure that the run folder includes at least one allelic ladder for genotyping, the run folder must be grouped by plate, not by run.

When you select Plate Name, an INVALID NAME message is displayed. Disregard the message. All samples from the plate are saved into one run folder.


Note: Sample name, run folder name, and path name, *combined*, can contain no more than 250 characters.

Setting Automated Processing

1. Click the **Automated Processing** tab.

Note: This tab is available only if you selected GeneMapperIDX-*<computer name>* in [step 2 on page 74](#).

2. Select **Only when the result group is complete** to analyze after all the samples that use the same results group have been processed and click **OK** to save the results group.



Select **Only when the result group is complete**

When you are finished, go to Procedure 4 to create a plate record for autoanalysis.

Procedure 4: Creating a Plate Record for Autoanalysis

A plate record specifies the instrument protocol, results group, and analysis parameters for a run of samples. For more information about plate records, see your instrument user guide.

IMPORTANT! Create new plate records after installing or reinstalling the GeneMapper® *ID-X* Software v1.0 or v1.1, and before each autoanalysis run. This procedure is provided in the *GeneMapper® ID-X Software Online Help*.

In the Navigation pane of the Data Collection Software, click the **Plate Manager** icon to display the Plate Editor.

IMPORTANT! Before you can create a new plate record, you must specify or import the size standards, panels, and analysis methods in the GeneMapper *ID-X* Software, then synchronize the Data Collection Software with the GeneMapper® *ID-X* Software (see [page 90](#)).

Defining New Plate Information

1. Click **New** to open the New Plate Dialog box.
2. Enter a name for the plate. The Plate Name is the name assigned to the run folder that resides in the shared folder where the .fsa sample files are stored.
3. (Optional) Enter a description for the plate.
4. Select **GeneMapperIDX-*<computer name>*** (the name of the computer that you intend to have run autoanalysis) in the Application drop-down list.
5. Select **96-well** or **384-well** in the Plate Type drop-down list.
6. Enter names for the owner and operator.
7. Click **OK**.

The 'New Plate Dialog' window contains the following fields and controls:

- 2: Name text input field.
- 3: Description text area.
- 4: Application dropdown menu (currently set to 'None').
- 5: Plate Type dropdown menu (currently set to '96-Well').
- 6: Owner Name and Operator Name text input fields.
- 7: OK and Cancel buttons at the bottom right.

Describing the Sample

For each row in the Plate Editor (see [page 82](#)), describe the sample:

1. Enter the name of the sample.

After you make selections for the first row, you can select the row, then press **Ctrl+D** to auto-populate (fill-down) the remaining rows.
2. (Optional) Enter a comment about the sample. This information is displayed in the Comment field of the Samples table in the GeneMapper® ID-X Software.
3. Enter a priority number to specify the sequence in which the samples are run.
4. Specify sample type:
 - a. Select **Allelic Ladder** as the sample type for the appropriate wells (you need at least one allelic ladder per plate to allow genotyping).
 - b. Select **Positive Control** or **Negative Control** for your control samples to allow an automated concordance check to ensure that the:
 - Positive control produces the expected profile.
 - Negative control does not contain peaks called above your peak amplitude threshold.

- c. Select **Sample** for all other samples.
5. Select a size standard from the drop-down list.

IMPORTANT! If the expected size standard, panel, and analysis method are not listed, you may need to synchronize the Data Collection Software and the GeneMapper® ID-X Software. See [page 90](#).

6. Select an applicable AmpFℓSTR® kit and panel from the drop-down list.
7. Select an analysis method from the drop-down list.
8. (Optional) Enter text in the user-defined text columns.
9. Select the results group you created using the procedures on [page 72](#).
10. Select the instrument protocol you created using the procedures on [page 70](#). Click **OK**.

Plate Editor

For more information on the Plate Editor, see your instrument user guide. For information on performing autoanalysis and using the Autoanalysis Manager, see the *GeneMapper® ID-X Software Online Help*.

Recommended Data Processing Tasks

Optimizing Data Sharing

Overview If you store .fsa sample files on a GeneMapper® *ID-X* Software computer (ie, the computer that creates the project AND stores the sample files), only the computer that created the project can reanalyze the project. As a result, any other computer must reassociate the files to reanalyze the project.

Additionally, if you export/import reference data and projects following installation but do not reassociate .fsa sample files, you can view results and make or edit allele calls, but you cannot reanalyze the projects.

To reanalyze projects:

1. Move .fsa sample files to a non-GeneMapper® *ID-X* Software computer to optimize file sharing.
2. Reassociate the .fsa sample files with the imported projects (**Edit ► Define New Sample Path**) to view raw and EPT data and reanalyze the imported projects on your client computer.

[Figure 5](#) shows the recommended configuration for optimizing data sharing.




Figure 5 Autoanalysis configuration to optimize data sharing

The following sections describe data optimization in detail. For additional details on how to perform each procedure, see the *GeneMapper® ID-X Software Online Help*.


Setting Up a Shared Folder (Remote Autoanalysis Only)

Remote autoanalysis of data requires a shared folder in which to store .fsa sample files. To make the sample data files accessible to all client computers, set up the shared folder on a computer that is not running the GeneMapper® ID-X Software.

On the computer on which you intend to store the .fsa sample files:

1. Select **Start ▶ My Computer**, or double-click **My Computer** on your desktop, then double-click the drive on which you want the shared folder.
2. Select **File ▶ New ▶ Folder**.


3. Name the folder (for example, **Remote_Autoanalysis**).
4. Right-click the new folder, then select **Properties**.
5. In the Sharing tab, select **Share this folder**.
6. (Optional) Enter a comment for the folder.
7. Click **OK**.




Setting Permissions and Security for the Shared Folder (Windows XP Only)

If you use a Microsoft® Windows® XP® operating system, modify the permissions and security settings of the shared folder:


1. Right-click the shared folder, then select **Properties**.
2. Select the **Sharing** tab, then click **Permissions**.
3. In the Permissions for <shared folder name> dialog box, select the check box for **Full Control** in the Allow column.
4. Click **OK**.



5. Select the **Security** tab, then click **Add**.
6. In the Select Users, Computers, or Groups dialog box, enter **Guest** in the “Enter the object names to select” field.
7. Click **OK**.



8. In the Group or user names list box in the Security tab, configure guest access:
 - a. Select **Guest**.
 - b. Select the check box for **Full Control** in the Allow column.
 - c. Click **OK**.



The shared folder is now configured for use.

Importing and Reassociating Data Files

Depending on your current network configuration, you can import and reassociate the files in two ways:

- **If the .fsa sample files are stored on a non-GeneMapper® ID-X Software computer other than a user's local computer** – The software saves the full networked path location of the .fsa sample files (for example: //server/folder/folder*.fsa). Because the path location is available to GMIDX, any computer in the network can access the project and reanalyze it without reassociating the files.

- If the .fsa sample files are stored on a user's local computer (that is, the computer that creates the project AND stores the sample files) - The software saves the path location of the .fsa sample files, but not the name of the computer on which the files reside (for example: C:/folder/folder*.fsa). The local computer can reanalyze the files because GMIDX looks on the local computer for the .fsa sample files. However, the name of the computer is not available to other computers running GMIDX so they cannot reanalyze projects without first reassociating files (**Edit ► Define New Sample Path**) (see [Figure 6](#) for file storage comparisons).





Figure 6 Storing .fsa sample files: non-Genemapper® ID-X Software computer vs. local computer

Synchronizing GMID-X Software with the Data Collection Software

After you install the GeneMapper® *ID-X* Software, and any time you create or change analysis methods, panels, or size standards in the GeneMapper® *ID-X* Software, you must synchronize the new data with the Data Collection software before you can select the updated items in the Data Collection Software plate record.

To synchronize the GeneMapper® *ID-X* Software with the Data Collection Software:

1. Select **File ▶ Exit** to close the GeneMapper® *ID-X* Software.
2. Verify that the Data Collection Software is running on the Data Collection computer. If not, start the Data Collection Software. For information on starting the Data Collection Software, see [“Starting the Data Collection Software on the 31xx Data Collection Computer” on page 26](#).
3. On your desktop:
Double-click  **GeneMapper *ID-X* v1.0 (or v1.1)**,
or
Select **Start ▶ All Programs ▶ Applied Biosystems ▶ GeneMapper ▶ GeneMapper *ID-X* v1.0 (or v1.1)** and log in.



Installing the Generic Updater Software

A

This appendix covers:

- Installing the Generic Updater Software 92
- Removing the Generic Updater Software 94

Installing the Generic Updater Software

The Generic Updater software is an add-on to the 3730 and 3130 Data Collection 3.0 software. Use the Generic Updater to define custom plates and write sample files to be used with the GeneMapper® *ID-X* Software v1.0 or v1.1.

If you have either 3730 or 3130 Data Collection 3.0 software installed on your computer, and you do not intend to install the GeneMapper® *ID-X* Software v1.0 or v1.1 on the same computer for either local or remote autoanalysis, you can install the Generic Updater software on the Data Collection computer.

The generic Plate Editor, residing in the Plate Manager function, has columns for GMID-X plate data that can be edited as text fields. The software provides you with the option to add generic or non-automated GMID-X plate definitions and sample files to your Data Collection software installation.

To review the requirements for installing the Generic Updater software, see [“Installation Options” on page 3](#).

1. Insert one of the following CDs into the Data Collection computer CD drive:
 - *GeneMapper® ID-X Software Version 1.0 (or Version 1.1) Full Install* CDs - Disk 1
 - *GeneMapper® ID-X Software Version 1.0 (or Version 1.1) Client Install* CD
2. Ensure that the Data Collection software is running (see [page 26](#)).
3. On the splash screen, click **Software Extras**, then on the next screen click **Generic Updater**.
4. On the following screen, click **Install Generic Updater**.
5. Click **Browse**, then select a location for the Generic Updater zip files and click **Unzip** to extract the files.
6. Remove the *GeneMapper® ID-X Software v1.0 (or v1.1)* CD from the CD drive.

7. Navigate to the extracted folder on your hard drive. Open the Deployment-Generic folder, then double-click **RegisterGMID-XGeneric.exe**.
8. Click **OK** when you see the message “Application GeneMapper® successfully registered in JNDI.”
9. Stop the Data Collection software and restart your computer.

Start the Data Collection software and open the Viewer window, then use the Plate Manager to create new GeneMapper *ID-X* Generic plate definitions.

For instructions on creating and using GeneMapper® *ID-X* Generic plates, see the User Guide for your specific instrument.

Resolving Installation Problems

If **RegisterGMID-XGeneric.exe** reports an installation problem:

1. Close the window and check the **JNDIRegistrationGMID-XGeneric** log file in the extracted folder.
2. Check your Data Collection installation to verify that the Data Collection software is running.
3. Try running **RegisterGMID-XGeneric.exe** again.

Removing the Generic Updater Software

IMPORTANT! If you decide to install the GeneMapper *ID-X* Software directly onto the Data Collection computer or as a remote autoanalysis configuration, you must first remove the Generic Updater Software from the Data Collection computer.

To uninstall the Generic Updater Software:

1. Ensure that the Data Collection Software is running (see [page 26](#)).
2. Navigate to the extracted folder on your hard drive.
3. Open the Deployment-Generic folder, then double-click **UnRegisterGMID-XGeneric.exe** to remove the GeneMapperIDX-Generic option from the Data Collection Software.
4. Click **OK** when you see the message “Application GeneMapper® successfully unregistered in JNDI.”
5. Verify that the Generic Updater Software was removed:
 - a. Open the Data Collection Viewer window.
 - b. Open the Plate Manager.
 - c. Verify that the GeneMapperIDX-Generic option is not available as an application choice.
6. Delete the Generic Updater Software files from your hard drive.

Resolving Uninstallation Problems

If **UnRegisterGMID-XGeneric.exe** reports a problem and the Generic Updater software is still displayed in the Data Collection Viewer Plate Manager:

1. Check the **JNDIRegistrationGMID-XGeneric** log file in the extracted folder.
2. Check your Data Collection installation to verify that the Data Collection software is running.
3. Run **UnRegisterGMID-XGeneric.exe** again.



Appendix A Installing the Generic Updater Software

Removing the Generic Updater Software



Troubleshooting the Installation

This appendix covers:

| | |
|-----------------------------------|-----|
| ■ Troubleshooting Tasks | 98 |
| ■ Troubleshooting Tools | 105 |

B



Troubleshooting Tasks

Disabling Firewall Protections

During a full GeneMapper® *ID-X* Software installation, the security firewall settings may prevent the database from installing. You may see Windows security alerts that prompt you to “keep blocking” or “unblock” the security firewalls. Alternatively, it may appear that the installation was complete but when you start the software, the database may not launch.

1. Before installing GMID-X v1.0 or v1.1, disable the firewall settings: go to **Start ▶ Control Panel ▶ Windows Firewall** to disable the settings.
2. At or during installation, if you are prompted with Window firewall alerts, click **Unblock** and proceed with installation.

If the software application doesn't launch:

- Uninstall the software (see [“Uninstalling the GMID-X Software” on page 100](#)).
 - Disable the firewall settings by going to **Start ▶ Control Panel ▶ Windows Firewall**,
or
Disable the firewall settings for your particular Windows operating system. Go to **Start ▶ Control Panel ▶ Security Center** for Windows online support.
3. Reinstall the software for your specific ID-X configuration.
 4. Reset firewall.

Canceling an Installation

DO NOT cancel an installation before the installation process is finished. If you do so, or if the installation otherwise stalls, you may need technical support to perform an uninstall. To uninstall your installation, see [“Uninstalling the GMID-X Software” on page 100](#). If uninstallation fails, contact your local Applied Biosystems technical support organization.

General Installation Solutions

Note: You may need to send the GMIDX install and MSI log files to Applied Biosystems Technical Support if the troubleshooting procedures in this appendix do not solve the installation problem.

If you have trouble installing the GeneMapper® *ID-X* Software v1.0 or v1.1 on your computer:

| Task | See page |
|--|---------------------|
| 1. Obtain the installation log file. | 99 |
| 2. Check the database installation. | 100 |
| 3. If necessary, uninstall the software. | 100 |
| 4. If necessary, have a knowledgeable technician run the clean-up utility. | 102 |

Obtaining the Installation Log File

To obtain the installation log file, go to:

`x:{installdir}\AppliedBiosystems\GeneMapper ID-X` where *x* is the drive on which the GeneMapper *ID-X* Software is installed and *{installdir}* is the user-selected installation path.

- You will find 2 logs for the full installation:
 - GMIDXInstall.log
 - GMIDX_MSI.log
- You will find 2 logs for the client installation:
 - GMIDXClientInstall.log
 - GMIDXClient_MSI.log



Checking the Database Installation

1. Select **Start ▶ All Programs ▶ Accessories ▶ Command Prompt**.
2. Enter **sqlplus**, then press **Enter**.
3. Enter the user name, **system**, then press **Enter**.
4. Enter the password, **manager**, then press **Enter**.
 - If a “Connected to:” message is displayed, the database is installed correctly.
 - If no “Connected to:” message is displayed, the database is not installed correctly. Complete [Table 6, “Troubleshooting Checklist,” on page 105](#), then contact Applied Biosystems Technical Support.

Uninstalling the GMID-X Software

To uninstall the GeneMapper® *ID-X* Software and the multi-user database:

1. If you are uninstalling a *co-installation*, *Data Collection client* installation, or a *remote autoanalysis* configuration, start the Data Collection software (see [page 26](#)).
2. Close the GeneMapper® *ID-X* Software and all other applications, then restart the computer before proceeding with the uninstall.

IMPORTANT! If you do not restart the computer before proceeding with the uninstall and/or if the uninstall does not complete successfully the first time, perform this procedure a second time. If the uninstall still does not work, have a knowledgeable technician review and proceed with [“Running the Clean-up Utility” on page 102](#) or contact Applied Biosystems Technical Support to proceed.


3. Select **Start ▶ Control Panel**.

4. Double-click **Add or Remove Programs**.

Note: If a message indicates that you need administrator rights or privileges, log off the computer, then log on again as a user with administrator privileges.

Note: If you receive an error message, dismiss the message and proceed to [“Running the Clean-up Utility” on page 102](#).


5. Select the GeneMapper® *ID-X* Software.
6. Click **Change/Remove**. The InstallShield Wizard opens.
7. In the Welcome page, select **Remove**, then click **Next**.
8. At the prompt, click **OK** to verify the uninstall.
 - If you are uninstalling a *Co-installation* **or** *Data Collection client* installation, the message below appears. Be sure the Data Collection Service Console is open, then click **OK** to start the uninstall.



- If you are uninstalling a *remote autoanalysis* configuration, no message is displayed; however, the Data Collection Software must be running and the Service Console open on the remote Data Collection computer to uninstall a remote autoanalysis configuration.

9. When the Maintenance Complete window opens, finish the uninstallation process:

IMPORTANT! If you uninstall a *co-installation* or *Data Collection client* installation, stop the Data Collection Software before restarting your computer by clicking **Stop All** in the Service Console.



- a. Select **Yes, I want to restart my computer now.**
- b. Click **Finish.**

Running the Clean-up Utility

Run the clean-up utility under only the following conditions:

- You have attempted the uninstall procedure using the **Add/Remove** Program uninstall option,
and
- The GMID-X v1.0 or v1.1 software files were not removed or were only partially removed.

Have only a knowledgeable technician run the *msicuu2.exe* utility to clean up the remaining files.

To perform the clean-up procedure:

1. Insert the GeneMapper® *ID-X* Software CD, then go to **DVD/CD-RW Drive ▶ Utilities**.
2. Run *msicuu2.exe* to install the utility.
3. Launch the utility from the desktop icon or from **Start ▶ Programs ▶ Windows** to install.
4. Select the GeneMapper *ID-X* application from the list provided, then click **Remove**.
 - a. Check to ensure that the *gmidx* install folder and *oraclegmidx* folder under your {installdir} are deleted.
 - b. Check to ensure that the Oracle folder in the Program Files folder is deleted.
5. Delete variable strings:

IMPORTANT! Steps 5 and 6 must be performed *carefully* and *correctly*. If they are not performed correctly, other applications could stop working, and the stability of your system could be at risk. If you do not feel comfortable performing these steps, please contact your system administrator to perform them for you.

- a. Right-click on **My Computer**, then select **Properties**.
- b. Click the **Advanced** tab.
- c. Click the **Environment Variables** button.
- d. In the *System variables* Variable pane, delete the words *ORACLE_HOME* and *ORACLE_HOSTNAME*.
- e. In the *System variables* Variable pane, select the *Path* variable and click the **Edit** button.
- f. Remove the path to *oraclegmidx* from the path environment variable.
- g. Click **OK** to exit the Environment Variables dialog box and **OK** to exit the System Properties dialog box.



6. Go to **Start ▶ Run** and enter **regedit**. In the left-pane of the Registry Editor, open the following folders sequentially:
HKEY_LOCAL_MACHINE ▶ SYSTEM
▶ CurrentControlSet ▶ Services.
 - a. In the Services folder remove all Oracle folders containing **IFA** letters only.
 - OracleIFAGMIDXCMan
 - OracleIFADataGatherer
 - OracleIFAGMIDXTNSListener
 - OracleIFAAgent
 - OracleIFAClientCache
 - OracleIFAGMIDXCMAAdmin
 - OracleServiceIFA
 - OracleJobSchedulerIFA
 - b. Log off and reboot the system.

Troubleshooting Tools

Complete the Troubleshooting Checklist and Interdepartment Forms before contacting Applied Biosystems Technical Support for assistance.

Checklist

Table 6 Troubleshooting Checklist

| Check | Information for Technical Support |
|-------|--|
| | Summarize the problem: |
| | Can you repeat the problem? If yes, list the steps that you perform: 1. 2. 3. 4. 5. 6. 7. |
| | Applied Biosystems personnel that you have contacted: <input type="checkbox"/> Field Applications Specialist <input type="checkbox"/> Field Service Engineer <input type="checkbox"/> Technical Support <input type="checkbox"/> Sales Representative <input type="checkbox"/> Order Administration <input type="checkbox"/> Other |



Appendix B Troubleshooting the Installation

Troubleshooting Tools

| Check | Information for Technical Support |
|-------|--|
| | Computer specifications: <ul style="list-style-type: none">• Operating system:• Version:• Processor:• Memory:• Hard disk space:• Hard disk configuration: |
| | Software installed: <ul style="list-style-type: none"><input type="checkbox"/> Data Collection Software version:<input type="checkbox"/> Status of Data Collection services:<input type="checkbox"/> GeneMapper <i>ID-X</i> Software version:<input type="checkbox"/> Other Applied Biosystems software: |
| | Computer login information: <ul style="list-style-type: none">• User privileges:• Local or networked domain: |
| | Software configuration installed: <ul style="list-style-type: none"><input type="checkbox"/> Co-installation<input type="checkbox"/> Remote autoanalysis<input type="checkbox"/> Stand-alone<input type="checkbox"/> Multi-user database |
| | Instrument and instrument computer information: <ul style="list-style-type: none">• Model:• Data Collection Software version:• Status of Data Collection services:• Other Applied Biosystems software:• Capillary length:• Capillary lot number:• Run module:• Dye set: |
| | Chemistry kit or reagent, with version number: |

| Check | Information for Technical Support |
|-------|--|
| | <p>Be prepared to send to Technical Support:</p> <ul style="list-style-type: none">• Installation log file (page 99)• Exported panels• Exported bins• Exported size standard definition• Exported analysis method• Sample (.fsa) files• GMIDXInstall.log or GMIDXClientInstall.log• PanellImportLog.txt• Printed results |

Interdepartmental Forms

Answering questions in the following forms helps manage complex troubleshooting procedures. Your answers will be shared among support groups, product groups and software development teams at Applied Biosystems.

To address any issues that may arise when installing or running the software, please answer the applicable list of questions, and collect the necessary files before contacting Technical Support.



For Installation Issues

Collect the following files and complete the questions below:

For full installation:

- GMIDXInstall.log
- GMIDX_MSI.log

For client installation:

- GMIDXClientInstall.log
- GMIDXClient_MSI.log

Do the computer and operating system meet the required computer specifications? (See [Chapter 2, "Installation Requirements,"](#) on [page 7.](#))
List any potential discrepancies.

Did you receive any error messages or alerts? If so, what did they say? (if possible, take screen captures.) At what point during the installation process did these errors/alerts occur?

How many installation attempts have you tried?

Are you logged in to the local computer?

Is the TCP/IP configured?

Do you have administrative privileges?

Is this setup stand-alone, networked, or connected to a Genetic Analyzer?



**For Post-
Installation
Related Issues**

| |
|---|
| <p>Do the computer and operating system meet the required computer specifications? (See Chapter 2, “Installation Requirements,” on page 7.) List any potential discrepancies.</p> |
| <p>Problem description:</p> |
| <p>Did you receive any error messages or alerts? If so, what did they say? (If possible, take screen captures.) At what point did these errors/alerts occur?</p> |
| <p>Is the problem reproducible?</p> |
| <p>Steps taken prior to encountering or to reproduce the problem:</p> <ol style="list-style-type: none">1.2.3.4.5. |





Appendix B Troubleshooting the Installation

Troubleshooting Tools



Glossary of Terms

This appendix covers:

| | |
|---------------------------------|-----|
| ■ Definition of Terms | 112 |
|---------------------------------|-----|



Definition of Terms

| | |
|---|---|
| Client Installation | A computer installed with the GeneMapper® <i>ID-X</i> Software only (software with no database). A client installation computer connects to a full installation computer to access the multi-user database. A client can be installed on a Data Collection (DC) computer, or it can be installed on a non-Data Collection (non-DC) computer in either a remote autoanalysis or stand-alone configuration. |
| Co-Installation (31xx) | A full installation (database and GeneMapper® <i>ID-X</i> Software) on a Data Collection computer that can be set to communicate locally with the Data Collection software for the collection and processing of .fsa sample files, and the analysis of data. Not recommended in a multiuser database configuration because the multi-user host is limited to a one-client-connection. |
| Data Collection Client Installation | An installation (GeneMapper® <i>ID-X</i> Software only) on a Data Collection computer that can be set to communicate locally with the Data Collection software for the collection and processing of .fsa sample files, and the analysis of data. Connects to a multi-user database host computer to access, analyze, and share data and analysis settings. |
| 310 Data Collection (DC) Computer | A stand-alone computer with Data Collection software that controls the Genetic Analyzer instrument and generates .fsa sample files used for analysis by the host or client computer GeneMapper® <i>ID-X</i> Software. |
| 31xx or 3730 Data Collection (DC) Computer | A computer with Data Collection software that controls the Genetic Analyzer instrument and generates .fsa sample files used for analysis by the host or client computer GeneMapper® <i>ID-X</i> Software. Can support a co-installation or client installation. |
| Data Collection Software | The software on the Data Collection computer that controls the Genetic Analyzer instrument and generates the .fsa sample files used for analysis. |



Full Installation A computer installed with both the multi-user database and GeneMapper *ID-X* Software.

A full installation computer can be part of a multi-user database configuration, or it can be a data collection (DC) or non-data-collection (non-DC) computer having no interaction with other computers. A full installation can occur in any of the following configurations, depending on the function of the computer:

- Co-installation
- Remote Autoanalysis
- Stand-alone

Genetic Analyzer The capillary electrophoresis instrument, controlled by the Data Collection software, which generates .fsa sample files.

GMID Abbreviation for GeneMapper® *ID* Software.

GMID-X Abbreviation for GeneMapper® *ID-X* Software.

Local Area Network A communications network that serves users within a confined geographical area.

Local Autoanalysis A co-installation (database and GeneMapper® *ID-X* Software) or client installation (GeneMapper® *ID-X* Software only) that is set up on a Data Collection computer to locally communicate with the Data Collection software for the collection and processing of .fsa sample files, and the analysis of data.

Multi-User Database Host Computer A computer on which data is stored on a central, multi-user database. Commonly referred to as the host computer, from which other client computers (computers without a database) can access, analyze, and share data. Can be installed for a remote autoanalysis or stand-alone configuration.

Non-Data Collection (Non-DC) Computer A non-data collection computer that may or may not be connected to a local area network (LAN). Can be a client installation or full installation, with a stand-alone or remote autoanalysis configuration.





**Stand-alone
Configuration**

A full installation (database and GeneMapper® *ID-X* Software) or client installation (GeneMapper® *ID-X* Software with no database) on a non-Data Collection computer that does not remotely communicate with the Data Collection computer. Can be the database host or client computer in a multi-user database configuration.

**Upgrade
Installation**

A software upgrade from GeneMapper® *ID-X* Software v1.0 or v1.0.1 to GeneMapper® *ID-X* Software v1.1. The upgrade can be performed on a client or full installation computer that already has the v1.0 or v1.0.1 software installed. However, to upgrade on a client computer, you must also upgrade the computer that hosts the full multi-user database installation.

**Remote
Autoanalysis**

A full or client installation on a non-Data Collection computer that is set up remotely to communicate with the Data Collection computer (data collection software) for the automatic acquisition and analysis of data. Can be the database host or client computer in a multi-user database configuration.

Numerics

31XX Data Collection computer
 client installation 46

A

administrator privileges 34, 44
 custom profiles 69
 logging on 25, 34, 44, 101
allelic ladder, Results Group 72
analysis parameters, sample run 80
Applied Biosystems
 contacting xi
 customer feedback on documentation x
 Information Development department x
 Technical Support xi
autoanalysis
 autoanalysis manager definition 64
 autoanalysis manager for remote
 autoanalysis 68
 configuration components 64
 instrument protocol 70
 process 62
 running 82
 setup, following installation 68
 user accounts 69

B

backing up data 16

C

changing installation configurations 67
client installation 8
 computer requirements 11

 connecting to new host 56
 DOS window 51
 installation time 12
 installing client 46
 local area network 44
 log file 99
 monitor requirements 11
 on 31XX Data Collection computer 46
 on Data Collection computer 4, 66
 on non-DC computer 4, 46
 operating system requirements 11
 optimizing data sharing 83
 troubleshooting 98

co-installation

 on Data Collection computer 66
 uninstalling GMID-X v1.0 software 100
 upgrade on Data Collection computer 4

computer requirements

 administrator privileges 25, 34, 44, 101
 client installation 11
 Ethernet, full installation 10
 full installation 9
 monitor, full installation 9
 operating system, full installation 10
 user accounts 34, 44
 verifying setup 13
 verifying user accounts 34, 44

computer setup

 verifying 13

connecting to new host, client 56

conventions

 for describing menu commands viii
 IMPORTANT! ix
 Notes ix
 user attention words ix

- creating new password for GMID-X v1.0 55
- creating user accounts 59, 69
- custom profiles, setting privileges 69
- customer feedback, on Applied Biosystems documents x

D

- Data Collection client installation
 - uninstalling GMID-X v1.0 software 100, 101
- Data Collection computer
 - client installation 4, 66
 - co-installation 66
 - co-installation upgrade 4
 - computer name 29, 50
 - enabling autoanalysis 26
 - Generic Updater Software 92
 - installation on 4
 - local autoanalysis 38
 - synchronizing software 90
- Data Collection Software
 - enabling autoanalysis 26
 - fsa sample files 38, 48
 - Results Group 72
 - Results Group and Plate list 72
 - service console 27, 41, 51, 102
 - starting manually 27
 - starting on 31XX data collection
 - computer 26
 - stopping 41, 51, 102
 - verifying error-free startup 28
- data options, instruments 14
- data sharing
 - optimizing 83, 84
 - remote autoanalysis 83
- destination folder, selecting 39, 40, 50
- disclaimer, license ii
- documentation, related x
- domain name
 - user accounts 34, 44
 - verifying 34, 44

- DOS window
 - client installation 51
 - during full installation 41
- dye set 70

E

- enabling autoanalysis
 - co-installation 26
 - remote 26
- Ethernet requirements, full installation 10
- exporting data 16
 - automatically 18
 - error messages 20
 - manually 18, 23
 - on non-DC computer 4
 - projects 59

F

- folder privileges
 - client installation 52
 - resetting permissions 52
- fsa sample files
 - Data Collection Software 38, 48
 - in Results Group 72
 - reassociating 83, 89
 - remote autoanalysis 38, 49
 - shared folder 80, 84
 - sharing 77, 83, 84
 - storing 77, 83, 84
- full installation
 - computer requirements 9
 - DOS window 41
 - Ethernet requirements 10
 - installation time 10
 - log file 99
 - monitor requirements 9
 - new 8
 - on non-DC computer 4, 21, 23
 - troubleshooting 98

G

GeneMapper® Software

- login 55
- logout 58
- registration code 54

Generic Updater Software

- installing 92
- on Data Collection computer 92
- uninstalling 94

genetic analyzer

- data options 14
- instrument compatibility 14
- operating system compatibility 65
- software compatibility 65

GMID export utility 16

GMID-X login 55

GMID-X logout 58

GMID-X v1.1

- upgrade installation 31

H

host computer name

- enabling autoanalysis 26

I

importing projects 59

Information Development department, contacting x

installation problems, troubleshooting 99

installation time

- client installation 12
- full installation on data collection
computers 10
- full installation on non-data collection
computer 10

installer

- client 8, 52, 92
- resetting permissions 52
- specifying folder privileges 52
- upgrade 8

installing software

- remote autoanalysis, non-DC
computer 39

instrument compatibility

- operating system 65
- software 65

instrument protocol 80

- creating 70
- dye set 70
- for autoanalysis 70
- name 70
- protocol name 70
- run module 70
- run type 70

L

license disclaimer ii

local area network

- client installation 44
- login requirements 44

local autoanalysis

- 31XX Data Collection computer 38
- Data Collection computer 38

log file

- client installation 99
- full installation 99
- obtaining 99
- troubleshooting 99

login

- GeneMapper® Software 55
- GMID-X 55
- password 55

logout

- GeneMapper® Software 58
- GMID-X 58

M

minimum computer requirements

- software installation 13
- System Properties dialog box 13

monitor requirements
 client installation 11
 full installation 9

MSDSs, obtaining xi

multi-user database
 autoanalysis components 64
 dating sharing 83
 obtaining host name 30
 optimizing data sharing 84

N

non-DC computer
 client installation 4, 21, 23, 46
 exporting data 4
 full installation 4, 21, 23
 installation on 4, 21, 23
 NTFS permissions, setting manually 52

O

operating system requirements
 client installation 11
 full installation 10
 optimizing data sharing for autoanalysis 84

P

password
 reusing 31
 password, creating
 for GMID-X v1.0 55
 performing autoanalysis 82
 permissions, setting 52
 plate dialog 80
 Plate Editor, for autoanalysis 80
 plate name 80
 plate records, specifying
 analysis parameters, for sample runs 80
 for autoanalysis 80
 instrument protocol 80
 Results Group 80

plate records, synchronizing software 80, 90

R

reassociating fsa sample files 83, 89
 registration code, software 54
 remote analysis
 optimizing data sharing 83
 remote autoanalysis
 Autoanalysis Manager 68
 changing installation configurations 67
 Data Collection computer name 29, 50
 data sharing 83
 installing on non-DC computer 39
 saving fsa sample files 38, 49
 uninstalling GMID-X v1.0 software 100
 resetting folder permissions 52
 Results Group and Plate list
 allelic ladder 72
 Data Collection Software 72
 Results Group, creating for autoanalysis 72, 80
 reusing existing password 31
 run folder 80
 run module 70
 run type 70

S

sample run, analysis parameters 80
 service console 41, 51, 102
 starting Data Collection Software 27
 software installation
 minimum computer requirements 13
 software registration code 54
 synchronizing software 80, 82, 90

T

- technical support
 - checklist [105](#)
 - contacting [xi](#)
 - login,client installation [44](#)
 - needed information [105](#)
- training, information on [xi](#)
- troubleshooting
 - disabling firewall settings [98](#)
 - obtaining log file [99](#)
 - technical support checklist [105](#)
 - uninstalling software [98](#)

U

- uninstalling Generic Updater software [94](#)
- uninstalling GMID v3.x software [16](#), [24](#)
- uninstalling GMID-X v1.0 software
 - changing installation configurations [67](#)
 - co-installation [100](#), [101](#)
 - Data Collection client installation [100](#), [101](#)
 - multi-user database [100](#)
 - remote autoanalysis [100](#)
- upgrade installation [8](#)
 - GMID-X v1.1 [31](#)
 - reusing password [31](#)
- user accounts
 - connecting to new host [56](#), [69](#)
 - creating [59](#), [69](#)
 - for autoanalysis [69](#)
 - full installations [34](#), [44](#)
 - multi-user database [56](#)
 - verifying domain name [34](#), [44](#)
- user attention words, definitions [ix](#)
- user guide
 - attention words [ix](#)
 - audience [viii](#)
 - text conventions [viii](#)

Worldwide Sales and Support

Applied Biosystems vast distribution and service network, composed of highly trained support and applications personnel, reaches 150 countries on six continents. For sales office locations and technical support, please call our local office or refer to our Web site at **www.appliedbiosystems.com**.

Applera is committed to providing the world's leading technology and information for life scientists. Applera Corporation consists of the Applied Biosystems and Celera Genomics businesses.

Headquarters

850 Lincoln Centre Drive
Foster City, CA 94404 USA
Phone: +1 650.638.5800
Toll Free (In North America): +1 800.345.5224
Fax: +1 650.638.5884

07/2008